

Sep 20, 2023

Rapid and robust cloning of sgRNA expression plasmids

DOI

dx.doi.org/10.17504/protocols.io.4r3l229ojl1y/v1

Ghanem Kassem¹, Michael Böttcher¹

¹Medical Faculty of the Martin Luther University Halle-Wittenberg

Atlas of Variant Effects ...



Ghanem Kassem

Universitätsklinikum Halle (Saale)

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l229ojl1y/v1>

Protocol Citation: Ghanem Kassem, Michael Böttcher 2023. Rapid and robust cloning of sgRNA expression plasmids.
protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l229ojl1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 12, 2023

Last Modified: September 20, 2023

Protocol Integer ID: 87667

Keywords: sgRNA Cloning, T4 Ligation, Heat shock Transformation, robust cloning of sgrna expression plasmid, sgrna expression plasmid, spcas9 sgrna, sgrnas expression vector, sgrna, step procedure for cloning, ligating oligonucleotide, cloning, robust cloning, digested plasmid, heat shock transformation procedure to clone, t4 ligation, using oligo annealing

Abstract

In this lab protocol, we will outline a step-by-step procedure for cloning of spCas9 sgRNAs using oligo annealing and T4 ligation. This protocol will guide you through the crucial steps of designing, annealing, ligating oligonucleotides and digested plasmids, as well as heat shock transformation procedure to clone sgRNAs expression vectors.

Attachments



[sgLenti vector for C...](#)

16KB

Protocol materials

 DNA Gel Loading Dye (6X) **Thermo Fisher Scientific Catalog #R0611**

 Macherey-Nagel™ NucleoSpin™ Gel and PCR Clean-up Kit **Fischer Scientific Catalog #11992242**

 NucleoSpin Plasmid Transfection-grade, Mini kit for ultrapure plasmid DNA **Macherey-Nagel Catalog #740490.250**

 T4 DNA Ligase Buffer (10X) **Thermo Fisher Catalog #B69**

 AarI (2 U/μL) **Thermo Fisher Catalog #ER1581**

 T4 DNA Ligase (5 U/μL) **Thermo Fisher Catalog #EL0011**

Troubleshooting



Safety warnings

- ⚠ Wear gloves and a UV filter mask and ensure sleeves are completely covering wrists to avoid direct UV light exposure

Before start

Modify the ssOligo overhangs and Sanger sequencing primer according to the plasmid that you are using.



Oligo Design

- 1 Order **TOP** and **BOTTOM** strand oligos, where the TOP sequence is the desired 20 nt sgRNA 'spacer' sequence and the BOTTOM sequence is the reverse complement of the TOP strand. After annealing of both oligos, the 5' overhangs should result in compatible sticky ends for T4 ligation into the target vector. Shown are example sticky ends to clone into the vector sgLenti (Addgene #105996). Alternative sgRNA expression vectors may require different sticky ends.

5'-**TTGG**NNNNNNNNNNNNNNNNNNNNNN -3'
3'-NNNNNNNNNNNNNNNNNNNNNN**CAAA**-5'

TOP strand oligo: 5'-**TTGG**NNNNNNNNNNNNNNNNNNNNNN-3'

BOTTOM strand oligo: 5'-**AAAC**NNNNNNNNNNNNNNNNNNNNNN-3'

Oligo Annealing

1h 40m

- 2 Set up the annealing reaction:

10m

🧴 1 µL Top strand oligo [M] 100 micromolar (µM)

🧴 1 µL Bottom strand oligo [M] 100 micromolar (µM)

🧴 1 µL 🧬 T4 DNA Ligase Buffer (10X) **Thermo Fisher Catalog #B69**

🧴 7 µL ddH2O

- 3 Anneal oligos by running the following thermocycler program:

1h 30m

1. 🔥 95 °C for 5 min

2. Ramp down at 🔥 0.1 °C /sec from 🔥 95 °C to 🔥 25 °C

Plasmid Restriction Digest

1d

- 4 Set up AarI plasmid digest:

10m

🧴 1 µg vector for spCas9 sgRNA expression

🧴 2 µL 10x AarI Buffer

🧴 1 µL 🧬 AarI (2 U/µL) **Thermo Fisher Catalog #ER1581**

🧴 0.4 µL 50x Oligo

ddH2O to 🧴 20 µL

- 5 Incubate in a Thermocycler at 🔥 50 °C 🕒 Overnight

16h



Digested Plasmid Purification by Agarose Gel Electrophoresis

3h

6 Prepare  1 % Agarose gel

30m

7 Add  4 μL of 6x

10m

 DNA Gel Loading Dye (6X) **Thermo Fisher Scientific Catalog #R0611** to

 20 μL of the restriction digestion mix.

8 Load samples on the  1 % Agarose gel and run the electrophoresis at 120 V for 50 min

1h

9 Place the gel in a UV light box

10m

Safety information

Wear gloves and a UV filter mask and ensure sleeves are completely covering wrists to avoid direct UV light exposure

10 Cut the digested plasmid band from the gel and purify the DNA using the

1h

 Macherey-Nagel™ NucleoSpin™ Gel and PCR Clean-up Kit **Fischer Scientific Catalog #11992242**

11 Determine the concentration of the purified plasmid DNA (e.g. NanoDrop One)

10m

Ligation of Plasmid and Annealed Oligos

30m

12 Dilute the annealed oligos 1:200 in ddH₂O

10m

13 Set up the ligation reaction:

 50 ng purified vector

10m



🧪 1 μ L diluted oligos

🧪 1 μ L 10x T4 Ligation buffer

🧪 1 μ L 🔗 T4 DNA Ligase (5 U/ μ L) **Thermo Fisher Catalog #EL0011**

ddH₂O to 🧪 10 μ L

Note

As a control set up a ligation reaction without oligos

14 Incubate for 10 min at 🌡 Room temperature

10m

Heat Shock Transformation of the Ligation Product into Chemically Competent *E. Coli*

1d

15 Take DH5a chemocompetent cells from 🌡 -80 °C storage and thaw on ice for 30 min. Use 1 🧪 50 μ L aliquot for the ligation reaction and 1 🧪 50 μ L aliquot for the control reaction

30m

16 Add 🧪 2 μ L of the T4 Ligation reaction and the control reaction to the cells and mix gently by flicking the tube

17 Incubate on ice for 30 min

30m

18 Incubate the tube containing the cells in a water bath at 🌡 42 °C for 45 sec and place back on ice for 2 min

2m 45s

19 Add 🧪 0.5 mL of LB medium to the cells and incubate in a thermo-block at 🌡 37 °C with shaking at 250 rpm for 1 hr

1h

20 Plate 🧪 100 μ L of the cells on 🌡 Room temperature LB agar plates with [M] 100 μ g/mL Carbenicillin

10m



21 Incubate at 37 °C Overnight

16h



22 After the overnight incubation, count the colonies on the ligation reaction and the control plates. The ligation reaction plate should have a minimum of 10x more colonies.

10m

23 Seal the plates with parafilm and store at 4 °C for up to 2 weeks

Plasmid DNA Purification

1d

24 Pick up to 3 colonies from the plate using a pipette tip and inoculate 5 mL LB medium with 100 µg/mL Carbenicillin in 15 mL cell culture tubes

25 Incubate the cell culture tubes Overnight at 37 °C and 225 rpm in a shaking incubator

16h



26 Extract the plasmid DNA using the
 NucleoSpin Plasmid Transfection-grade, Mini kit for ultrapure plasmid DNA **Macherey-Nagel Catalog #740490.250**

1h 30m

27 Measure the concentration of the purified plasmid DNA (e.g. NanoDrop One)

28 Send the plasmids for Sanger sequencing with the following primer to confirm the correct sequence:
GGCTTGGATTCTATAACTTCGTATAGCA

Note

Provided primer is complementary to the 'sgLenti vector for CRISPR sgRNA expression' provided in the attachments and needs to be adjusted for alternative vectors.