

Jul 21, 2022

Version 5

Rapid extraction of total lipids from microalgae V.5

DOI

dx.doi.org/10.17504/protocols.io.dm6gpr9jdvzp/v5



Ying-Yu Hu¹, Zoe V. Finkel¹

¹Dalhousie University

Marine Microbial Macroecology Lab

Tech. support email: ruby.hu@dal.ca



Ying-Yu Hu

Dalhousie University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpr9jdvzp/v5>

Protocol Citation: Ying-Yu Hu, Zoe V. Finkel 2022. Rapid extraction of total lipids from microalgae. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpr9jdvzp/v5> Version created by **Ying-Yu Hu**

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 06, 2022

Last Modified: July 21, 2022

Protocol Integer ID: 66101

Keywords: lipids, microalgae, Folch solvent, rapid extraction of total lipid, total lipids from miroalgae, extracted lipid, rapid extraction, total lipid, extraction, microalgae, contamination from filtration system, filtration system, filtration, filter debris, folch solvent, mixture of water, centrifugation,

Abstract

In this protocol, total lipids from miroalgae is extracted with a mixture of water and Folch solvent (2:1 chloroform-methanol v/v). Filter and cell debris is commonly removed by filtration, which is laborious and time consuming. It is also the main reason to either cause sample loss and therefore underestimation; or contamination from filtration system and therefore overestimation. We now use centrifugation to remove filter debris, which stays in between the organic phase and water phase. The extracted lipids is dried under N₂ gas flow, stored under -80 °C for further measurement.

Citation

FOLCH J, LEES M, SLOANE STANLEY GH

. A simple method for the isolation and purification of total lipides from animal tissues.
J Biol Chem, 1957, 226, 497-509.

Citation

Axelsson M, Gentili F (2014). A single-step method for rapid extraction of total lipids from green microalgae..
PloS one.

<https://doi.org/10.1371/journal.pone.0089643>

LINK



Guidelines

Biomass requirement

Considering that

- (1) lipids are approximately 10~30% of microalgal dry mass
- (2) the linear range for colorimetric lipid analysis is 4.2 to 80 µg, the low limit of quantitation is 20 µg

The minimum requirement of sample volume for total lipids is calculated as following:

$$V_L = 20 / \text{Chl-a} / (17.3 / 1.1)$$

If both total lipids and phospholipids are expected to be measured, the minimum sample volume needs to be at least doubled.

Protocol materials

☒ Chloroform (HPLC grade) Merck MilliporeSigma (Sigma-Aldrich) Catalog #439142-4L

☒ Methanol Merck MilliporeSigma (Sigma-Aldrich) Catalog #34860

☒ Potassium chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #P3911-500G

Troubleshooting

Safety warnings



Safety information

Operate chloroform in fume hood.

Safety information

Follow the disposal guidelines regarding the halogenated organic waste.



Collect microalgae samples

1 Precombust GFF filter at 450 °C for 04:00:00

4h

2 Rinse forceps with 95% ethanol, air-dry.

Equipment

Filter forceps

NAME

blunt end, stainless steel

TYPE

Millipore

BRAND

XX6200006P

SKU

Note

Wipe-dry forceps can cause carbon contamination of samples.

3 Filter microalgae in liquid media onto precombusted GFF filters, using gentle vacuum pressure (130 mm Hg).

4 Rinse sample with filtered seawater

5 Place sample filters in cryogenic vials

6 Filter blank media (without cells) through precombusted GFF filter as blank.

7 Flash freeze filters and stored at -80 °C



8 Freeze dry before measurement.

Equipment

FreeZone® 2.5 L Benchtop Freeze Dryers

NAME

Labconco®

BRAND

700202000

SKU

9 Follow <Total particulate carbohydrate from microalgae> protocol to hydrolyze the sample.

Hydrolysis treatment can improve the extraction efficiency:

1. Hydrolysis releases bound lipids into easily extractable forms.
2. Acidified water fraction can facilitate separation of the lipid fraction from extraneous protein and other material.
3. Acid can charge phospholipid to optimize extraction.

Prepare glassware

10 Precombust the centrifuge tubes at  500 °C for  06:00:00

6h

11 Precombust pasteur pipets at  500 °C for  02:00:00

2h



Equipment

Disposable Soda-Lime Glass Pasteur Pipets	NAME
5 3/4"	TYPE
Fisherbrand	BRAND
13-678-6A	SKU

12 Rinse caps with 95% ethanol and air-dry prior to use

13 Latex bulbs are required for Pasteur pipets

Prepare reagent

14 Folch solvent (CHCl_3 : MeOH=2:1 v/v)

14.1 Mix two parts of chloroform and one part of methanol in a 1 L amber bottle. Log the volume of each solvent for double checking the ratio.

	A	B
	Chloroform (mL)	
	Methanol (mL)	



Chloroform (HPLC grade) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #439142-4L**



Methanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #34860**

14.2 Attach dispensette to the bottle, mix well.



Equipment

Bottle-top dispenser

NAME

BrandTech Dispensette® S

BRAND

4731330

SKU

14.3 Label bottle with MSDS label.

15 KCl solution ([M] 0.88 %)

Note

1.5 mL per sample

15.1 Weigh the pyrex media bottle and tare.

Equipment

PYREX® Media Bottles

NAME

Corning®

BRAND

1395-100

SKU

15.2 Directly weigh  0.44 g KCl in the bottle.

 Potassium chloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P3911-500G**



15.3 Top bottle with MilliQ water to  50 g

	A	B
	KCl (g)	
	Final (g)	

Prepare for extraction

16 If lipids samples are not processed for carbohydrate, transfer freeze dried samples and blanks into muffled centrifuge tubes

Equipment

Disposable Glass Screw-Cap Centrifuge Tubes^{NAME}

10 mL^{TYPE}

Corning®^{BRAND}

99502-10^{SKU}

Equipment

Polypropylene Screw Caps^{NAME}

Linerless, 15-415^{TYPE}

Kimble Chase^{BRAND}

73805-15415^{SKU}



16.1 Add  100 μ L MilliQ directly onto the sample.

16.2 Freeze at  -80 $^{\circ}$ C  00:10:00

10m

16.3 Remove vials from freezer.

16.4 Purge the dispensette, fill the tubing with solvent before dispensing solvent into sample tube.

16.5 Dispense  2.0 mL Folch solvent into sample tube.

17 If lipids samples have been hydrolyzed for carbohydrate and solvent has already been added, go to the vortex step directly.

Extract lipids

18 Vortex  00:28:00 by using a tube insert.

28m

Equipment

VWR ANALOG VORTEX MIXER

NAME

VWR

BRAND

10153-838

SKU

With tube insert

SPECIFICATIONS

19 Sonicate  00:02:00

20 Vortex  00:30:00 by using a tube insert.

30m

- 21 Prepare one set of precombusted tubes (#T1), label the tubes, cap is not required.

Equipment

Disposable Glass Screw-Cap Centrifuge Tubes ^{NAME}

10 mL ^{TYPE}

Corning® ^{BRAND}

99502-10 ^{SKU}

- 22 Place one pasteur pipet into each tube (#T1)

Equipment

Disposable Pasteur Pipet ^{NAME}

9 inch ^{TYPE}

VWR ^{BRAND}

14672-380 ^{SKU}

- 23 Prepare another set of precombusted tubes (#T2) for extract. Cap the tube to avoid contamination.

Equipment	
Disposable Glass Screw-Cap Centrifuge Tubes	NAME
10 mL	TYPE
Corning®	BRAND
99502-10	SKU

Equipment	
Polypropylene Screw Caps	NAME
Linerless, 15-415	TYPE
Kimble Chase	BRAND
73805-15415	SKU

24

Centrifuge at

 3200 rpm, Room temperature, 00:05:00

5m



Equipment

General-purpose benchtop centrifuge

NAME

IEC CENTRA CL2

TYPE

Thermo

BRAND

00427 0F

SKU

25 Use the pasteur pipet to transfer organic layer to #T2.

Note

If some debris precipitate at the bottom, gentle blow air bubble through pasteur pipet to relocate it to water layer. Avoid transferring water layer and debris into #T2

26 Add  1 mL Folch solvent to the residue.

27 Sonicate  00:02:00

28 Vortex by using a tube insert for about  00:10:00

10m

29 Centrifuge at  3200 rpm, Room temperature, 00:05:00

Equipment

General-purpose benchtop centrifuge	NAME
IEC CENTRA CL2	TYPE
Thermo	BRAND
00427 0F	SKU

30 Use the pasteur pipet to transfer organic layer to #T2.

31 Add  1 mL Folch solvent to the residue.

32 Sonicate  00:02:00

33 Vortex by using a tube insert for about  00:10:00

34 Centrifuge at  3200 rpm, Room temperature, 00:05:00

Equipment

General-purpose benchtop centrifuge	NAME
IEC CENTRA CL2	TYPE
Thermo	BRAND
00427 0F	SKU



35 Use the pasteur pipet to transfer organic layer to #T2.

36 Add  1 mL Folch solvent to the residue.

37 Sonicate  00:02:00

38 Vortex by using a tube insert for about  00:10:00

Collect lipids extract

1h 30m

39 Centrifuge at  3200 rpm, Room temperature, 00:05:00

5m

40 Use the pasteur pipet to transfer organic layer to #T2.

41 Add  1.5 mL KCl solution to T2, vortex and then centrifuge at
 3200 rpm, Room temperature, 00:05:00

5m

Note

Volume of Folch solvent to KCl is about 4 to 1.

42 Turn on heat block to  37 °C , use a thermometer to monitor the actual temperature.

**Equipment****LSE digital dry bath heater**

NAME

Corning

BRAND

6885-DB

SKU

Equipment**Blocks for Corning® LSE Digital Dry Bath Heaters**

NAME

Corning

BRAND

480124

SKU

- 43 Place tubes (#T2) in the heater.

Note

Organic layer turns foggy when temperature is lower than  37 °C

- 44 Use a new pasteur pipet to transfer the lower organic phase to a clear 12 mL storage vial. Do not disturb the cell debris in between the two phases.

Note

Clear vial helps to check if there is water drops or impurities in lipids extract after dried.

Equipment

Glass Vials PTFE/SILiCone SEPTA Clear	NAME
16 mL	TYPE
Thermo Scientific	BRAND
B7990-4	SKU

Equipment

Screw Vial Convenience Kit, 12mL solid top PTFE cap	NAME
Thermo Scientific	BRAND
B7800-12A	SKU

45 Dry organic phase extract at  37 °C under a stream of N₂ gas (<2 psi) for about  00:30:00 .

30m

	A	B	C
		Time	Gas cylinder pressure
Start			
End			



Equipment

Reacti-Vap Evaporator

NAME

Thermo Scientific

BRAND

TS-18825

SKU

Purify lipids extract

30m

- 46 The lipids extract might still have water residue (which can't be dried by nitrogen gas) or water soluble impurities.

Redissolve it with  5 mL chloroform by using glass serological pipet, transfer certain amount of chloroform dissolved extract for lipids measurement (based on the estimation, 20 to 100 ug) into a new vial. Log the actual volume transferred.

Equipment

Safetypette

NAME

Jencons

BRAND

75856-442

SKU

- 47 Dry extract at  37 °C under a stream of N₂ gas (<2 psi) (Generally 2 mL/5 min).

	A	B	C
		Time	Gas cylinder pressure



	A	B	C
Start			
End			

48 Store dried extract and excess extract (in chloroform) at  -80 °C .

Citations

FOLCH J, LEES M, SLOANE STANLEY GH. A simple method for the isolation and purification of total lipides from animal tissues

Axelsson M, Gentili F. A single-step method for rapid extraction of total lipids from green microalgae.

<https://doi.org/10.1371/journal.pone.0089643>