

Sep 09, 2025

## Rapid, in-situ fungal DNA extraction ("Toothpick Method")

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7n5n8lwz/v1>

Harte Singer<sup>1,2,3,4</sup>

<sup>1</sup>FUNDIS; <sup>2</sup>North American Mycological Association; <sup>3</sup>Dikarya LLC; <sup>4</sup>California State University, East Bay



Harte Singer

FUNDIS, North American Mycological Association, Dikarya LLC,...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.q26g7n5n8lwz/v1>

**Protocol Citation:** Harte Singer 2025. Rapid, in-situ fungal DNA extraction ("Toothpick Method"). protocols.io <https://dx.doi.org/10.17504/protocols.io.q26g7n5n8lwz/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** September 08, 2025



**Last Modified:** September 09, 2025

**Protocol Integer ID:** 226755

**Keywords:** fresh fungal tissue for dna barcoding, situ fungal dna extraction, fungal dna extraction, sampling fresh fungal tissue, dna barcoding, dna, tissue for extended storage

## Abstract

An explanation of a simple, affordable, and versatile protocol for sampling fresh fungal tissue for DNA barcoding. Can be use to stabilize tissue for extended storage without refrigeration.

## Materials

Individually wrapped wooden toothpicks (cellophane-wrapped is ideal)

### XNA Extraction and Dilution buffers

Individual 0.2 mL swing-top PCR tubes [Amazon](#) \$55

PCR Tube Racks [Amazon](#)

Single-channel 20-200 uL micropipette [Amazon](#) \$25 OR

8-channel 50-300 uL and 5-50 uL micropipettes [Amazon](#) and [Amazon](#) \$250 for both

100 uL filtered pipette tips [Amazon](#) \$40 (As long as they are universal and not for the LTS system, they should work)

Reagent reservoirs [Amazon](#) \$50



## Preparing PCR tubes for extraction

10m 3s

- 1 Pour an appropriate volume of XNA ([dx.doi.org/10.17504/protocols.io.14egnyq8pv5d/v1](https://dx.doi.org/10.17504/protocols.io.14egnyq8pv5d/v1)) into a brand-new reagent reservoir. Each tube will need  30  $\mu\text{L}$  so use that to determine the amount you need. Add a little extra to account for pipetting error and variability.
- 2 Using a 5-50  $\mu\text{L}$  8-channel micropipette, add  30  $\mu\text{L}$  of XNA buffer into 0.2  $\mu\text{L}$  individual attached-cap PCR tubes. You can also use a p200 single-channel micropipette, however it is much slower. Close the caps after adding buffer.
- 3 Label or mark the PCR tubes before DNA extraction so that you can easily reference which sample is in which tube. For a guide on how to sample tissue and a suggestion on how to label, see this protocol which is designed around using dried tissue [dx.doi.org/10.17504/protocols.io.eq2lyw93rvx9/v2](https://dx.doi.org/10.17504/protocols.io.eq2lyw93rvx9/v2)
- 4 Open the cap of a tube of buffer, insert the toothpick tip that now has fungal tissue on it and swirl around the buffer for  00:00:03 . Remove and discard the toothpick, close the cap of the PCR tube, and mark the location of the sample in your notes.
- 5 Remove the wrapper of a fresh toothpick by pushing the toothpick through the wrapper. Avoid touching the exposed toothpick or otherwise introducing contamination.

Either dip, or scratch the fertile surface of the mushroom with the toothpick. The objective here is to get some tissue stuck to the rough, wood surface of the toothpick, but to avoid getting a chunk of tissue. You should see some schmutz on the tip, but not an actual chunk of material. This is important! Too much tissue can inhibit PCR. See the link in the previous step for suggestions on how to sample tissue.
- 6 As long as the tissue is covered with buffer, these samples are essentially shelf-stable. I have stored samples like this for at least 8 weeks at room temperature without any noticeable reduction in PCR activity. It is recommended that these be stored at  -20  $^{\circ}\text{C}$  until ready to proceed to the next step.

3s

II



- 7 Once all of the tubes are filled, spin down in a microcentrifuge for  00:00:10 then place in a PCR thermal cycler and incubate at  95 °C for  00:10:00 followed by a ramp-down to  12 °C .

10m 10s 
- 8 Spin tubes down in a microcentrifuge for  00:00:10 to bring all of the liquid and sample to the bottom of the tube.

10s
- 9 Using a brand new reagent reservoir and a 50-300 uL 8-channel micropipette (or p200 single-channel) add  80 µL of Low TE (Recommended) or molecular-grade water to each tube. Be especially careful to open the tubes without splashing any liquid. Hover the pipette tips over the open tubes and use one hand to brace the pipette in a proper position. Using a smooth and precise motion, expel the buffer into the open PCR tubes. The goal is to have the buffer come out of the pipette tips without splashing or dribbling too much. If some buffer doesn't make it into the PCR tube, it is fine. This dilution step is optional, but recommended.

 