

Jan 11, 2026

Rapid Isolation of Lysosomes (LysolP) from mouse brain tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzj414lx1/v1>

Ali Ghoochani^{1,2,3,4}, Monther Abu-Remaileh^{1,2,3,4}

¹Department of Chemical Engineering, Stanford University, Stanford, CA 94305, USA;

²Department of Genetics, Stanford University, Stanford, CA 94305, USA;

³The Institute for Chemistry, Engineering and Medicine for Human Health (Sarafan ChEM-H), Stanford University, Stanford, CA 94305, USA;

⁴Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA

asap



Ali Ghoochani

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzj414lx1/v1>

Protocol Citation: Ali Ghoochani, Monther Abu-Remaileh 2026. Rapid Isolation of Lysosomes (LysolP) from mouse brain tissue. protocols.io <https://dx.doi.org/10.17504/protocols.io.rm7vzj414lx1/v1>



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 13, 2024

Last Modified: January 12, 2026

Protocol Integer ID: 109774

Keywords: LysolP, lysosomes isolation, mouse brain lysosomes, lysosomes from fresh mouse brain tissue, lysoip from mouse brain tissue, rapid isolation of lysosome, lysosome, assessing lysosomal protein profile, lysosomal protein profile, preparing lysoip, mouse brain tissue this protocol, proteomic analysis, fresh mouse brain tissue, lysoip, purification method, mouse brain tissue, using immunoprecipitation, purification, recycling of biomolecule, protein

Abstract

This protocol outlines a method for the rapid isolation of lysosomes from fresh mouse brain tissue. Lysosomes serve as metabolic hubs for the degradation and recycling of biomolecules. Despite their critical cellular functions, quantitatively assessing lysosomal protein profiles has been challenging. To address this, we developed a rapid harvesting and purification method using immunoprecipitation (LysolP). This protocol provides detailed instructions for preparing LysolP from mouse brain tissue for proteomic analysis.

Materials

1x extraction buffer :

- 50 mM HEPES, pH 7.4
- 40 mM NaCl
- 2 mM EDTA
- 1% Triton X-100 (Sigma)
- 1.5 mM sodium orthovanadate (NaVO₄)
- 30 mM sodium fluoride (NaF)
- 10 mM sodium pyrophosphate (Na₄P₂O₇)
- 10 mM sodium β-glycerophosphate
- Complete EDTA-free Protease Inhibitor Cocktail

Materials

1 Materials and Reagents

1. **Anti-HA magnetic beads:** Thermo Fisher Scientific
2. **Optima LC/MS water**
3. **Phosphate-buffered saline (PBS):** pH ~7.4
4. **Complete EDTA-free Protease Inhibitor Cocktail:** Roche
5. **Magnetic separation equipment:** DynaMag Spin Magnet
6. **1x extraction buffer (Material section)**
7. **Centrifuge**
8. **Homogenizer:** 2 ml homogenizer

Equipment Preparation

- 2 Wash the glass vessel homogenizer with MilliQ Water, 10 times each. wash the tissue grinder homogenizer thoroughly with DI Water and MilliQ Water, especially the gap between the white parts, don't touch the part that goes into the glass vessel. Then dry upside-down using paper towels. Carefully place the glass vessels against something to prevent falling down. Minimize any contact between the grinder and anything else.
 - **Pre-chill all equipment and buffers** used for homogenization and immunoprecipitation to 4°C.
 - **Set centrifuge temperature to 4°C** for all centrifugation steps.
- 3 **Prepare microcentrifuge tubes as follows on a metal rack on ice (for each sample, from left to right):**
 - ① 2 mL tube for posthomogenization
 - ② 1.5 mL tube for whole brain fraction
 - ③ 1.5 mL tube for beads
 - ④ 1.5 mL tube for 3rd wash
 - ⑤ 1.5 mL tube for final lysolP samples

Preparation of Anti-HA beads

- 4 Pool all required beads volumes together ( 100 μL / sample, e.g.  800 μL total for 8 samples, extra is not needed).
- 4.1 Wash 3 times with the same volume ( 100 μL / sample, e.g.  800 μL total for 8 samples) of ice-cold PBS containing Complete EDTA-free Protease Inhibitor Cocktail. After each wash, place the tube on the magnet, aspirate the supernatant, and add the same volume of ice-cold PBS containing Complete EDTA-free Protease Inhibitor Cocktail.
- 4.2 After this wash, Resuspend beads with the same volume ( 100 μL / sample, e.g.  800 μL total for 8 samples) of ice-cold PBS containing Complete EDTA-free Protease Inhibitor Cocktail.
- 4.3 Aliquot 100 μL into each ③ 1.5 mL tube for beads.

Isolation of Lysosomes from Mouse Brain

- 5 Tissue Collection and Preparation:
 - 5.1 Process each mouse brain separately to ensure rapid isolation of lysosomes.
 - 5.2 Euthanize the mouse and transfer the whole brain into 2 ml homogenizer (on ice)

Homogenization


- 6 **Critical: All steps of the procedure must be performed on ice** 
- 6.1 Immediately add  950 μL ice-cold PBS containing Complete EDTA-free Protease Inhibitor Cocktail into 2 ml homogenizer.
- 6.2 Homogenize with 25 gentle strokes (on ice).
- 6.3 Transfer the homogenized tissue to a 2 ml tube.

6.4 Reserve  25 μL of the homogenate (2.5% of the total tissue) for processing as the whole brain fraction.

6.5 **Recommended:** Add 1 ml of ice-cold PBS with Protease Inhibitor Cocktail to the 2ml tube (Step 6.3) to minimize pellet contamination after centrifuge (next step).

6.6 Centrifuge at  1000 x g, 00:02:00 at 4°C.

2m



Immunoprecipitation of Lysosomes

15m 25s

7 Transfer the supernatant containing cellular organelles to a 1.5 ml tube containing  100 μL of anti-HA magnetic beads

7.1 Incubate on a gentle rotator shaker for  00:15:00 at  4 °C

15m

7.2 Place the tube on magnet. Count at least  00:00:25 to allow for beads to be pulled by magnets.

25s

Note: it is important to keep this count the same between each wash and each sample



for consistency i.e. 25 seconds each time .

- 7.3 Wash the beads fraction three times with  1 mL PBS containing Protease Inhibitor Cocktail using a DynaMag Spin Magnet.

Note:

- During the first wash, make sure to aspirate any liquid trapped on the inner side of the cap.
- Pipet up and down 3 times and keep consistent each wash, each sample.

- 7.4 After the second wash, resuspend beads with PBS containing Protease Inhibitor Cocktail and then transfer it to the clean ④ 1.5 ml tube for third wash (this step helps give cleaner results).

- 7.5 Place the tube on magnet wait  00:00:25 to allow beads to be pulled by magnets. Then, aspirate the supernatant.

25s

Lysosomal Protein Extraction

30m

- 8 Resuspend the beads with  80 µL of ice-cold **1x extraction buffer**.

- 8.1 Incubate at least for  00:15:00 at  4 °C

15m

- 8.2 After  00:15:00 finishing the last IP, place the IP samples on the magnet to separate the beads. **Collect the supernatant and transfer it** to a new 1.5 ml tube. Store the IP samples at  -80 °C until further processing.

15m

- 8.3 For whole brain fraction add  160 µL of ice-cold 1x extraction buffer.

25m

- Incubate at least for  00:15:00 at  4 °C
- Centrifuge  15000 rpm, 00:10:00  4 °C
- Transfer the supernatant to a new tube. Store the samples at  -80 °C until further processing.

Quality Control



- 9
 - Add loading dye to the sample before running it on SDS-PAGE. Load 1% of the whole brain lysate and 10% of LysolP.
 - Detect lysosomal markers (LAMP1, Cathepsin B), Golgi marker (Golgin-97), endoplasmic reticulum marker (Calreticulin), and mitochondrial marker (VDAC) to confirm the enrichment of lysosomes.