



Dec 29, 2025

Version 2

Rapid On-Site Detection of Crop RNA Viruses Using CRISPR/Cas13a V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk9wb5v5r/v2>

Hagit Hak¹, Steffen Ostendorp², Anton Reza², Shany Ishgur Greenberg³, Gur Pines⁴, Julia Kehr², Ziv Spiegelman¹

¹Department of Plant Pathology and Weed Research, Agricultural Research Organization - Volcani Institute, Rishon LeZion, Israel;

²Department of Biology, Institute of Plant Science and Microbiology, Universität Hamburg, Hamburg, Germany;

³Extension Service, Ministry of Agriculture and Rural Development, Beit Dagan 50250, Israel;

⁴Department of Entomology, Agricultural Research Organization - Volcani Institute, Rishon LeZion, Israel

**Hagit Hak**

ARO, Volcani Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)OPEN  ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk9wb5v5r/v2>

Protocol Citation: Hagit Hak, Steffen Ostendorp, Anton Reza, Shany Ishgur Greenberg, Gur Pines, Julia Kehr, Ziv Spiegelman 2025. Rapid On-Site Detection of Crop RNA Viruses Using CRISPR/Cas13a. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.j8nlk9wb5v5r/v2> Version created by [Hagit Hak](#)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 29, 2025

Last Modified: December 29, 2025

Protocol Integer ID: 235993

Keywords: CGMMV, Diagnostics, On-site, Tomato brown rugose fruit virus, ToBRFV, TuMV, site detection of crop rna virus, site detection of plant virus, crop rna virus, cas13a crispr, plant virus, using crispr, crispr, nucleic acid, virus, nucleic acids sequence, site detection, rna, rnasealert, concentration of rnasealert

Funders Acknowledgements:

Israeli Ministry of Agriculture and Rural Development

Grant ID: 20-02-0187

National Center for Genome Editing Applications and Technologies in Agriculture

Grant ID: 20-01-0209

European Research Council (ERC) under the European Union's Horizon 2020 research and innovation programme

Grant ID: 810131

Deutsche Forschungsgemeinschaft

Grant ID: DFG; GA No. 433194101, Research Unit 5116

Abstract

CRISPR/Cas technology is an emerging tool for the identification of nucleic acids sequences. Here, we present a user-friendly, extraction-free, rapid protocol for specific on-site detection of plant viruses using CRISPR/Cas13a.

Version 2: corrected the concentration of RNaseAlert in step 9 (from 200μM to 200nM).

Materials

Before starting:

- crRNAs – synthetic RNA oligos
- LbuCas13a
- RNaseAlert® substrate (IDT, Coralville, USA)
- Lysis buffer (15% polyethylene glycol 4000, 20 mM NaOH)
- **Reaction buffer:**

	A	B
	407 Hepes pH 6.8	20 mM
	KCl	50 mM
	MgCl ₂	5 mM
	Glycerol	5%

- A healthy plant to serve as control.

Equipment:

- P51™ Molecular Fluorescence Viewer from MiniPCR Bio (Cambridge, MA, USA) Mobile phone with camera

Nuclease assay:

Mix the following

	A	B
	Reaction buffer	5 µL
	RNPs	1 µL (100 nM final volume)
	RNaseAlert® substrate (IDT, Coralville, USA)	1 µL (200 µM final volume)
	Crude RNA	3 µL



crRNAs

- 1 Design 23-long crRNAs that complement the RNA sequence of the virus of interest.
- 2 Order the crRNAs fused to the 3' end of the LbuCas13a stem (5'GACCACCCCAAAAAUGAAGGGGACUAAAAC 3').
- 3 Resuspend crRNA oligos with molecular grade water to [M] 1 micromolar (μM) .

RNP Assembly

30m

- 4 Mix crRNA with LbuCas13a at a 1:1 ratio. Incubate at [T] Room temperature for [D] 00:30:00 .
 
- 5 RNPs can be stored at [T] -20 °C for up to at least 6 months.

On-site: Crude RNA Extraction

5m 5s

- 6 Punch a 1 cm leaf disc and place in [V] 300 μL of lysis buffer.
- 7 Incubate at [T] Room temperature for [D] 00:05:00 .

- 8 Shake vigorously for [D] 00:00:05 .


Note

Use immediately or place [T] On ice .

On site: Nuclease Assay

10m



9 Mix the following in a PCR tube:

	A	B
	Reaction buffer	5 μ L
	RNPs	1 μ L (100 nM final volume)
	RNaseAlert substrate (IDT, Coralville, USA)	1 μ L (200 nM final volume)
	Crude RNA	3 μ L

10 Incubate at  Room temperature for  00:10:00 .

10m



11 View fluorescence using P51™ Molecular Fluorescence Viewer. Since exposure to light interferes with fluorescence, it's best to place the Fluorescence Viewer in a dark box (we use a small shoebox).

12 Take a photo with any mobile phone camera. It is better to set the brightness to a medium level.

13 Compare fluorescence levels of tested plants with that of the control.

14 **Optional:** Fluorescence quantification using ImageJ (<https://imagej.net/ij/>).

Select an area of the liquid in the tube. Be careful not to include reflections of the LED lights. Measure grayscale by clicking Analyze → Measure. Compare Mean value with that of a control sample.

