

May 08, 2023

Rapid Ribosome (Polysome) Profiling

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7ykmqgwz/v1>



Jason D Limberis¹, Yuxiang Chen²

¹University of California, San Francisco; ²Stanford University, Stanford



Jason D Limberis

University of California, San Francisco

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7ykmqgwz/v1>

Protocol Citation: Jason D Limberis, Yuxiang Chen 2023. Rapid Ribosome (Polysome) Profiling. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7ykmqgwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

While we have not tested this protocol, it is theoretically sound.

Created: April 10, 2023

Last Modified: May 08, 2023

Protocol Integer ID: 80279

Keywords: profiling ribosome profiling, traditional ribosome profiling method, rapid ribosome, ribosome, positions of ribosome, mrna fragment, protected mrna fragment, translation rate, translated open reading frame, genome, transcript, translation, sequencing, open reading frame

Abstract

Ribosome profiling is a powerful technique used to study translation at a genome-wide level. It involves the sequencing of ribosome-protected mRNA fragments to determine the positions of ribosomes on transcripts. This information can be used to infer translation rates and identify translated open reading frames. While traditional ribosome profiling methods can be time-consuming and expensive, our method is rapid and cost-effective, and the inclusion of UMIs allows for precise quantification.

Attachments



[riboSeek.png](#)

244KB

Materials

A	B
Name	Sequence
RiboS_linker	rCAAGCAGAAGACGGCATACGAGAT
RiboS_linker_primer	ATCTCGTATGCCGTCTTCTGCTTG
RiboS_IlluminaAdapt_TSO_UMI_RNA_F	ACACTCTTTCCCTACACGACGCTCTTCCGATCTNNNNNN GATrGrGrG
RiboS_blocking_oligo	CTACCCCAAGCAG
Index 1 Read	GATCGGAAGAGCACACGTCTGAACTCCAGTCAC[i7]ATCTC GTATGCCGTCTTCTGCTTG
Index 2 Read	AATGATACGGCGACCACCGAGATCTACAC[i5]ACACTCTTT CCCTACACGACGCTCTTCCGATCT

Oligonucleotides. "r" indicated RNA base.

☒ Monarch RNA Cleanup Kit (10 µg) - 100 preps **New England Biolabs Catalog #T2030L**

☒ Agencourt RNAClean XP Magnetic Beads **Beckman Coulter Catalog #A63987**

☒ RtcB Ligase - 25 rxns **New England Biolabs Catalog #M0458S**

☒ Superase-In RNase Inhibitor **ThermoFisher Catalog #AM2694**

☒ Template Switching RT Enzyme Mix – 100 rxns **New England Biolabs Catalog #M0466L**

Optional

Urea page

TBE

Agarose

Before start

Read through the protocol and ensure you have the correct reagents depending on your chosen approach



RNA isolation and nuclease footprinting

- 1 Perform RNA isolation and nuclease footprinting as appropriate for your cell type (see steps 1-17 [here](#) for cultured mammalian cells)
- 2 Optional but recommended preselection of small RNA fragments using either a urea PAGE gel or bead and membrane selection.

STEP CASE

15% UREA-TEB PAGE

45 steps

Recipe for casting own gel:

Urea 19.2g

40% Acrylamide/Bis (19:1) 15ml

10X TBE 4ml

25% APS 200ul

TEMED 20ul

Adding nuclease free water to 40ml

- 3 Add the appropriate volume of 2X RNA loading buffer to each RNA sample
- 4 Denature the samples for 80 °C 00:01:00 followed by ≥ 00:02:00 on On ice 3m
- 5 Load the samples on the polyacrylamide gel with TBE running buffer
- 6 Pre-run the gel in 1X TBE buffer at 100V for 10min

Optional: preselection of small RNA fragments

5m

- 7 Separate by electrophoresis for 01:30:00 at 100V 1h 30m
- 8 Stain the gel for 00:03:00 with 1X SYBR Gold in 1X TBE running buffer on a shaker 3m
- 9 Visualize the gel and excise desired region (20-40nts)



- 10 Transfer the excised gel slice to a  1.5 mL microcentrifuge tube and weigh it.
- 10.1 Pro Tip: Smashing the gel in the tube increases RNA recovery.
- 11 Add 4 volumes of Monarch RNA Cleanup Binding Buffer to the tube with the slice
- 12 Incubate the sample between  55 °C , gently mixing periodically until the gel slice is completely dissolved (generally <  00:10:00) 10m
- 13 Add two volumes of absolute ethanol to the sample and mix well by pipetting up and down
- 14 Incubate the sample between  55 °C or an additional  00:05:00 5m
- 15 Insert an RNA cleanup column into a collection tube, load the sample onto the column and close the cap
- 15.1 Spin for  00:01:00 , then discard flow-through 1m
- 16 Re-insert the column into the collection tube
- 16.1 Add  500 µL RNA Cleanup Wash Buffer
- 16.2 Spin for  00:01:00 , then discard flow-through 1m
- 16.3  15 for a total of two washes.
- 17 Transfer the column to a clean 1.5 ml microfuge tube



18 Add  10 μ L nucleae-free water directly to the membrane and incubate at  Room temperature for  00:05:00


5m

19 Spin for  00:01:00

1m

Overview

20



Linker ligation

21

Add 1 μL of 1M 100 micromolar (μM) RiboS_linker (linker RNA oligo) to 10 μL of RNA, denature it for 00:01:00 at 70 $^{\circ}\text{C}$, and then cool it to

1m



Room temperature

22 Set up the ligation reaction below and incubate for 01:30:00 at 37 °C 1h 30m

	A	B
Component	Volume (μl)	
RNA and linker	11	
RtcB Reaction Buffer (10X)	2	
SUPERase-In (20 U/μl)	1	
1 mM GTP	2	
10 mM MnCl ₂	2	
RtcB RNA Ligase	1.5	
H ₂ O	0	

RNA cleanup

23 Add 100 μL RNA Cleanup Binding Buffer

24 Add 150 μL absolute ethanol and mix by pipetting

25 Insert an RNA cleanup column into a collection tube, load the sample onto the column and close the cap

25.1 Spin for 00:01:00 , then discard flow-through

26 Re-insert the column into the collection tube

26.1 Add 500 μL RNA Cleanup Wash Buffer

26.2 Spin for 00:01:00 , then discard flow-through

26.3  6 for a total of two washes.

27 Transfer the column to a clean 1.5 ml microfuge tube

28 Add  10 μL nucleae-free water directly to the membrane and incubate at  Room temperature for  00:05:00

29 Spin for  00:01:00

RT (TSO) reaction

2h 37m 30s

30

6m

A	B
Component	Volume (μl)
RNA	9
RiboS_linker_primer (100uM; reverse transcription primer)	2
dNTP mix (10 mM each)	2

Mix thoroughly by gently pipetting up and down at least 10 times, then centrifuge briefly to collect the solution to the bottom of the tube

Denature for  00:05:00 at  70 $^{\circ}\text{C}$ in a thermal cycler and then place on ice for $\geq$  00:01:00

31 Vortex the Template Switching RT Buffer briefly followed by a quick spin
Combine the following components in a reaction tube

A	B
Component	Volume (μl)
Template Switching RT Buffer	5
RiboS_IlluminaAdapt_TSO_UMI_RNA_F (75uM)	1
Template Switching RT Enzyme Mix	1



RT reaction mix

Mix thoroughly by gently pipetting up and down at least 10 times, then centrifuge briefly to collect the solution to the bottom of the tube

- 32 Combine  7 μL RT reaction mix (above) with  13 μL of the annealed mix, mix well by gently pipetting up and down at least 10 times, then centrifuge briefly to collect the solution to the bottom of the tube

Incubate as below

	A	B
	Temp (C)	Time (m)
	42	90
	85	5
	4	1

Index PCR

2h 37m 30s

33

	A	B
	Component	Volume (μL)
	10 μM Forward Index Primer	2.5
	10 μM Reverse Index Primer	2.5
	RiboS_blocking_oligo	2
	2X Phusion Master Mix	25
	TSO reaction	18

PCR components

	A	B	C
	Cycles	Temp (C)	Time (s)
	1	98	30
	10	98	5



	A	B	C
		60	5
1		72	10

PCR cycle parameters (for selection of <40bp inserts)

Fragment is 150bp excluding insert at this stage

AATGATACGGCGACCACCGAGATCTACAC[i5]ACACTCTTTCCCTACACGACGCTCTTCCGATCTNNNNNNG

ATGGG[INSERT]CAAGCAGAAGACGGCATACGAGAT[i7]GTGACTGGAGTTCAGACGTGTGCTCTTCCGATC

34 Add  90 μ L AMPure XP beads to the PCR product

8m

Incubate at  Room temperature for  00:05:00

Place on a magnetic rack

Aspirate supernatant

Add  200 μ L [M] 70 % (v/v) ethanol

Wait for  00:00:30

Aspirate and discard the supernatant

Add  200 μ L [M] 70 % (v/v) ethanol

Wait for  00:00:30

Aspirate and discard the supernatant

Resuspend beads in  10 μ L of H₂O

Incubate for  00:02:00

Transfer to a clean PCR tube

Optional: Agarose gel size selection

35 Run a 1% agarose gel, excise the desired band, and perform a gel cleanup.

Optional: Depletion of Abundant Sequences by Hybridization

36 To perform DASH to deplete unwanted sequences (e.g., rRNA) follow the protocol [here](#).

Quantification

37 Run a tapestation or bioanalyzer chip



Do qubit or qPCR using an Illumina library quant kit to quantify the library and pool the samples (see [Illumina Sequencing Coverage Calculator](#) for pooling information)

Analysis

- 38 Analysis can be done using common pipelines with the addition of [umi_tools](#) to count unique reads

Protocol references

Ingolia, N., Brar, G., Rouskin, S.*et al.* The ribosome profiling strategy for monitoring translation *in vivo* by deep sequencing of ribosome-protected mRNA fragments. *Nat Protoc* **7**, 1534–1550 (2012).

<https://doi.org/10.1038/nprot.2012.086>