

Jun 13, 2019

RbCl Uber-Competent Cells

DOI

dx.doi.org/10.17504/protocols.io.3hygj7w

Sricharan Kadimi¹

¹University of Connecticut

UConn iGEM



Sricharan Kadimi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3hygj7w>

Protocol Citation: Sricharan Kadimi 2019. RbCl Uber-Competent Cells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3hygj7w>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our group and it is working, though open to further improvement/modification

Created: May 30, 2019

Last Modified: June 13, 2019

Protocol Integer ID: 23832

Keywords: Competent Cells, competent cell, good competent cell, good competent cells in summer, rbcl uber, cell, protocol, similar to the hanahan, hanahan, mit igem team

Abstract

A protocol kindly provided to our team by the advisor to the 2018 MIT iGEM team (Dr. Brian Teague). It is modified slightly from that version (made by Stefan Maas in 1997) to fit our needs. Stefan Maas' protocol is similar to the Hanahan, D. protocol from 1983.

It worked pretty well to make good competent cells in Summer 2018 and we just used it again in Summer 2019.

Materials

MATERIALS

⊗ MOPS P212121

⊗ Manganese(II) chloride tetrahydrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #M3634

⊗ Glycerol Bio Basic Inc. Catalog #GB0232.SIZE.500ml

⊗ Rubidium chloride Bio Basic Inc. Catalog #RB0668.SIZE.25g

⊗ Potassium acetate Merck Millipore (EMD Millipore) Catalog #1.04820.1000

⊗ Calcium Chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #C4904

⊗ Potassium hydroxide Merck MilliporeSigma (Sigma-Aldrich) Catalog #1050121000

⊗ Acetic acid Merck MilliporeSigma (Sigma-Aldrich) Catalog #695092

Make sure you also have all the glassware/plasticware you need (listed in the first step), as well as:

- SOB media + 20mM MgCl₂
- LB plates without antibiotic (2)
- LB media
- Plastic tube racks with foam inside
- Labeled boxes for storing 1.5mL tubes

Troubleshooting



Day 0

1 Autoclave the following materials:

2h

	Amount	Material	Comments
	2	2L Erlenmeyer flasks	For culturing cells, second one is if you mess up the first time
	2	500 mL bottles	For TfbI
	2	100 mL bottles	For TfbII
	2	Stir bars	For stirring buffers
	~200	1.5mL tubes in a glass beaker	For aliquoting competent cells (the whole beaker should just be covered)



			in foil)
	1	500 mL graduated cylinder	For measuring media
	1	250 mL graduated cylinder	For QSing TfbI
	1	100 mL graduated cylinder	For QSing TfbII
	4	200 mL centrifuge bottles	For centrifuge steps, have different orange caps from other bottles


Make sure all bottles have their caps on loosely!

For best results, autoclave glassware 3/4-full with DI (to remove residual detergent that could interfere with bacterial growth or competency). Plasticware should be rinsed with DI 2-3x and then autoclaved on dry cycle. Any containers without caps should be covered securely in aluminum foil. Autoclave tape should be applied if necessary to foil or outside container walls. See the Autoclave protocol for more details.

- 2 Take a tube of commercially competent cells and put it on ice (Thermo-Fisher TOP10 preferred, we think). Make sure to use a tube of commercially competent cells to start, or you may have poorer results. Work near a flame to ensure sterility.

20m



Using a sterile pipet tip or a flamed and cooled inoculation loop, take a drop of the thawed competent cells and streak gently onto an LB plate. Put the plate upside-down (agar-side up) in the  37 °C incubator for 16-18 hours. Discard the tube of competent cells in bacterial waste (or use for a transformation if possible).



See the Streaking Cells protocol for more details.

- 2.1 If you work fast, you might be able to pick a chip of ice from the tube quickly, then return the tube to the -80°C freezer. However, we don't recommend this unless you're experienced, since once the competent cells start warming over -80°C, they lose competency.

Day 1

- 3 Make sure you take the plates from the incubator in the morning and leave them in the cold room until you are ready to do your overnight culture.

20m



Work near a flame to ensure sterility. Using a serological pipet, fill 2-3 14mL round-bottomed culture tubes with  5 mL of LB media. For each culture tube, pick a colony from your plate with a p200 tip and drop it into the media. Put the tubes on a rack and put the rack in the  37 °C shaker at 250rpm for 16-18 hours.

(You should be able to safely use your SOB + 20mM MgCl₂ for this step, that might actually be better. However, we've only tried this with LB).

Do NOT add any antibiotic to the culture media - the cells are not antibiotic-resistant and should not be either!

See the Overnight Culture protocol for more details.

Day 2


- 4 Work near a flame to ensure sterility. Add  800 µL fresh overnight culture to  400 mL medium in 2L flask.

5m

(We used LB but 2xYT and NZDT work, SOB + 20mM MgCl₂ is optimal).


Flask should be 5x the volume of the medium and covered with foil.



- 5 Shake the culture at 250 rpm at  37 °C for 3-4 hours until at A_{600} of 0.5-0.7 (for cuvettes of length 1 cm, $A_{600}=OD_{600}$).

Use the spectrophotometer and cuvettes at our bench (in Quartzzy hopefully) for checking the OD (it has a cuvette length of 1cm). Make sure to blank with media. Swirl the culture before taking a 1mL aliquot. Do not return the aliquot to the flask. Work near a flame for sterility.


Note



This is the step where you are most likely to mess up! Start checking your culture after about  02:30:00 and then start checking every 10-15 minutes after that (depending on how much the culture grew). Note that while the culture initially takes a while to get to 0.1 or 0.2 OD_{600} (lag phase), the OD will start shooting up during the log (exponential phase).

- 6 While the cells are shaking, load foam racks and plastic tube racks with open tubes in the cold room. (Put some paper towels over them to prevent dust or anything else from getting into them).

15m

Also leave the centrifuge bottles, a box of p200 pipette tips, and several 10mL and 25mL serological pipettes in the cold room to chill.

Also turn on the large centrifuge next to the 2nd floor cold room and set it to  4 °C ,

 4000 rpm ,  00:10:00 so it can pre-cool.

- 7 Also while the cells are shaking, make fresh TfbI and TfbII. Ideally, these should be made fresh, but if necessary, they can be made the day before as well.

1h 30m

	Com pon ent	TfbI (150 mL)	TfbII (20 mL)
	Rubi dium chlo ride	1.8 g	0.02 4 g
	Man gane se (II) chlo ride tetra	2.30 g	non e

	hydrate		
	Calcium chloride	0.17g	0.17g
	Potassium chloride	0.44g	none
	MOPS	none	0.042g
	Glycerol	22.4 mL (28.4g)	3.0 mL (3.78g)

Add all dry reagents to the bottle, add glycerol (usually easier just to weigh it) ~3/4 of the autoclaved MilliQ water (in the past, we just used DI, but autoclaved MilliQ is best) needed to QS.

Then adjust pH:

- For TfbI adjust to pH 5.8 with 0.2M acetic acid.

Note


You may need to let the TfbI stir a long while to get the pH up to 5.8 or higher before you pH adjust back down!

Add slowly! If you add too much then the solution will turn yellow-green and you will have to remake it.

- For TfbII adjust to pH 6.5 with 1M KOH.

Then transfer to graduated cylinder to QS to final volume with autoclaved MilliQ water.

For TfbI and TfbII, filter sterilize with a  0.4 μm filter (not sure about the pore size, but it's what we used).

8 Chill cell culture in ice-water bath for 10-15 minutes.  00:10:00

10m

**Note**




All steps from here on should be done in the cold room. Make sure to keep the cells as cold as possible for the remainder of the protocol - so all containers, buffers, tips, etc should be cold as well!

- 9 Separate culture evenly between two chilled 250 mL centrifuge bottles.

10m

Note

Make sure that the cap of the centrifuge bottle screws on correctly and tightly! There are some dud caps that could cause your cell culture to spill out into the centrifuge chamber if used. Also check if there might be any autoclaved O-rings that go with the centrifuge bottles, because that would also help prevent spills.

Centrifuge at  4000 rpm (about 3000x g) at  4 °C for  00:10:00 .

- 10 In the cold room, decant medium and drain on a paper towel to allow all traces of medium to drain out.

5m


Tap the bottles on paper towels but be careful not to lose any cells (the pellet may start to slide down the side).

- 11 Gently resuspend each cell pellet in 4-8 mL ice-cold Tfb1 (use serological pipette).

2h 10m

"Gently" means lightly swirling the flask until all cells are resuspended. Don't shake, poke, or pipet.




DO NOT disturb the cell pellet by pipetting back and forth or shaking vigorously.

This may take a while, but make sure the pellets are fully resuspended (uniform solution). Pour one of the resuspended pellets into the bottle with the other resuspended pellet and bring the total volume to  120 mL (use serological pipette).

Incubate on ice in the cold room for  02:00:00 .


**Note**

This is usually a good time to go get lunch :)


- 12 Centrifuge again at  4000 rpm (about 3000x g) at  4 °C for  00:10:00 .

10m

Make sure to label the tube containing the cells. Then, fill another centrifuge bottle with water as a balance (and label that differently).

- 13 In the cold room, decant the buffer and discard. Gently resuspend cells in  16 mL ice cold TfbII (this could take 5-10 minutes, similar to before). Aliquot cells into chilled 1.5mL tubes.


25m

Aliquots should be 50-200µL (usually  100 µL). Chill your pipette tips in the cold room before aliquotting and swirl cells regularly during aliquotting. Try to work quickly!

- 14 Close the caps of all the 1.5mL tubes in the foam rack. Fill a plastic ice bucket with liquid nitrogen (ask for a PI or grad student's help here). Then place the foam rack into the liquid nitrogen for snap freezing. You will know the tubes are frozen when most of the popping/crackling sounds subside.

15m

Snap freezing ideally should be done in liquid nitrogen (as detailed above). However, you can also use dry ice if necessary.

Make sure you have some (pre-labeled!) empty freezer boxes to transfer the 1.5mL tubes to after snap freezing. Then put the boxes in the  -80 °C freezer for storage.

If necessary, transfer more of the 1.5mL tubes to the emptied foam rack and repeat until all cell aliquots are snap frozen.

Testing Competency

- 15 Make sure to test your cells for transformation efficiency using the iGEM Comptent Cell Test Kit or similar. (BBa_J04450).

Refer to the RbCl Cells Transformation protocol for more details.