

May 08, 2017

Real-time PCR

DOI

<https://dx.doi.org/10.17504/protocols.io.huzb6x6>



Ruiyun Bi

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.huzb6x6>

Protocol Citation: Ruiyun Bi: Real-time PCR. protocols.io <https://dx.doi.org/10.17504/protocols.io.huzb6x6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 08, 2017

Last Modified: November 09, 2017

Protocol Integer ID: 5753

Keywords: real time pcr, time pcr

Total RNA isolation

- 1 The bilateral TGs from four rats per group were dissected and pooled for total RNA isolation using TRIzol reagent (Invitrogen, Carlsbad, CA, USA).

Reverse transcription

- 2 Reverse transcription PCR was conducted with an iScript cDNA synthesis kit (Bio-Rad) in 20 µl reaction volume containing 1 µg of total RNA incubated at 25°C for 5 min, transcribed at 42°C for 30 min, and terminated by heating at 85°C for 5 min. The synthesized cDNA was stored at -20°C until use.

Real-time PCR

- 3 Real-time PCR was performed with Power SYBR Green PCR Master Mix (Applied Biosystems) using a 7500 real-time PCR System (Applied Biosystems). The reactions were run in duplicate with 1 µl of cDNA template in a 20 µl reaction volume with the program running at 50°C for 2 min and 95°C for 10 min, followed by 40 cycles of 94°C for 15 s and 60°C for 1 min. The amplification specificity was confirmed by melting curve. The mRNA level of the target gene was acquired from the value of threshold cycle (Ct) as a relative level to that of β -actin through the formula $2^{-\Delta Ct}$ ($\Delta Ct = \beta\text{-actin Ct} - \text{gene of interest Ct}$). The efficiency of the primers was confirmed by sequencing the conventional PCR products before applying for real-time PCR.