


Jan 04, 2023

Recombinant expression and purification of HIV-1 RT

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4jzyjlo5/v1>



Javiera A Avilés¹, Tamara Matute², Isaac Núñez², Maira Rivera², Javiera Reyes², Jenny Molloy³, Cesar A Ramirez-Sarmiento², Fernan Federici¹

¹Millennium Institute for Integrative Biology (iBio), Santiago, Chile;

²Institute for Biological and Medical Engineering, Pontificia Universidad Católica de Chile;

³Open Bioeconomy Lab, University of Cambridge

Reclone.org (The Reagent Collaboration Network)

Tech. support email: protocols@recode.org

[Click here to message tech. support](#)



Cesar A Ramirez-Sarmiento

Pontificia Universidad Catolica de Chile

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4jzyjlo5/v1>



Protocol Citation: Javiera A Avilés, Tamara Matute, Isaac Núñez, Maira Rivera, Javiera Reyes, Jenny Molloy, Cesar A Ramirez-Sarmiento, Fernan Federici 2023. Recombinant expression and purification of HIV-1 RT. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4jzyjlo5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 04, 2023

Last Modified: January 04, 2023

Protocol Integer ID: 74741

Keywords: recombinant expression, dialysis, purification, sequence of the plasmid, plasmid, fast buffer exchange

Funders Acknowledgements:

ANID Millennium Science Initiative Program

Grant ID: ICN17_022

ANID CONCYTEC

Grant ID: covbio0012

Abstract

This protocol has been optimized for the recombinant expression of HIV-1 RT. The sequence of the plasmid encoding HIV-1 RT can be found on reclone.org.

The goal of this protocol was to eliminate the use of large volumes for dialysis and its fast buffer exchange into storage conditions.



Materials

MATERIALS

☒ Sodium phosphate monobasic monohydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9638**

☒ PMSF **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7626**

☒ Sodium phosphate dibasic **Merck MilliporeSigma (Sigma-Aldrich) Catalog #7558-79-4**

☒ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**

☒ Imidazole **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I5513**

☒ NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**

☒ HisTrap FF Crude Column **GE Healthcare Catalog #17528601**

☒ Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516**

☒ Dextrose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9434**

☒ Tween-20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416**

☒ EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #ED2SS**

☒ Trizma Base **Catalog #93362**

☒ Lysozyme **Thermo Fisher Scientific Catalog #89833**

☒ Amicon Ultra-15 Centrifugal Filter Unit **Merck Millipore (EMD Millipore) Catalog #UFC910024**

Buffer A, pH 8.0

[M] 50 millimolar (mM) NaPO₄, pH 8.0

[M] 50 millimolar (mM) dextrose

[M] 300 millimolar (mM) NaCl

[M] 1 millimolar (mM) EDTA

[M] 0.1 % volume Tween-20

[M] 40 millimolar (mM) Imidazole, pH 8.0

Buffer B, pH 8.0

[M] 50 millimolar (mM) Tris-HCl, pH 8.0

[M] 300 millimolar (mM) NaCl

[M] 1 millimolar (mM) EDTA

[M] 0.1 % volume Tween-20

[M] 10 % volume Glycerol

[M] 40 millimolar (mM) Imidazole, pH 8.0



Buffer C, pH 8.0

- [M] 50 millimolar (mM) Tris-HCl, pH 8.0
- [M] 300 millimolar (mM) NaCl
- [M] 1 millimolar (mM) EDTA
- [M] 0.1 % volume Tween-20
- [M] 10 % volume Glycerol
- [M] 150 millimolar (mM) Imidazole, pH 8.0

Buffer D, pH 8.0

- [M] 50 millimolar (mM) Tris-HCl, pH 8.0
- [M] 300 millimolar (mM) NaCl
- [M] 1 millimolar (mM) EDTA
- [M] 0.1 % volume Tween-20

Storage Conditions

- [M] 25 millimolar (mM) Tris-HCl, pH 8.0
- [M] 150 millimolar (mM) NaCl
- [M] 0.5 millimolar (mM) EDTA
- [M] 5 millimolar (mM) DTT
- [M] 0.05 % volume Tween-20
- [M] 50 % volume Glycerol



DAY 1 – Plasmid transformation

1d

- 1 Transform  100 ng of plasmid containing HIV-1 RT into *E. coli* BL21(DE3) competent cells using either heat shock or electroporation. 2h
- 2 Spread transformed cells in LB Agar plates supplemented with  0.05 mg/mL **Kan**. Grow plate overnight at  37 °C . 12h

DAY 2 – Preinoculum

1d

- 3 Select a single colony from the LB agar plate to prepare a preinoculum in  10 mL LB media supplemented with  0.05 mg/mL **Kan**. Grow overnight at  250 rpm, 37°C . 1d

DAY 3 – Protein Overexpression

1d

- 4 Use the full volume of the preinoculum to inoculate  1 L of LB media supplemented with  0.05 mg/mL **Kan** (1% inoculation). Grow at  250 rpm, 37°C until reaching an optical density at 600 nm (OD₆₀₀) = 0.8. 4h
- 5 Upon reaching OD₆₀₀ = 0.8, add  0.5 millimolar (mM) IPTG and incubate for 2h at  220 rpm, 37°C . 2h

DAY 4 – Protein Purification by IMAC

6h

- 6 Centrifuge the cell culture  4000 x g, 4°C, 00:20:00 .Then, resuspend the cell pellet in  50 mL of **Buffer A** freshly supplemented with  0.5 millimolar (mM) PMSF and  0.2 mg/mL lysozyme. 30m
- 7 Incubate the resuspended cells  80 rpm, Room temperature , 00:30:00 . 30m
- 8 Sonicate on ice for  00:04:00 using cycles of  00:00:06 ON and  00:00:06 OFF at 40% amplitude (Qsonica Q125, 125W). 30m
- 9 Centrifuge the unclarified lysate  20000 x g, 4°C, 00:20:00 and collect the supernatant. You might want to collect a small sample for SDS-PAGE afterwards. 30m



- 10 On a **1 mL HisTrap column (Ge Healthcare)** preequilibrated with 10 column volumes (c.v.) (here, 10 mL) of **Buffer A**, load the supernant. Wash with 10 c.v. of **Buffer A**. Then, wash with 10 c.v. of **Buffer B**, and elute with 5 c.v. of **Buffer C**, collecting the eluted fractions every  0.5 mL in 1.5 ml tubes. 1h
- 11 To quickly pool the fractions containing the protein of interest, prepare a 96-well plate or 1.5 mL tubes with  40 μL of Bradford reagent and  160 μL of distilled water. Then, add  10 μL of each protein fraction and compare against a blank reference sample corresponding to  10 μL of **Buffer C**. You can determine your protein-containing fractions either by absorbance at 595 nm on a plate reader or visually by comparing the blue coloration of each fraction against the blank reference. Pool your fractions and collect a  10 μL sample for SDS-PAGE 5m
- 12 To decrease the imidazole concentration, perform a buffer exchange step with an Amicon Ultra-15 concentrator (Merck Millipore). Centrifuge  3000 x g, 10°C, 00:10:00, discard the flowthrough, add Buffer D to decrease the imidazole concentration and repeat this step, until the imidazole concentration reaches < 30 mM. 40m
- 13 Recover the concentrated protein and determine its concentration using the Bradford assay. Also, collect a  10 μL sample for SDS-PAGE. 5m
- 14 For storage, supplement your pooled fraction with  10 millimolar (mM) DTT. Then, add glycerol up to  50 % volume to reach **Storage Conditions**. Do consider that a final protein concentration of  1 mg/mL is appropriate for subsequent experiments. 5m
- 15 Generate  200 μL aliquots of the enzyme and store it at  -20 °C until required. 30m

IMAC SDS-PAGE Result

10m

16

Expected result

