


Sep 05, 2020

Recombinant protein expression and purification of codon-optimized Bst-LF polymerase

DOI

<https://dx.doi.org/10.17504/protocols.io.bksrkwd6>



Maira Rivera¹, Severine Cazaux¹, Ariel Cerda², Anibal Arce Medina², Isaac Núñez¹, Tamara Matute¹, Alex Brown³, Javier Gasulla⁴, Fernan Federici², Cesar A Ramirez-Sarmiento¹

¹Institute for Biological and Medical Engineering, Pontificia Universidad Católica de Chile;

²Millennium Institute for Integrative Biology (iBio), Santiago, Chile; ³National Jewish Health;

⁴FCEN-Universidad de Buenos Aires, Argentina

Maira Rivera: Millennium Institute for Integrative Biology (iBio), Santiago, Chile;

Reclone.org (The Reage...

XPRIZE Rapid Covid Tes...



Cesar A Ramirez-Sarmiento

Pontificia Universidad Catolica de Chile

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bksrkwd6>

Protocol Citation: Maira Rivera, Severine Cazaux, Ariel Cerda, Anibal Arce Medina, Isaac Núñez, Tamara Matute, Alex Brown, Javier Gasulla, Fernan Federici, Cesar A Ramirez-Sarmiento 2020. Recombinant protein expression and purification of codon-optimized Bst-LF polymerase. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bksrkwd6>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol in our workspace and it is working. Typical protein yield is about 20 mg per liter of cell culture. Functionally tested in a LAMP assay using both commercial and in-house buffers equivalent to Bst2.0 from New England Biolabs and an initial protein concentration of 0.25 ug/ml (8-fold dilution from the initial 2 mg/ml stock). LAMP and RT-LAMP Assay protocols soon to be published.

Created: September 04, 2020

Last Modified: September 05, 2020

Protocol Integer ID: 41521

Keywords: RT-LAMP, isothermal amplification, COVID-19, SARS-CoV-2, If polymerase this protocol, If polymerase, purification of codon, If enzyme, recombinant protein expression, recombinant expression, purification, optimized bst, plasmid, codon, dialysis through the use,

Abstract

This protocol has been optimized for the recombinant expression of a codon-optimized Bst-LF polymerase.

The goal of this protocol was to eliminate the use of large volumes for dialysis through the use of concentrators for buffer exchange before storage conditions.

The plasmid encoding the codon-optimized Bst-LF enzyme used here can be found on reclone.org



Materials

MATERIALS

- ☒ Amicon Ultra-15 Centrifugal Filter Unit **Merck Millipore (EMD Millipore) Catalog #UFC910024**
- ☒ Sodium phosphate monobasic monohydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9638**
- ☒ PMSF **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7626**
- ☒ Sodium phosphate dibasic **Merck MilliporeSigma (Sigma-Aldrich) Catalog #7558-79-4**
- ☒ Imidazole **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I5513**
- ☒ NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**
- ☒ HisTrap FF Crude Column **GE Healthcare Catalog #17528601**
- ☒ Lysozyme **Thermo Fisher Scientific Catalog #89833**
- ☒ Glycerol **Merck Millipore (EMD Millipore) Catalog #104092**
- ☒ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #DTT-RO**
- ☒ Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #X100-100ML**
- ☒ Trizma® base **Merck MilliporeSigma (Sigma-Aldrich) Catalog #93362**
- ☒ EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #ED2SS**
- ☒ KCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9541**

Buffer A, pH 8.0

- [M] 50 millimolar (mM) NaPO₄, pH 8.0
- [M] 300 millimolar (mM) NaCl
- [M] 30 millimolar (mM) Imidazole, pH 8.0

Buffer B, pH 8.0

- [M] 25 millimolar (mM) Tris-HCl, pH 8.0
- [M] 200 millimolar (mM) KCl
- [M] 30 millimolar (mM) Imidazole, pH 8.0

Buffer C, pH 8.0

- [M] 25 millimolar (mM) Tris-HCl, pH 8.0
- [M] 100 millimolar (mM) KCl
- [M] 150 millimolar (mM) Imidazole, pH 8.0

Buffer D, pH 8.0

- [M] 25 millimolar (mM) Tris-HCl, pH 8.0



[M] 100 millimolar (mM) KCl

Storage conditions

[M] 12.5 millimolar (mM) Tris-HCl, pH 8.0

[M] 50 millimolar (mM) KCl

[M] 0.1 millimolar (mM) EDTA

[M] 0.1 % volume Triton X-100

[M] 50 % volume Glycerol

[M] 1 millimolar (mM) DTT



DAY 1 – Plasmid transformation

1d

- 1 Transform  100 ng of pET15b plasmid containing codon-optimized Bst-LF polymerase into *E. coli* C41 competent cells using either heat shock or electroporation. 2h
- 2 Spread transformed cells in LB Agar plates supplemented with  0.1 mg/mL Amp. Grow plate overnight at  37 °C . 12h

DAY 2 – Preinoculum

1d

- 3 Select a single colony from the LB agar plate to prepare a preinoculum in  10 mL LB media supplemented with  0.1 mg/mL Amp. Grow overnight at  250 rpm, 37°C . 1d

DAY 3 – Protein Overexpression

1d

- 4 Use the full volume of the preinoculum to inoculate  1 L of LB media supplemented with  0.1 mg/mL Amp (1% inoculation). Grow at  200 rpm, 37°C until reaching an optical density at 600 nm (OD_{600}) = 0.8. 4h
- 5 Upon reaching OD_{600} = 0.8, add IPTG to a final concentration of  0.5 millimolar (mM) and incubate overnight at  200 rpm, 18°C . 16h

DAY 4 – Protein Purification and Storage

6h

- 6 Centrifuge the cell culture  4000 x g, 4°C, 00:30:00 .Then, resuspend the cell pellet in  40 mL of **Buffer A** freshly supplemented with  1.0 millimolar (mM) PMSF and  0.2 mg/mL lysozyme. 30m
- 7 Incubate the resuspended cells  80 rpm, Room temperature , 00:30:00 . 30m
- 8 Sonicate on ice for  00:08:00 using cycles of  00:00:01 ON and  00:00:04 OFF at 40% amplitude (Qsonica Q125, 125W). 30m
- 9 On an ultracentrifugation tube, incubate the unclarified lysate at  65 °C for  00:25:00 to precipitate most of *E. coli* proteins, and then place on ice for 1h




00:05:00 . Centrifuge 20000 x g, 4°C, 00:30:00 and collect the supernatant.

You might want to collect a small sample for SDS-PAGE afterwards.

- 10 On a **5 mL HisTrap column (GE Healthcare)** pre-equilibrated with 10 column volumes (c.v.) (here, 50 mL) of **Buffer A**, load the supernatant. Wash with 10-20 c.v. of **Buffer B**. Then, elute with 5 c.v. of **Buffer C**, collecting the eluted fractions every 1 mL in 1.5 ml tubes. 1h
- 11 To quickly pool the fractions containing the protein of interest, prepare a 96-well plate or 1.5 mL tubes with 40 μL of 5X Bradford reagent and 160 μL of distilled water. Then, add 10 μL of each protein fraction and compare against a blank reference sample corresponding to 10 μL of **Buffer B**. You can determine your protein-containing fractions either by absorbance at 595 nm on a plate reader or visually by comparing the blue coloration of each fraction against the blank reference. Pool your fractions and collect a 10 μL sample for SDS-PAGE. 5m
- 12 To decrease the imidazole concentration, perform a buffer exchange step with an Amicon Ultra-15 concentrator (Merck Millipore). Centrifuge 3000 x g, 10°C, 00:10:00 , discard the flowthrough, add **Buffer D** to decrease the imidazole concentration and repeat this step, until the imidazole concentration reaches < 30 mM. 1h
- 13 Recover the concentrated protein and determine its concentration using the Bradford assay. Then, supplement with [IM] 0.2 millimolar (mM) EDTA, [IM] 0.2 % volume Triton X-100, [IM] 2 millimolar (mM) DTT and add glycerol up to [IM] 50 % volume to reach **Storage Conditions**. Do consider that a final protein concentration $\leq$ [IM] 2 mg/mL is appropriate for subsequent experiments. 5m
- 14 Generate 200 μL aliquots of the enzyme and store it at -20 °C until required. 30m

IMAC SDS-PAGE Result

15



SDS-PAGE at 10% PA of eluted fractions from IMAC purification of Bst-LF. The high MW contaminant is due to disulfide bond formation and eliminated upon addition of DTT (data not shown).