


Jan 04, 2023

Recombinant protein expression and purification of fuGFP

DOI

dx.doi.org/10.17504/protocols.io.e6nvwje79lmk/v1



Javiera A Avilés¹, Tamara Matute², Isaac Núñez², Maira Rivera², Javiera Reyes², Anibal Arce Medina¹, Cesar A Ramirez-Sarmiento², Fernan Federici¹

¹Millennium Institute for Integrative Biology (iBio), Santiago, Chile;

²Institute for Biological and Medical Engineering, Pontificia Universidad Católica de Chile

Reclone.org (The Reagent Collaboration Network)

Tech. support email: protocols@recode.org

[Click here to message tech. support](#)



Cesar A Ramirez-Sarmiento

Pontificia Universidad Catolica de Chile

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwje79lmk/v1>



Protocol Citation: Javiera A Avilés, Tamara Matute, Isaac Núñez, Maira Rivera, Javiera Reyes, Anibal Arce Medina, Cesar A Ramirez-Sarmiento, Fernan Federici 2023. Recombinant protein expression and purification of fuGFP. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwje79lmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 04, 2023

Last Modified: January 04, 2023

Protocol Integer ID: 74738

Keywords: RT-LAMP, isothermal amplification, COVID-19, SARS-CoV-2, recombinant expression of fugfp, encoding fugfp, fugfp, fugfp this protocol, fluorescent protein, recombinant protein expression, properties of fluorescent protein, purified protein, protein, recombinant expression, purification, plasmid, open pti vector, purification of fugfp,

Funders Acknowledgements:

ANID Millennium Science Initiative Program

Grant ID: ICN17_022

ANID CONCYTEC

Grant ID: covbio0012

Abstract

This protocol has been optimized for the recombinant expression of fuGFP encoded in an open pTi vector. The plasmid encoding fuGFP used here can be found on reclone.org. The purified protein can be used for teaching about the properties of fluorescent proteins.



Materials

MATERIALS

☒ Sodium phosphate monobasic monohydrate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9638**

☒ PMSF **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7626**

☒ Sodium phosphate dibasic **Merck MilliporeSigma (Sigma-Aldrich) Catalog #7558-79-4**

☒ Imidazole **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I5513**

☒ NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**

☒ HisTrap FF Crude Column **GE Healthcare Catalog #17528601**

☒ Lysozyme **Thermo Fisher Scientific Catalog #89833**

Buffer A, pH 8.0

[M] 50 millimolar (mM) NaPO₄, pH 8.0

[M] 300 millimolar (mM) NaCl

[M] 30 millimolar (mM) Imidazole, pH 8.0

Buffer B, pH 8.0

[M] 25 millimolar (mM) Tris-HCl, pH 8.0

[M] 200 millimolar (mM) NaCl

[M] 30 millimolar (mM) Imidazole, pH 8.0

Buffer C, pH 8.0

[M] 25 millimolar (mM) Tris-HCl, pH 8.0

[M] 100 millimolar (mM) NaCl

[M] 300 millimolar (mM) Imidazole, pH 8.0

Troubleshooting



DAY 1 – Plasmid transformation

1d

- 1 Transform  100 ng of the open pTi plasmid containing fuGFP into *E. coli* BL21 (DE3) competent cells using either heat shock or electroporation. 2h
- 2 Spread transformed cells in LB Agar plates supplemented with  0.05 mg/mL Kan. Grow plate overnight at  37 °C . 12h

DAY 2 – Preinoculum

1d

- 3 Select a single colony from the LB agar plate to prepare a preinoculum in  10 mL LB media supplemented with  0.05 mg/mL Kan. Grow overnight at  250 rpm, 37°C . 1d

DAY 3 – Protein Overexpression

1d

- 4 Use the full volume of the preinoculum to inoculate  1 L of LB media supplemented with  0.05 mg/mL Kan (1% inoculation). Grow at  200 rpm, 37°C until reaching an optical density at 600 nm (OD_{600}) = 0.8. 4h
- 5 Upon reaching OD_{600} = 0.8, add IPTG to a final concentration of  0.5 millimolar (mM) and incubate overnight at  180 rpm, 18°C 16h

DAY 4 – Protein Purification by IMAC

6h

- 6 Centrifuge the cell culture  4000 x g, 4°C, 00:20:00 .Then, resuspend the cell pellet in  40 mL of **Buffer A** freshly supplemented with  0.5 millimolar (mM) PMSF and  0.2 mg/mL lysozyme. 20m
- 7 Incubate the resuspended cells  80 rpm, Room temperature , 00:30:00 . 30m
- 8 Sonicate on ice for  00:08:00 using cycles of  00:00:01 ON and  00:00:01 OFF at 40% amplitude (Qsonica Q125, 125W). 8m 2s


- 9 Centrifuge the unclarified lysate 20000 x g, 4°C, 00:20:00 and collect the supernatant. You might want to collect a small sample for SDS-PAGE afterwards. 20m
- 10 On a **1 mL HisTrap column (GE Healthcare)** pre-equilibrated with 10 column volumes (c.v.) (here, 10 mL) of **Buffer A**, load the supernatant. Wash with 20 c.v. of **Buffer B**. Then, elute with 5 c.v. of **Buffer C**, collecting the eluted fractions every 1 mL in 1.5 ml tubes. 1h
- 11 To quickly pool the fractions containing the protein of interest, prepare a 96-well plate or 1.5 mL tubes with 40 µL of 5X Bradford reagent and 160 µL of distilled water. Then, add 10 µL of each protein fraction and compare against a blank reference sample corresponding to 10 µL of **Buffer C**. You can determine your protein-containing fractions either by absorbance at 595 nm on a plate reader or visually by comparing the blue coloration of each fraction against the blank reference. Pool your fractions and collect a 10 µL sample for SDS-PAGE. 5m
- 12 For storage, we suggest to do a dialysis against **Buffer A**, and store at 4° C.

IMAC SDS-PAGE Result

13



Eluted fractions from immobilized metal affinity chromatography (IMAC) after recombinant protein purification of fuGFP using open pTi vector.



Pooled fractions from immobilized metal affinity chromatography (IMAC) after recombinant protein purification of fuGFP using open pTi vector under blue light.