

Mar 20, 2024

Recombinant α S purification Protocol

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7popkgwz/v1>

Jane Balster¹

¹ASAP - Team Lee

Team Lee



Jane Balster

ASAP - Team Lee

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7popkgwz/v1>

Protocol Citation: Jane Balster 2024. Recombinant α S purification Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7popkgwz/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 21, 2024

Last Modified: May 31, 2024

Protocol Integer ID: 95779

Keywords: ASAPCRN, recombinant α S purification protocol, recombinant α S purification, recombinant, protocol this protocol detail, protocol detail, protocol

Abstract

This protocol details the recombinant α S purification.

Attachments



Generation of recomb...

22KB

Materials

Lysis Buffer:

	A	B
	NaOH	40 mM
	Tris pH 8.0	20 mM
	EDTA	1 mM
	Triton X-100	0.10%

Buffer A:

	Thris pH 8.0	25 mM
	NaCl	20 mM
	EDTA	1 mM

Buffer B (filtered) :

	Tris pH 8.0	25 mM
	NaCl	1M
	EDTA	1mM

Buffer C (filtered) :

	Tris pH 8	25mM
	NaCl	1M
	EDTA	1mM



Culture Growth and Induction

1d 14h 5m

- 1 **Day 1:** Inoculate  5 mL of LB/Amp [ 100 μ L :  1 μ L of  100 μ L Amp stock (freezer) for every  1 mL of culture] with 1 colony from LB/Amp plate and put in a shaker at  37 °C w/  250 rpm  Overnight .

8h



- Modification:

1. Inoculate  10 mL of LB/Amp with 1 colony (allows more rapid growth).

- If no viable plate:

1. Take a stab from glycerol stock ( -80 °C Freezer) and grow  Overnight .
2. Streak LB/Amp agar plate with  Overnight culture.

- 2 **Day 2:** Inoculate  500 mL of LB/Amp ( 100 μ L :  1 μ L of  100 μ L stock for every  1 mL of culture) with  5 mL of  Overnight culture and put in shaker at  37 °C until O.D = 0.5 – 0.6 (~1.5h -  03:00:00)

3h



- Modification:

1. Inoculate  1000 mL of LB/Amp with  10 mL of  Overnight culture and put in shaker.
2. Nano-drop:

- Use  2 μ L drop.
- Use LB or LB/Amp as a blank.

3. Formula to determine the time it will take the culture to hit an OD of 0.5 (assuming 30 min doubling time):



$$OD_f = OD_i \left[2^{(t/30min)} \right]$$

For example: for a culture with an $OD_i = 0.052$

$$0.5 = 0.052 \left[2^{(t/30min)} \right]$$

$$9.61 = \left[2^{(t/30min)} \right]$$

$$\ln[9.61] = t/30[\ln(2)]$$

$$t=97.3 \text{ minutes}$$

- 3 Induce with IPTG ( 400 μL :  6.4 μL of [M] 0.5 Mass Percent) for every  10 mL of culture) for 4h-  06:00:00 at  37 $^{\circ}\text{C}$.

6h



- 4 Spin down culture at  4600 x g, 4 $^{\circ}\text{C}$, 00:20:00 in  500 mL buckets (max volume  350 mL) in JLA-10.5 rotor (blue rotor, standing centrifuge).

20m



Note

First Stopping Point: Pellet can be frozen at  -20 $^{\circ}\text{C}$. For >  500 mL culture, split pellet into  500 mL equivalent fractions

- 5 **Day 3:** Resuspend in  50 mL of lysis buffer. Add 1 tablet of protease inhibitor cocktail (NOT MINI) and  50 μL of  500 μL PMSF ( 500 μL stock in floor  -20 $^{\circ}\text{C}$; Final conc.  500 μL ; incubate at  37 $^{\circ}\text{C}$ without shaking for 35min-  00:40:00 .

40m



Lysis Buffer:

	A	B
	NaOH	40 mM
	Tris pH 8.0	20 mM
	EDTA	1 mM



	A	B
	Triton X-100	0.1%

- 6 Add 500 μL of [M] 1 Mass Percent MgCl_2 , 500 μL of [M] 1 Mass Percent CaCl_2 , and 20 μL of DNase from Roche (200 μL total) and incubate at 37 $^\circ\text{C}$ with 250 rpm shaking for 01:00:00 hour.



1h

- 7 Add 1 mL [M] 0.25 Mass Percent EDTA, mix well, and remove cellular debris by centrifugation at 16900 x g, 00:15:00 (JA 25.50 rotor).



15m

- 8 Add 0.116 g of ammonium sulfate per mL of supernatant and stir at 4 $^\circ\text{C}$ for 01:00:00 . Centrifuge at 20000 x g for 00:30:00 , use new tubes. (~ 5.8 g ; Toss pellet).



1h 30m

- 9 Add 0.244 g of ammonium sulfate per mL of supernatant and stir at 4 $^\circ\text{C}$ for 01:00:00 (or Overnight). Centrifuge at 20000 x g for 00:30:00 , use new tubes. (~ 12.9 g ; **Keep pellet).



2h 30m

- Quality Control Checkpoint: Pellet should be at bottom of tube. Excessive proteolysis is indicated by floating pellet or significant smearing of pellet along the side of the tube.

- 10 **Day 4:** Resolubilize in 25 mL s Buffer A and add PMSF to 1 μL . Stir in refrigerator for ~ 00:30:00 .

30m

**Buffer A:**

	Thris pH 8.0	25 mM
	NaCl	20 mM
	EDTA	1 mM

- 11 Dialyze against 1 L of Buffer A with 0.5 μL PMSF Overnight .

8h





12 **Day 5:** Filter through  0.22 μL filter, Millex, Millipore.

13 LPLC: Run over Anion Exchange column HiTrap Q FF (program: SC IExAlphasS, Buffer A, filtered Buffer B), aS will start to come off at ~20% Buffer B). Keep samples at  4 $^{\circ}\text{C}$.



SEE DETAILED QFF PROTOCOL

Buffer B (filtered):

	A	B
	Tris pH 8.0	25 mM
	NaCl	1 M
	EDTA	1 mM

14 Concentrate αS positive fractions with Amicon Ultracel-3KD MWCO, Millipore, and exchange into filtered Buffer C. Perform during Day 6 equilibration.

20m

- Spin the samples at  4500 x g for ~  00:20:00 till final sample volume $\leq$  1 mL (lower volume better e.g.  750 μL).



Buffer C (filtered):

	A	B
	Tris pH 8	25 mM
	NaCl	1 M
	EDTA	1 mM

15 Day 6: LPLC: Run over size exclusion column in Buffer C. aS will start to come off at ~  80 mL . **SEE DETAILED SEC PROTOCOL – Freeze tubes until Conc./BCA.**



- 16 Concentrate aS positive fractions with Amicon Ultracel-3KD MWCO, Millipore, and exchange into sterile PBS.
- BCA samples and aliquot at desired concentration/volume (recommended at $5\ \mu\text{L}$).
 - Freeze aliquots at $-80\ ^\circ\text{C}$