



Jul 07, 2023

Regular maintenance of human pluripotent stem cells

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l69obdlqe/v1>

Narayana Yadavalli^{1,2,3,4,5}, Shawn M. Ferguson^{1,2,3,4,6,5}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Berrak Ugur

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l69obdlqe/v1>

Protocol Citation: Narayana Yadavalli, Shawn M. Ferguson 2023. Regular maintenance of human pluripotent stem cells. protocols.io <https://dx.doi.org/10.17504/protocols.io.bp2l69obdlqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 25, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82448

Keywords: hiPSCs, E8 media, Matrigel, EDTA, ASAPCRN, regular maintenance of human pluripotent stem cell, human pluripotent stem cell, passaging human ipsc, human ipsc, cell

Funders Acknowledgements:

ASAP

Grant ID: 000580

Abstract

This protocol describes the regular maintenance and passaging human iPSCs.

Attachments



[mqjubj3up.docx](#)

16KB

Materials

Reagents required

1. E8 Media (Gibco)
2. 0.5M EDTA (Gibco), working solution 0.5mM EDTA prepared in sterile PBS solution
3. Y-27632 Rock inhibitor (Tocris Bioscience)
4. Matrigel (corning)



Day 1

- 1 Pre coat a 6 well dish with Matrigel matrix for  24:00:00 or 1 hour. 1d
- 2 Thaw the frozen iPSCs by placing the vial in  37 °C water bath for  00:02:00 . 2m
- 3 After thawing, spray the tube with 70% ethanol and place in biosafety cabinet.
- 4 Aspirate the cells into 15 ml falcon tube and add  4 mL cold E8 media containing rock inhibitor. 
- 5 Spin down the tube at  0.3 rcf, 00:03:00 . 3m

- 6 Aspirate the supernatant and resuspend in fresh E8 media containing rock inhibitor.
- 7 Remove Matrigel matrix and dispense appropriate number of cells onto Matrigel coated dishes.

Day 2

- 8 Remove the E8 media containing rock inhibitor and feed cells with E8 media without rock inhibitor.
- 9 Feed every day with fresh media till plate gets 60-70% confluency.
- 10 Cells were passaged every 5 days with E8.

iPSC passaging with 0.5 mM EDTA solution: Day 1



11 Pre coat a 6 well dish with Matrigel matrix for  24:00:00 or 1 hour.

1d

12 Remove the culture media.

13 Rinse 1X with sterile PBS.



14 Add  1 mL of EDTA solution.



15 Place in the incubator for 5-7 minutes.



16 Bring the plate into biosafety cabinet and add  2 mL E8 media containing rock inhibitor.



17 Dispense all the cell suspension into 15 ml falcon tube.

18 Spin down the tube at  0.3 rcf, 00:03:00 .

3m



19 Aspirate the supernatant and resuspend in fresh E8 media containing rock inhibitor.

20 Remove Matrigel matrix and dispense appropriate number of cells onto Matrigel coated dishes.

iPSC passaging with 0.5 mM EDTA solution: Day 2

21 Remove the E8 media containing rock inhibitor and feed cells with E8 media without rock inhibitor.

22 Cells were passaged every 5 days with E8.