

Oct 16, 2019

Removal of gDNA from totalRNA

DOI

dx.doi.org/10.17504/protocols.io.8bfhsjn

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8bfhsjn>

Protocol Citation: Igem Dusseldorf 2019. Removal of gDNA from totalRNA. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8bfhsjn>

Manuscript citation:

https://assets.thermofisher.com/TFS-Assets/LSG/manuals/MAN0012000_DNase_I_RNasefree_1UuL_UG.pdf

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: October 16, 2019

Last Modified: October 16, 2019

Protocol Integer ID: 28743

Keywords: removal of gdna, removal of genomic dna, rna, genomic dna, gdna, dna, removal

Abstract

Removal of genomic DNA from RNA preparations

Guidelines

Always work with gloves and safety gear and work on ice.

Materials

- RNase-free H₂O
- DNaseI-buffer (10X)
- DNaseI
- 50 mM EDTA
- NaOAc (3M & pH 5,3) and 100% EtOH

Troubleshooting



Digestion of DNA


1

Calculate the volume required for 1 µg RNA

	Com pon ent	amo unt
	DNa sel buff er	1 µl
	total RNA	1 µg
	DNa sel	1 µl
	H2O	to 10 µl


2 pipet components into a RNase-free 1.5 ml tube (following order: H2O, buffer, total RNA, DNaseI)

3 incubate reaction at 37°C

 00:30:00

Extraction of DNaseI-digested RNA

4 add 1 µl 50 mM EDTA and incubate at 65°C


 00:10:00

5 add 1/10 Volume NaOAc and 3 Volumes 100% EtOH

6 freeze at -20°C over night or at -80°C for at least 1h



7 centrifuge at 13000 rpm at 4°C


 00:30:00

8 Discard supernatant

RNA pellet may be invisible, make sure to align tube lid so you remember where the pellet is located

9 add 500 µl 70% EtOH

10 centrifuge at 13000 rpm at 4°C

 00:10:00

11 discard supernatant, remove as much ethanol as possible

12 dry tubes at RT under the hood (~5min), do not overdry!

13 resuspend pellet in 20 µl RNase-free water

14 measure concentration

Quality control:

15 perform PCR on your DNA-free RNA. Include a positive control using gDNA as a template and primer that bind to it.

There should be no visible bands except for positive control!