

Mar 22, 2019

## Removal of genomic DNA from RNA preparations (Thermo Scientific )- (M4455 Version)



Forked from [Removal of genomic DNA from RNA preparations \(Thermo Scientific \)](#).

DOI

[dx.doi.org/10.17504/protocols.io.ze2f3ge](https://dx.doi.org/10.17504/protocols.io.ze2f3ge)

Maximilian Dietsch<sup>1</sup>

<sup>1</sup>Heinrich-Heine Universität Düsseldorf

Axmann Lab

M4455 - Synthetische B...



Sebastian ST Triesch

Institute for Synthetic Microbiology, HHU Düsseldorf

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.ze2f3ge](https://dx.doi.org/10.17504/protocols.io.ze2f3ge)

**Protocol Citation:** Maximilian Dietsch 2019. Removal of genomic DNA from RNA preparations (Thermo Scientific )- (M4455 Version). [protocols.io](https://dx.doi.org/10.17504/protocols.io.ze2f3ge) <https://dx.doi.org/10.17504/protocols.io.ze2f3ge>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 22, 2019

**Last Modified:** March 22, 2019

**Protocol Integer ID:** 21690

## Abstract

Removal of genomic DNA from RNA preparations

## Guidelines

allways waer gloves and work on ice



- 1

Add to an RNase free tube:	
RNA	1 µg
10X reaction buffer with MgCl <sub>2</sub>	1 µl
DNase I, RNase-free	1 µl (1U)
Water	to 10 µl
- 2 Incubate at 37 °C for 30 min

37 °C

00:30:00
- 3 Add EDTA, Water and PCI and vortex thoroughly.

1 µL EDTA

80 µL Water

100 µL PCI (phenol chloroform isoamyl alcohol)
- 4 Centrifuge for 10 min at 10000 rpm and 4 °C
- 5 transfer the upper phase into a fresh tube and add 3 volumen EtOH/ 3M Natrumacetat (30:1, ph 5.2 )
- 6 precipitate RNA over night at -20 °C
- 7 Centrifuge 30 min at 13000 rpm and 4 °C

4 °C

00:30:00
- 8 Discard supernatant and wash pellet with 75% EtOH ( do not resuspend the pellet)
- 9 Centrifuge 10 min at 13000 rpm and 4 °C

4 °C

00:10:00
- 10 go to step #8 repeat washing step



 [go to step #9](#) Centrifuge

- 11 Discard supernatant and dry pellet for 10 - 15 min
  
- 12 resuspend pellet with 30  $\mu$ l H<sub>2</sub>O