



Apr 18, 2023

Version 2

Reproductive Tissue Collection (Mammals): Post-mortem Sampling V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7yrb3gwz/v2>

sanaz.arenivas¹, Barbara Durrant², comizzolip³, mhouch², Jacquelyn Mountcastle⁴, Budhan Pukazhenth⁵, phil.purdy⁶

¹Viagen Pets, Cedar Park, TX, USA;

²Conservation Science Wildlife Health, San Diego Zoo Wildlife Alliance, Escondido, CA, USA;

³Smithsonian Conservation Biology Institute, National Zoological Park, Front Royal, VA, USA;

⁴Revive & Restore, Sausalito, CA, USA;

⁵Smithsonian National Zoological Park and Conservation Biology Institute, Front Royal, VA, USA;

⁶USDA, ARS, NLRP, National Animal Germplasm Program, Fort Collins, CO, USA



Revive and Restore

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7yrb3gwz/v2>



Protocol Citation: sanaz.arenivas , Barbara Durrant, comizzolip , mhouck , Jacquelyn Mountcastle, Budhan Pukazhenth, phil.purdy 2023. Reproductive Tissue Collection (Mammals): Post-mortem Sampling. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.q26g7yrb3gwz/v2> Version created by **Revive and Restore**

Manuscript citation:

Brahmasani, S.R., Sontakke, S., Ghosh, S., Mallapur, G., Yadav, S., Vasudevan, K. (2021). GRB Manual- An Introduction to Genetic Resource Banks for Wildlife Conservation.

https://cza.nic.in/uploads/documents/publications/english/An%20Introduction%20to%20Genetic%20Resource%20Banks%20for%20Wildlife%20Conservation_digital%20copy.pdf

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: March 21, 2023

Last Modified: April 18, 2023

Protocol Integer ID: 79152

Keywords: biobanking, cryopreservation, cryobanking, reproductive tissue, gametes, germplasm, testes, ovaries, reproductive tissue collection, collected tissue, living reproductive tissue, term cryopreservation, tissue, mammal, gamete collection, collection

Abstract

Version date: April 2023

The following protocol illustrates how to collect and ship living reproductive tissue from a deceased wild or captive mammal for long-term cryopreservation. Collected tissues can either be immediately cryobanked at -196°C or processed for gamete collection and later cryobanked at -196°C.

To print this protocol: Click the arrow to the right of the "MORE" button at the top of this page and select the option to print or download as PDF version.

Guidelines

The collection of samples should be opportunistic and follow all applicable regulations. Reproductive tissues collected from deceased animals should be collected **as close to the time of death as possible** to avoid tissue decay. Temperature, body condition, and time are factors that accelerate or decelerate tissue decomposition, dictating the chances of collecting living cells.

If the animal dies suddenly or if the animal must be euthanized, the carcass or any tissues **must not be frozen**. The carcass can be kept in a cool area and reproductive tissues must be dissected immediately and shipped the same day as collection to a biobanking facility.

Testes maximum time after collection before arrival at facility: **24 hours, at 4-15°C**

If testicular tissue is warm at time of collection, it should be kept at room temperature during the collection process. If the tissue is cold, it must be kept cold throughout the process to avoid rewarming and recooling. Testicular tissues must be kept at 4°C during shipment.

Ovaries maximum time after collection before arrival at facility: **4 hours, at 15-20°C**

Ovaries held at low temperatures (less than 15°C) for more than 24 hours will begin to rapidly degrade the oocytes inside. Ovarian tissue must be kept at 15-20°C at all times. Contact receiving facility to discuss transport of samples prior to collection.

For any animal carcass found in the wild, time, ambient temperature, and storage methods are critical factors that can impact how quickly a sample must be collected. If the time of death is unknown, tissues can be harvested as long as it has not been frozen, become necrotic, or started decomposing. In such cases, tissues should still be shipped immediately.

For urgent questions regarding protocol steps or collection guidelines, please contact the Revive & Restore Biobanking Team at informedbiobanking@reviverestore.org



Materials

User-supplied Materials:

Sample containers

Phosphate buffered saline (PBS) with antibiotics (Penicillin, Streptomycin, Gentamicin)

Sterile scalpel

Sterile scissors

Disposable gloves

Gauze

Rubbing alcohol

Field notebook

Pencil/Pen

Permanent/alcohol-resistant marker

Newspaper or bubble wrap

Ice packs (*testicular tissue only*)

Styrofoam shipping box

Sample metadata form

Pre-paid FedEx shipping label

Safety warnings

- ⚠ Steps for personal safety must be considered before going into the field. Some wildlife diseases are transmissible to humans. Refer to your agency's health and safety guidelines for personal protective equipment (PPE). At a minimum, field personnel should wear disposable gloves and a fresh pair should be used between handling different specimens to avoid cross contamination. Sampling instruments and equipment should be thoroughly cleaned and disinfected or disposed of after use.



Before start

Tissue Handling

Collected specimens are stored in vials containing a nutrient-rich media to keep the cells alive and a mixture of antibiotics to prevent bacterial growth. All samples must be collected under aseptic conditions to avoid contamination. Using sterile tweezers, scalpels, and scissors, as well as cleaning the sampling site with alcohol or surgical scrub (e.g. chlorhexidine) will decrease the chances of contamination. Collect testes and/or ovaries in a sterile container large enough to hold the entire organ, such as a large 50ml centrifuge tube, a clean Tupperware, or ziplock bag.

Cold Storage

Avoiding temperature fluctuations is very important for preserving the tissue samples. Use ice packs for transport of *testicular tissues* **only** from the field and ensure refrigeration is available immediately after returning from the field site. Do not use wet ice or other commercial therapeutic packs. Chill ice packs in the freezer the night before collection. Frozen ice packs will remain cold in an insulated container for up to 24 hours. If the field site is more than 24 hours from refrigeration, plan to bring a portable electric refrigerator. Vials containing tissue should not have direct contact with the ice packs. **Never freeze tissue biopsies** before shipping to a biobanking facility.

Shipping

All testicular tissue samples must be shipped to a biobanking laboratory via overnight express. Shipment Monday through Wednesday will guarantee arrival at the receiving facility the next day. **Do not ship samples on Friday.** Immediate transport of ovaries to the receiving facility is highly advised. Overnight shipment of ovarian tissue is not recommend but may be done if it is the only option available. Consult with the receiving facility for any questions regarding shipping timelines before collecting or sending samples.

Preparation

- 1 Pre-chill icepacks in the freezer the day before planning to collect and ship samples (*testicular tissues only*).
- 2 Record all information about the species and individual, including a picture of the animal for identification and GPS location where the animal was found.
- 3 Proper protective equipment must be worn (gloves, etc). Sterility must be maintained as much as possible.
- 4 Use a sterile container large enough to hold the entire organ, such as a large 50ml centrifuge tube, a clean Tupperware with a tight-sealing lid, or a durable ziplock bag. Gather containers and have them readily available.

STEP CASE

Ovary Collection

9 steps

The figure below demonstrates tissue decay of a reproductive tract between time of death and removal from the body.

(A) Two hours after death (four-horned antelope). This is an ideal sample for biobanking.;

(B) Thirteen hours after death (gaur). Visible black patches indicate necrotic change (cell death);

(C) Eighteen hours after death (marmoset). Further black discoloration indicates necrotic change.




Photo credit: Brahmasani et al., 2021

- 5 Wet the entire sampling area of the carcass in rubbing alcohol and blot with sterile gauze.
- 6 With sterile scissors or scalpel, dissect open the pelvic cavity. The reproductive tract will lie under the large intestine/rectum.



- 7 Collect the entire reproductive tract, including both ovaries, the uterine body, and oviduct. For larger animals, collect only the ovaries.

Tissue Preservation

- 8 Thoroughly soak gauze in phosphate buffered saline (PBS) and antibiotics and gently wrap around the tissue. Ensure that the gauze fully covers the tissue, keeping it moist. Place in sample container and close the cap tightly. If the specimen is too large for the container, you may use a clean ziplock bag.
- 9 Using an alcohol-resistant marker, label the container with an identifier that matches **exactly** the recorded metadata. Check to make sure that each container is easy to identify with the information provided.
Required information to record:
 - Scientific name of animal
 - Sex of individual denoted as ♀ (female)
 - Date of tissue collection
 - Tissue type
 - Any other identification number of individual
- 10 Ovarian tissue must be kept at 15-20°C. **Do not refrigerate or freeze samples.**

*****Notify the receiving facility immediately.*****

Transport

- 11 Place the containers with samples in an insulated shipping box. Add additional 3 inches of bubble wrap or newspaper between the samples to ensure that the samples do not move or become damaged during transport.
- 12 Enclose all paperwork and sample/animal information in the insulated shipping container.
- 13 Contact the receiving facility to arrange immediate pick up or shipment.



Protocol references

Brahmasani, S.R., Sontakke, S., Ghosh, S., Mallapur, G., Yadav, S., Vasudevan, K. (2021). GRB Manual- An Introduction to Genetic Resource Banks for Wildlife Conservation.

Gaur, A., Umapathy, G., Vasudevan, K., Sontakke, S., Rao, S., Goel, S., Kumar, D., Gupta, B., Singh, D.N., Kumar, A. (2017). Manual for Biological Sample Collection and Preservation for Genetic, Reproductive and Disease Analyses.

Trans Ova Genetics. (2022). Trans Ova Livestock Cloning. <https://transova.com/service/cloning-services/>

ViaGen Pets & Equine. (2021). Emergency Pets Protocol – U.S. Clients. <https://www.viagenpets.com/emergency-pets/>

White, C. L., & Dusek, R. J. (2015). Wildlife specimen collection, preservation, and shipment. In J. C. Franson, M. Friend, S. E. J. Gibbs, & M. A. Wild (Eds.), Techniques and Methods. US Geological Survey.
<https://doi.org/10.3133/tm15c4>