

Jun 05, 2023

Retina and RPE/Choroid RNA Extraction Protocol

DOI

dx.doi.org/10.17504/protocols.io.5qpvorjzdv4o/v1

Chottiwatt Jittprasong¹

¹City University of Hong Kong



Chottiwatt Jittprasong

City University of Hong Kong

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvorjzdv4o/v1>

Protocol Citation: Chottiwatt Jittprasong 2023. Retina and RPE/Choroid RNA Extraction Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvorjzdv4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
In Use

Created: June 04, 2023



Last Modified: June 05, 2023

Protocol Integer ID: 82864

Keywords: choroid rna extraction protocol this protocol, choroid rna extraction protocol, extraction of rna, retinal pigment epithelium, retina, rna, trizol as the main reagent, utilizing trizol, rpe, cell lysi, extraction, reagent

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol involves the extraction of RNA from retinal and retinal pigment epithelium (RPE) cells in the subject's retina, utilizing TRIzol as the main reagent for cell lysis.

Troubleshooting

Safety warnings

- ❗ TRIzol is extremely corrosive and toxic. Exercise extreme caution while handling TRIzol (e.g. handling within the fume hood).

Before start

Prior to the initiation of this protocol, extraction of the retinal sample from the subject's eyeball is required.



RNA Extraction

1h 48m

- 1 Precool the centrifuge to 4 °C
- 2 Add 1000 µL of **TRIZOL** per retina and pipette up and down 30 times to homogenize the mixture.

Safety information

This step should be done in fume hood due to toxicity of TRIZOL.

- 3 Incubate for 00:05:00 On ice 5m
- 4 Add 200 µL of **Chloroform** per 1000 µL of TRIZOL used for lysis and mix thoroughly.

- 5 Incubate for 00:03:00 On ice 3m

- 6 Centrifuge the sample at 12000 x g, 4°C, 00:15:00 15m

- 7 Transfer the aqueous phase containing the RNA to a new tube.

Note

Avoid penetration of the interphase layer.

- 8 Add 500 µL of **Isopropanol** to the aqueous phase, per 1000 µL of TRIZOL used.

- 9 Add 4 µL of **Glycoblue** and mix well.



10 Incubate at -4 °C for 01:00:00

1h



11 Centrifuge at 30000 x g, 4°C, 00:15:00

15m



Expected result

Blue pellet should be present at the bottom of the tube.

12 Discard supernatant carefully.



13 Resuspend the pellet in 1000 µL of On ice 75% ethanol per 1000 µL of TRIzol. Vortex and centrifuge at 7500 x g, 4°C, 00:05:00

5m



14 Repeat **Step 13**.

15 Air dry the RNA pellet for 00:05:00

5m

16 Resuspend the pellet in 20-50 µL of RNase free water.



Note

40 µL of RNase free water is recommended.

17 Quantify the sample by NanoDrop

