

Nov 20, 2018

Version 1

## rev-ChIP V.1

DOI

[dx.doi.org/10.17504/protocols.io.vp5e5q6](https://dx.doi.org/10.17504/protocols.io.vp5e5q6)

Lorane Texari<sup>1</sup>, Carlos Guzman<sup>1</sup>, Sven Heinz<sup>1</sup>

<sup>1</sup>University of California, San Diego



Carlos Guzman

UCSD

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.vp5e5q6>

**Protocol Citation:** Lorane Texari, Carlos Guzman, Sven Heinz 2018. rev-ChIP. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.vp5e5q6>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol in our group and it is working.

**Created:** November 18, 2018



**Last Modified:** November 20, 2018

**Protocol Integer ID:** 17885

**Keywords:** ChIP-seq, genomic regulatory network, characterization of genomic regulatory network, precise regulation of transcriptional program, in vivo genome, wide location of transcription factor, genome, genomic, seq experiment, transcription factor, chromatin accessibility, nascent transcription, transcriptional program, histone modification, chromatin immunoprecipitation, transcription, novel chip, regulatory network, chip, antibody, seq

## Abstract

Understanding the precise regulation of transcriptional programs in human health and disease requires the accurate identification and characterization of genomic regulatory networks. Next-generation sequencing (NGS) technologies are powerful, and widely applied tools to map the in vivo genome-wide location of transcription factors (TFs), histone modifications, chromatin accessibility, and nascent transcription that make up these regulatory networks. While chromatin immunoprecipitation followed by sequencing (ChIP-seq) is one of the oldest, and most-utilized experimental techniques to study the location and abundance of TFs, experiments still frequently require optimization to reproducibly yield good data with high signal-to-noise ratios due to the massive variability between possible antibody-antigen combinations and commercial reagents.

To overcome these obstacles, we systematically carried out well over 500 ChIP-seq experiments designed to test every aspect of typical ChIP-seq experiments and developed rev-ChIP, a novel ChIP-seq method that is optimized for scalability, robustness, low-input, speed, cost efficiency and data quality. We find that rev-ChIP can be scaled to work for cell numbers ranging from millions to under a thousand, and from a single sample to 500 samples a week in a non-automated fashion with minimal hands-on time. Additionally, rev-ChIP has been tested on a variety of sample types ranging from cell lines to sorted primary cells and solid tissues.

## Troubleshooting



## Lysis and Sonication

- 1 Thaw cell pellet on ice and resuspend cells in  $500\ \mu\text{L}$  of lysis buffer.

### Note

1We do not recommend using  $500\ \mu\text{L}$  of lysis buffer or tip sonication when sonicating less than 500K cells. In this case we suggest using Covaris (  $130\ \mu\text{L}$  ) or PIXUL (  $60\ \mu\text{L}$  ).

- 2 Sonicate samples for 7 cycles.

### Note

This step is dependent on crosslinking method, and cell line or tissue type and should be optimized.

### Note

Double crosslinked DNA is harder to sonicate and requires more rounds of sonication.

### Note











Make sure not to oversonicate your samples and keep them constantly cold.

### Expected result

Chromatin size range of 200-500.



## Input Cleanup


- 3 Take  10  $\mu\text{L}$  of each sample (for input) and put them into new PCR strip.
- 4 Dilute lysis buffer (LB3 - add 55ul 10% Triton X-100, Metivier - add 750ul Metivier Dilution Buffer, RIPA - None). If splitting lysates for IP, make sure diluted samples are well mixed. Split lysates as required if needed.
- 5 Add  68  $\mu\text{L}$  of Elution Buffer (10mM Tris pH8, 0.5% SDS, 5mM EDTA, 280mM NaCl) +  1  $\mu\text{L}$  RNase A to each input sample and incubate for  00:15:00 at  37 °C .
- 6 Add  1  $\mu\text{L}$  55 of Proteinase K to each input sample and incubate at  55 °C for  01:00:00 and then at  65 °C for  00:30:00 .

## Immunoprecipitation

- 7 Prepare Dynabeads A/G: capture Dynabeads on magnet, remove supernatant and resuspend in equal volume of appropriate lysis buffer.

### Note

For LB3 use LB3 + 1/9th volume of 10% Triton X-100, for Metivier use an equal volume of Metivier Dilution Buffer and for RIPA use an equal volume of RIPA buffer.






- 8 Add the appropriate volume of beads + antibody to each ChIP sample.
- 9 Incubate IP overnight on wheel at  4 °C (rotating at 8rpm).

### Note

Optionally, you can incubate IP at  4 °C for  01:00:00 in most cases.

## Input Cleanup



- 10 Create mastermix of  2  $\mu\text{L}$  1 SpeedBeads +  120  $\mu\text{L}$  120% PEG8000/1.5M NaCl (8.5% PEG, 1M NaCl), mix thoroughly, and add  122  $\mu\text{L}$  of mastermix to each input sample.
- 11 Incubate at RT for  00:10:00 .
- 12 Wash 2x with  200  $\mu\text{L}$  of 80% EtOH on magnet (move strip side to side 10x slow, 10x fast).

**Note**

We suggest using a repeater pipette here to speed things up.


**Note**

If magnet not available we suggest vortexing beads.

- 13 Air-dry until cracks appear on bead pellet.

**Note**

Air-drying should take approximately  00:14:00 .

- 14 Elute in  15  $\mu\text{L}$  TT (0.05% Tween 20, 10mM Tris pH 8.0, cold).

**Note**

We suggest using a repeater pipette here to speed things up.



#### Note

Can store at -20 °C and stop for the day.

## Washes

- 15 Add 100x PIC to WBI/III and TET ( 10  $\mu\text{L}$  PIC per 1 mL WBI/WBIII/TET). Do 3x washes with 180  $\mu\text{L}$  WBI + PIC, 3x washes with 180  $\mu\text{L}$  WBIII + PIC, and 2x washes with 185  $\mu\text{L}$  cold TET + PIC.

#### Note

We suggest using a multichannel pipette here to speed things up.

#### Note



















The addition of PIC (Protease Inhibitor Cocktail) is not required, but highly recommended.

- 16 Resuspend beads in 25  $\mu\text{L}$  of cold TT using repeater pipette.

## Library Preparation

- 17 Collect beads and take 2  $\mu\text{L}$  of each input supernatant (1-2ul for 500K cells) that will be library prepped and add 23  $\mu\text{L}$  of TT to each input taken.








- 18 Create a mastermix of  1.5  $\mu\text{L}$  of Enzyme Mix End Prep +  3.5  $\mu\text{L}$  of End Prep Reaction Buffer per sample, mix well and add  5  $\mu\text{L}$  of mastermix to each sample. Incubate for  00:30:00 at  20  $^{\circ}\text{C}$  and then  00:30:00 at  65  $^{\circ}\text{C}$  .
- 19 Add  1  $\mu\text{L}$  of Bioo ChIP Adaptors (10.625uM) to each sample.
- 20 Create a mastermix of  15  $\mu\text{L}$  Ligation Master Mix +  0.5  $\mu\text{L}$  of Ligation Enhancer per sample, mix well, and add  15.5  $\mu\text{L}$  of mastermix to each sample and incubate for  00:15:00 at  20  $^{\circ}\text{C}$  .
- 21 Create a STOP solution mastermix of  4  $\mu\text{L}$  10% SDS +  3  $\mu\text{L}$  0.5M EDTA +  20  $\mu\text{L}$  water per sample and add  27  $\mu\text{L}$  of mastermix to each sample.
- 22 Add  4.5  $\mu\text{L}$  of 5M NaCl to each sample.

**Note**


We recommend using a multichannel to speed things up.

## Proteinase K and Reverse Crosslinking



- 23 Add  1  $\mu\text{L}$  of Proteinase K to each sample using multichannel and incubate for  01:00:00 at  55 °C and then  00:30:00 at  65 °C .





#### Note

Alternatively, FA-fixed can be incubated for  00:30:00 at  55 °C .

#### Note

Can leave at  4 °C overnight.

## Cleanup


- 24 Create mastermix of  2  $\mu\text{L}$  SpeedBeads +  61  $\mu\text{L}$  20% PEG8000/1.5M NaCl (8.5% PEG, 1M NaCl), mix thoroughly, and add  63  $\mu\text{L}$  of mastermix into new tube strips, then collect and transfer the supernatant of samples into these new tubestrips with  63  $\mu\text{L}$  of speedbeads + peg.

#### Note

DO NOT MIX PROTEIN A/G BEADS WITH SPEEDBEADS

- 25 Incubate at RT for  00:10:00 .



- 26 Wash 2x with  200  $\mu\text{L}$  of 80% EtOH on magnet (move strip side to side 10x slow, 10x fast).


#### Note

We recommend using a repeater pipette to speed things up.

- 27 Air-dry until cracks appear on bead pellet.

#### Note

Air-drying should take approximately  00:14:00 .

- 28 Elute in  25  $\mu\text{L}$  TT (0.05% Tween 20, 10mM Tris pH 8.0, cold). Collect beads on magnet and transfer supernatant into new PCR strips.

#### Note

We recommend using a multichannel to speed things up.

## Library Prep Amplification PCR






- 29 Do PCR.

	Mastermix Library PCR		
	<b>25.5ul MM + 24.5ul sample</b>	<b>1x</b>	<b>20x</b>
	Sample	24.5	
	(Blue Cap) NEBNext Ultra II Q5 2x MM	25	500
	100uM Solexa 1GA	0.25	5
	100uM Solexa 1GB	0.25	5



	PCR Program		
		98°C	30 seconds
	R	98°C	10 seconds
	60°C	15 seconds	
	72°C	30 seconds	
		72°C	1 minute
		4°C	ON

## Final Cleanup

- 30 Make mastermix of  2  $\mu\text{L}$  SpeedBeads +  38.5  $\mu\text{L}$  of per sample and add  40.5  $\mu\text{L}$  of mastermix to each sample.
- 31 Incubate at RT for  00:10:00 .
- 32 Wash 2x with  200  $\mu\text{L}$  of 80% EtOH on magnet (move strip side to side 10x slow, 10x fast).

### Note

We recommend using a repeater pipette to speed things up.

- 33 Air-dry until cracks appear on bead pellet.

### Note

Air-drying should take approximately  00:14:00 .

- 34 Elute in  20  $\mu\text{L}$  of TT.