


Sep 08, 2025

REVA #3: Magnetic Resonance Imaging (MRI) of Embalmed Cadaver



DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvrw9ezlmk/v1>

Daniel A Herzka¹, Michael Markley², Noa B Nuzov³, Goksel Sali¹, Shruti Kumari¹, Nicole A Pelot⁴, Andrew J. Shoffstall^{3,5}, Chris Flask¹, Andrew R. Crofton^{6,7}, Ari Blitz¹

¹Department of Radiology, Case Western Reserve University and University Hospitals, Cleveland Medical Center, Cleveland, OH, USA, 44106;

²University Hospitals Cleveland Medical Center, Cleveland, OH, USA 44106;

³Department of Biomedical Engineering, Case Western Reserve University, Cleveland, OH, USA, 44106;

⁴Department of Biomedical Engineering, Duke University, Durham, NC, USA, 27708;

⁵APT Center, Louis Stokes Cleveland Department of Veterans Affairs Medical Center, Cleveland, OH;

⁶Department of Anatomy, Case Western Reserve University, Cleveland, OH;

⁷Department of Pathology and Cell Biology, University of South Florida, Tampa, FL

Daniel A Herzka: ORCID: 0000-0002-9400-7814;

Noa B Nuzov: ORCID: 0000-0001-8187-2115;

Goksel Sali: ORCID: 0009-0009-1638-8509;

Shruti Kumari: ORCID: 0000-0002-3077-1982;

Nicole A Pelot: ORCID: 0000-0003-2844-0190;

Andrew J. Shoffstall: ORCID: 0000-0002-0881-2180;

Andrew R. Crofton: ORCID: 0000-0002-1105-3971



Nicole A Pelot

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvrw9ezlmk/v1>

Protocol Citation: Daniel A Herzka, Michael Markley, Noa B Nuzov, Goksel Sali, Shruti Kumari, Nicole A Pelot, Andrew J. Shoffstall, Chris Flask, Andrew R. Crofton, Ari Blitz 2025. REVA #3: Magnetic Resonance Imaging (MRI) of Embalmed Cadaver. [protocols.io https://dx.doi.org/10.17504/protocols.io.6qpvrw9ezlmk/v1](https://dx.doi.org/10.17504/protocols.io.6qpvrw9ezlmk/v1)

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 15, 2025

Last Modified: September 08, 2025

Protocol Integer ID: 225834

Keywords: Magnetic resonance imaging, MRI, Human anatomy, Gross anatomy, Embalmed cadaver, embalmed cadaver this protocol, embalmed human cadaver, embalmed cadaver, human cadaver, mri, cadaver, magnetic resonance imaging, parameters for magnetic resonance imaging, imaging, reva, procedure

Funders Acknowledgements:

NIH SPARC REVA

Grant ID: 75N98022C00018

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol describes the procedures and parameters for magnetic resonance imaging (MRI) of embalmed human cadavers.

Guidelines

N/A

Materials

1. MRI-safe gurney
2. 3 T MRI Scanner (MAGNETOM Vida, Siemens Healthineers, Erlangen, Germany)
3. Transfer board ([Link](#))
 - a. Manufacturer: Allied Healthcare Products, sold on Medline. Name: XTRA Emergency Backboards.
4. MRI 20-channel head coil (Head/Neck 20 , Siemens Healthineers)
5. MRI 36-channel spine coil (Siemens Healthineers)
6. MRI 18-channel anterior body coils (x2) (Body Array, Siemens Healthineers)
7. Disinfectant wipes

Troubleshooting

Safety warnings

-  Personnel who will enter the MRI suite must be properly trained MR safety and approved by the local safety committee. Cadavers with significant metal implants (e.g., pacemakers, may not be optimal for this protocol due to potential image artifacts.

This protocol might include items and/or substances that may pose hazards (e.g., chemical, physical, biological, or otherwise) to your health upon use or exposure. Before engaging in the processes described in this protocol, familiarize yourself with and follow the safety data sheets, manufacturer safety recommendations, and local regulations.

Ethics statement

Be sure to seek approval for or an exemption from human subjects research from your local regulatory body(ies) as required by local and/or institutional regulations before initiating studies.

This study was determined to be exempt from IRB oversight by the Case Western Reserve University Institutional Review Board (IRB) because it involved de-identified cadaveric tissue and no protected health information was collected from the donors.

Before start

See the protocol for embalming and preparing a human cadaver for imaging (dx.doi.org/10.17504/protocols.io.kxygx4wm4l8j/v1), including placing it in a vacuum-sealed bag on a dissection cart for transport. On the day of the scan, before transporting the cadaver, ensure the vacuum seal bag is still tight and no air has entered the bag; if air has entered the vacuum seal bag, use an electric vacuum to remove it.



Section 1: Cadaver Placement in the MR Scanner

- 1 Move the cadaver (in a vacuum-sealed mattress bag) from the dissection cart used for transportation onto the MRI-safe gurney while outside the scan room. An MR-compatible patient transfer board can be used to facilitate transfer.
- 2 Place the cadaver on the MRI scanner table, positioning the head as far as possible inside of the head coil which should result in the shoulders coming in contact with the edge of the head coil.
- 3 The spine coil is already in place as it is built in with the table.
- 4 Place two (2) body coils over the chest and abdomen of the cadaver. Depending on the cadaver height, the rectangular coils can be placed with the longest edge in the left-right or head-foot direction as needed for coverage. Use straps to secure the coils to the table and minimize distance between coils and cadaver. The straps should also press the arms against the torso to minimize the required field of view during imaging. Coils should be abutted and there should be minimal overlap between coils. Ensure the cadaver is secured and that the arms are pressed against the torso.
- 5 Landmark the position of the cadaver using scanner controls using the nose as a point of reference.

Section 2: MRI Scan Session

- 6 Run a fast localizer scan spanning the entire volume to be imaged to visualize the general positioning of the cadaver. Reposition the cadaver as needed to straighten the spine and neck.
- 7 Run a higher resolution localizer, based on the HASTE sequence. Use the multi-station approach to ensure sufficient coverage from the top of the head to the pubic symphysis. Overlap each station by 15% to facilitate stitching of images from different stations into a single DICOM image with full coverage.
- 8 Run the Dixon-VIBE scan as defined in **Table 1**. Ensure inclusion of acromia for landmarking purposes.

	A	B	C	D
		CISS	Dixon-VIBE	VIBE

A	B	C	D
Orientation	Sagittal	Axial	Sagittal
Minimum TR	5.99 ms	10.00 ms	10.00 ms
TE	3.00 ms	Echo 1: 2.46 ms Echo 2: 3.69 ms	2.46 ms
Flip angle	35 deg	10 deg	12 deg
FOVx	256 mm	405 mm	323 mm
FOVy*	256 mm	278 mm	212 mm
FOVz*	224 mm	145.6 mm	128 mm
Matrix (X, Y *)	512, 512	576, 396	640, 420
Number of slices	448	208	256
Oversampling in z	0.0%	15.4%	12.5%
Oversampling in y	0.0%	10%	10%
Elliptical scanning window	on	off	on
Scan duration per station	35 min 57 sec	10 min 28 sec	14 min 23 sec
Number of stations**	4	7	3
Bandwidth	488 Hz/Px	430 Hz/Px	240 Hz/Px
Averages	1	1	1
Orientation	sagittal	axial	sagittal

A	B	C	D
Image resolution (acquired and reconstructed voxel size)***	Acquired: 0.5 mm isotropic Reconstructed: 0.5 mm isotropic	Acquired: 0.7 mm isotropic Reconstructed: 0.7 mm isotropic	Acquired: 0.5 mm * 0.6 mm * 0.6 mm Reconstructed: 0.5 mm isotropic
Dimensionality (2D vs 3D)	3D	3D	3D

Table 1. MRI scan parameters.

* Adapted for full anterior-posterior coverage dependent on body habitus; resolution is fixed at 0.5 mm.

** Adapted depending on the height of the subject.

*** 6/8 partial Fourier under sampling in both slice and phase for Dixon-VIBE

- 9 Run the CISS scan as defined in **Table 1**.
- 10 Run the VIBE scan as defined in **Table 1**. Ensure the number of stations and fields of view are exactly matched to CISS scan.
- 11 Stitch the images from sequential stations from scans in **Steps 8-10** into a single image using the vendor-provided stitching software.
- 12 Reformat all data (multiplanar reconstruction, MPR) to convert the sagittally acquired images into axial images.
- 13 Remove the cadaver from the scanner table to the MRI-compatible gurney using the transfer board.
- 14 Move the cadaver from the gurney to the dissection cart for removal from the MRI suite.
- 15 Use disinfectant wipes to clean the MR scanner table, control panels, and imaging coils.

Acknowledgements

The authors thank the donors and staff of the Case Western Reserve University Anatomical Gift Program. Without the selfless donations of the donors and the tireless efforts of the staff of this program, the methods described in this protocol and the insights and advancements that result from studies conducted according to this protocol would not be possible.