

Aug 09, 2023

Reverse-phase high pH fractionation (using Thermo Fisher, Cat# 84868)

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm32bnl3p/v1>

Louise Uoselis¹

¹Lazarou Lab, WEHI



Louise Uoselis

WEHI

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm32bnl3p/v1>

Protocol Citation: Louise Uoselis 2023. Reverse-phase high pH fractionation (using Thermo Fisher, Cat# 84868). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.yxmvm32bnl3p/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 09, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 86212

Keywords: ASAPCRN, phase peptide fractionation kit, phase high ph fractionation, high ph fractionation, using thermo fisher, phase, thermo fisher

Abstract

This protocol uses the Pierce™ High pH Reversed-Phase Peptide Fractionation Kit (Thermo Fisher, Cat# 84868)



Conditioning the columns

4m

- 1 Remove the white cap on the end of the column and place the column in a 2 mL collection tube.
- 2 Centrifuge at 5000 rcf, 00:02:00 at Room temperature , discard the liquid. 2m
- 3 Remove the screw cap and add 300 μ L acetonitrile (ACN) to the column (replacing the screw cap after).
- 4 Centrifuge at 5000 rcf, 00:02:00 at Room temperature , discard the liquid. 2m
- 5 Repeat steps 3 and 4 (for a total of 2 washes with ACN).
- 6 Repeat steps 3 – 5 with 0.1% TFA instead of ACN (total of 2x washes with 0.1% TFA).
- 7 The column is now ready to use.

Fractionating the samples

16m 20s

- 8 Add 300 μ L of 0.1% v/v trifluoroacetic acid to each sample.
- 9 Vortex for ~ 00:00:10 10s
- 10 Leave to incubate for 00:05:00 at Room temperature 5m
- 11 Vortex for ~ 00:00:10 10s



- 12 Sonicate in a waterbath sonicator for  00:05:00 in an ice slurry. 5m
- 13 Load each sample into a fractionation column, replace the cap and centrifuge at  3000 rcf, 00:02:00 (keep eluate as 'flow through' fraction). 2m
- 14 Place the column into a new tube, and load  300 μL of water, and centrifuge at  3000 rcf, 00:02:00 (keep eluate as 'wash' fraction). 2m
- 15 Place the column into a new tube, and load the TMT wash solution (5% ACN, 0.1% triethylamine (TEA)) .
- 16 Place the column into a new tube, and load  300 μL of the appropriate elution solution (see table below), and centrifuge at  3000 rcf, 00:02:00 to collect the fraction. 2m
- 17 Repeat step 5 for each step of the gradient fraction.
- 18 If you are concatenating the fractions, combine the fractions into the desired combinations
- 19 Lyophilise all samples until there are only a few μL left in the tube