

Oct 14, 2019

Reverse transcription of RNA to cDNA

DOI

dx.doi.org/10.17504/protocols.io.77chriw

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.77chriw>

Protocol Citation: Igem Dusseldorf 2019. Reverse transcription of RNA to cDNA. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.77chriw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 14, 2019



Last Modified: October 14, 2019

Protocol Integer ID: 28612

Keywords: pipet 100 ng rna, rna 40x yellow sample buffer, cdna 100 ng rna, μ l revertaid reverse transcriptase, ng template rna, μ l ribolock rnase, revertaid reverse transcriptase, transcription of rna, rnase, yellow sample buffer to each reaction cdna, transcription, pcr programm, rna, yellow sample buffer random hexamer, completion of rt reaction, pcr stripe, rt reaction, μ l 40x yellow sample buffer, mastermix of the following reagent, μ l mastermix to each reaction, ng for rt, pcr, reaction cdna, mm of each nucleotide, following reagent, μ l dntp mix, buffer dntp mix, rt, 5x rt, yellow sample buffer



Abstract

- 100 ng RNA
- 40x yellow sample buffer
- random hexamer primers
- 5x RT-buffer
- dNTP Mix (10 mM of each nucleotide)
- RevertAid Reverse Transcriptase
- *You will require 100 ng for RT and 100 ng for -RT control (include a -RT control for each sample!)*

Preparation of 1.33X yellow buffer

33,25 µL 40x yellow sample buffer + 966,75 µL H₂O

Start

In PCR stripes, pipet in the following order:

100 ng template RNA

1 µL random hexamer primers

RNase-free H₂O to 12,5 µL volume

Prepare a Mastermix of the following reagents:

4 µL 5x RT-buffer

0.5 µL RiboLock RNase-Inhibitor

2 µL dNTP Mix (10 mM of each nucleotide)

1 µL RevertAid Reverse Transcriptase

Add 7.5 µL Mastermix to each reaction

For your -RT control, pipet 100 ng RNA in PCR stripes, leave out buffer and everything, just add H₂O to a final volume of 20 µL. Mark as -RT so you can distinguish it from your actual cDNA!

PCR Programm

10 min 25 °C

60 min 42 °C

10 min 70 °C

After completion of RT reaction:

Add 60 µL 1,33 x yellow sample buffer to each reaction

cDNA can be stored at 4 °C for a short time, otherwise, freeze at -20 °C

Troubleshooting

