

Oct 16, 2019

Reverse transcription of RNA to cDNA

DOI

dx.doi.org/10.17504/protocols.io.8bjhskn

Igem Dusseldorf¹

¹Heinrich-Heine Universität Düsseldorf



Igem Dusseldorf

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.8bjhskn

Protocol Citation: Igem Dusseldorf 2019. Reverse transcription of RNA to cDNA. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8bjhskn>

Manuscript citation:

https://assets.thermofisher.com/TFS-Assets/LSG/manuals/MAN0012757_RevertAid_Reverse_Transcriptase_UG.pdf,
https://assets.thermofisher.com/TFS-Assets/LSG/manuals/MAN0012904_DyNAmo_ColorFlash_SYBR_Green_qPCR_F416L_UG.pdf

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 16, 2019

Last Modified: October 16, 2019

Protocol Integer ID: 28747

Abstract

Synthesizing cDNA (for qPCR)



Guidelines

Always work with gloves and safety gear and work on ice

Materials

- 100 ng RNA
- 40x yellow sample buffer (from DyNAmo Color Flash Kit)
- random hexamer primers
- 5x RT-buffer
- dNTP Mix (10 mM of each nucleotide)
- RevertAid Reverse Transcriptase
- **You will require 100 ng for RT and 100 ng for -RT control (include a -RT control for each sample!)**



Preparation of 1.33X yellow buffer

- 1 33,25 μL 40x yellow sample buffer + 966,75 μL H₂O

optional: This was used to be able to see the pipetting scheme for qPCR more easily

Start

- 2 In PCR stripes, pipet in the following order:

100 ng template RNA
1 μL random hexamer primers
RNase-free H₂O to 13 μL volume

- 3 Prepare a Mastermix of the following reagents:

4 μL 5x RT-buffer
2 μL dNTP Mix (10 mM of each nucleotide)
1 μL RevertAid Reverse Transcriptase

- 4 Add 7 μL Mastermix to each reaction

- 5 For your -RT control, pipet 100 ng RNA in PCR stripes, leave out buffer and everything, just add H₂O to a final volume of 20 μL . Mark as -RT so you can distinguish it from your actual cDNA!

PCR

- 6

PCR protocol:

10 min	25 °C
60 min	42 °C
10 min	70 °C

- 7 Add 60 μL 1,33 x yellow sample buffer to each reaction

optional: This was used to be able to see the pipetting scheme for qPCR more easily



8 cDNA can be stored at 4 °C for a short time, otherwise, freeze at -20 °C