


Nov 06, 2018

RNA extraction from adult Aiptasia

DOI

dx.doi.org/10.17504/protocols.io.sp5edq6



Kelsey Speer¹, Angela Poole², modified from Weis lab protocol by Mauricio Rodriguez-Lanetty³

¹University of Pennsylvania; ²Berry College; ³Oregon State University

Aiptasia Symbiodiniacea...



Kelsey Speer

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.sp5edq6

External link: http://people.oregonstate.edu/~weisv/assets/trizol_rneasyhybrid.pdf

Protocol Citation: Kelsey Speer, Angela Poole, modified from Weis lab protocol by Mauricio Rodriguez-Lanetty 2018. RNA extraction from adult Aiptasia. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.sp5edq6>

Manuscript citation:

Rodriguez-Lanetty, M., Phillips, W.S., Weis, V.M., 2006. Transcriptome analysis of a cnidarian-dinoflagellate mutualism reveals complex modulation of host gene expression. BMC Genomics 7, 23. doi:10.1186/1471-2164-7-23

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 17, 2018

Last Modified: February 21, 2019

Protocol Integer ID: 14813



Keywords: Aiptasia, mRNA

Abstract

A protocol for extracting RNA from Aiptasia, modified from the Weis lab protocol by Angela Poole and Mauricio Rodriguez-Lanetty. Note that the final product has significant DNA contamination and needs to be cleaned up via DNase treatment.

Materials

MATERIALS

☒ RNeasy Mini Kit **Qiagen Catalog #74104**

☒ TRIzol Reagent **Thermo Fisher Scientific Catalog #15596026**

☒ Micropestle

Safety warnings

! All Trizol steps should be performed in the hood. All Trizol waste should be treated as hazardous.

Before start

Make sure to setup a clean area in the fume hood, wipe down all pipettes with RNase Away/Zap, and cool down the centrifuge.



- 1 Homogenize starting material in Trizol using a micropestle. Trizol volume should be ~ 1 mL per 0.1g tissue. Once the homogenate is uniform, incubate at room temperature for 5 minutes.

Safety information

All Trizol steps (1-5) should be performed in a fume hood and all Trizol waste should be treated as hazardous.

- 2 Centrifuge samples for 10 minutes at 12,000 x g at 4°C. This step will pellet any debris and polysaccharides. Transfer the supernatant to a new tube and discard the pellet.
- 3 Add chloroform to the sample (0.2 mL per 1 mL of Trizol used) and shake vigorously for 20 seconds. Incubate samples at room temperature for 2-3 minutes.
- 4 Centrifuge samples for 18 minutes at 10,000 x g at 4°C.

Note

Note: past this step you should be careful of RNase contamination

- 5 Carefully remove the top aqueous (clear) layer and transfer to a new sterile, RNase-free tube.

Note

IMPORTANT: When you take your samples out of the centrifuge, you will see three layers: aqueous (clear, top), interphase (white, middle) and organic (red, bottom). Stay away from the white middle layer - that is where the DNases and RNases are. It is better to leave a little of the aqueous layer behind than to risk contamination.

- 6 **Slowly** add an equal volume of 100% RNase-free EtOH and mix gently by pipetting or inverting tube.
- 7 Load sample (up to 700 µL) into an RNeasy spin column (Qiagen) and spin at room temperature for 1 min at 8,000 x g. Discard flow through.

**Note**

If your total volume exceeds 700 μ L, repeat this step until all of the sample has passed over the column.

- 8 Add 700 μ L of Buffer RW1 to the column and spin at room temperature for 1 min at 8,000 x g. Discard flow through.
- 9 Transfer spin column to a new collection tube and add 500 μ L of Buffer RPE (check that ethanol has been added) and spin at room temperature for 1 min at 8,000 x g. Discard flow through.
- 10 Add 500 μ L of Buffer RPE and spin at room temperature for 2 min at 8,000 x g. Discard flow through.
- 11 Spin for an additional 1 minute at 8,000 x g to get rid of remaining buffer in column.
- 12 Transfer the spin column to a new, RNase-free 1.5 mL collection tube and pipette 30-50 μ L of RNase-free water over the center of the membrane (do not touch the membrane with your pipette tip!). Let the sample sit at room temperature for 1 minute and then spin for 1 minute at 8,000 x g to elute RNA.
- 13 Discard the spin column and store your RNA at -80°C until further use.