

May 14, 2025

RNA extraction from Drosophila embryos

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwq73dvmk/v1>

Artem Ilin¹

¹SU



Artem Ilin

SU

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwq73dvmk/v1>

Protocol Citation: Artem Ilin 2025. RNA extraction from Drosophila embryos. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.e6nvwq73dvmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 14, 2025



Last Modified: May 14, 2025

Protocol Integer ID: 218238

Keywords: rna extraction from drosophila embryo, extraction from drosophila embryo, drosophila embryo, rna extraction, simple protocol for rna extraction, rna, cDNA synthesis, embryo, extraction

Abstract

A simple protocol for RNA extraction and cDNA synthesis from Drosophila embryos

Materials

Invitrogen (Thermo Fischer)

TRIzol (cat number: 15596026)

Chloroform

Isopropanol

RNase-free water

Sigma-Aldrich

DNase I (AMPD1)

Applied Biosystems (Thermo Fischer)

High-Capacity RNA-to-cDNA Kit (Cat Number: 4387406)



RNA extraction

1h 3m 30s

- 1 Collect 100-200 staged embryos in a 1.5 ml eppendorf tube. Wash with PBT twice. Either centrifuge them at a low speed (≤ 500 g) for a minute or let them precipitate by themselves. 10m
- 2 Remove all the liquid, add  50 μL of Trizol. Homogenize the embryos using a plastic pestle. Depending on the amount of material, add  50-150 μL of Trizol trying to wash the pestle. The protocol proceeds with the assumption that the total volume of Trizol is now  200 μL . Incubate for  00:05:00 at room temperature. 10m
- 3 Add 1/5 volume of chloroform (for  200 μL of Trizol it would be  40 μL of chloroform). Vortex for  00:00:30 and incubate at RT for  00:03:00 3m 30s
- 4 Centrifuge for  00:15:00 at  12000 x g at  4 $^{\circ}\text{C}$ 15m
- 5 Transfer the aqueous phase to a new tube and add 0.8-1 volume of isopropanol. Incubate for  00:10:00 at RT if expected yield is high, or at  -20 $^{\circ}\text{C}$ for  02:00:00 if you want RNA to precipitate better. 10m
- 6 Centrifuge for  00:10:00 at  12000 x g at  4 $^{\circ}\text{C}$ 10m
- 7 Discard the supernatant and wash the pellet one time with  200 μL of cold 70% EtOH. Try to dislodge the pellet from the tube side but do not suck it into the pipette tip when mixing.
- 8 Centrifuge for  00:05:00 at  10000 x g at  4 $^{\circ}\text{C}$. Discard the supernatant. 5m
- 9 (Optional) Dry the pellet at  55 $^{\circ}\text{C}$ for  00:01:00
- 10 Resuspend the pellet in  10-20 μL of RNase-free water.



- 11 (Optional) Dissolve at 55 °C for 00:01:00
- 12 Measure the RNA concentration at NanoDrop spectrophotometer.

DNase treatment and reverse transcription

1h 35m

- 13 Use 1 µg of RNA from the previous step (the total volume shouldn't exceed 6.4 µL).
- 14 Add 0.8 µL of DNase buffer and 0.8 µL of DNase I from Sigma DNase I Kit to your RNA, if the volume is less than 8 µL, add RNase free water until it is equal to 8 µL. Mix by pipeting. Leave for 00:20:00 on RT. (Adding 1 µL of buffer and DNase would not result in the failure of the experiment, so one can opt to do that, if their pipette is not very precise). **20m**
- 15 Add 1 µL of Stop solution, mix by pipeting and put on 65 °C for 00:10:00. **10m**
- 16 Transfer everything to a 200 µl PCR tube. Add 10 µl of RT buffer and 1 µl of RT enzyme from High-Capacity RNA-to-cDNA Kit. Mix by pipeting and put in the PCR machine for RT reaction. Set the program like this: **1h 5m**
- 37 °C for 01:00:00
 - 95 °C for 00:05:00
- 17 Add RNase-free water to the final volume of 200 µl. Your cDNA is ready.