



Feb 14, 2025

RNA Extraction From Mouse Tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly618pgx9/v1>

Yujie Fan¹, Chelsie Steele¹

¹Caltech



Chelsie Steele

Caltech

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly618pgx9/v1>

Protocol Citation: Yujie Fan, Chelsie Steele 2025. RNA Extraction From Mouse Tissue . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.eq2ly618pgx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 14, 2025



Last Modified: February 17, 2025

Protocol Integer ID: 118323

Keywords: rna from mouse tissue, rna extraction, rna, mouse tissue, extraction, tissue, mouse, extract

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020495

Abstract

how to extract RNA from mouse tissue

Materials

RNASE-ZAP

pentobarbital sodium and phenytoin sodium

DNA LO BIND 1.5 mL tubes

Qiagen Rneasy Mini Columns Protocol (Cat No 74104, Qiagen).

DPBS

Preparation before Perfusion

- 1 Before beginning perfusion make sure all surfaces and tools are sprayed down with 70% ethanol followed by RNASE-ZAP.
- 2 Ensure perfusion pump is PFA free, DPBS is sterile filtered and on ice, and the surface used to perfuse the mouse is also sterilized with 70% ethanol followed by RNASE-ZAP
- 3 Get a bucket of dry ice to snap freeze the extracted tissue
- 4 Get two cooling blocks - should be at least 0 degrees. Put Petri dishes full of sterile DPBS on top of the cooling blocks.
- 5 Prepare a 50 ml tube with RNASE-ZAP. Put your tools in this tube while you wait to begin the perfusion.
- 6 Label DNA LO BIND 1.5 mL tubes with the appropriate label for the mouse and tissue you want to convert into RNA.

Mouse Perfusion with only DPBS

- 7 Begin the perfusion as normal - Anesthetize mouse using pentobarbital sodium and phenytoin sodium. A 30g mouse would need roughly 100 ul of this pentobarbital sodium and phenytoin sodium solution.
- 8 Ensure deep anesthesia with paw pinch. Mouse should not respond.
- 9 Using four 26g needles, secure mouse paws to styrofoam.
- 10 Make an incision along the abdomen, extending up to the rib cage.
- 11 Expose chest cavity and secure with hemostats clamped to the xiphoid process.
- 12 Insert 22g needle into the left ventricle and hold via hemostat.



- 13 Make small incision in the right ventricle.
- 14 Immediately start peristaltic pump connected to sterile filtered DPBS, flowing at a rate of 10 mL/min
- 15 Once the perfusate runs clear, stop pump

Extracting tissue

- 16 Remove the tissue needed from the animal and place into the DPBS on the cooling core. For example: if we want to extract RNA from the duodenum, the animals intestinal track would be removed from the carcass and placed into the DPBS on the cooling core
- 17 Make small precise cuts to remove roughly 30mg of tissue and place into DNA LO BIND 1.5 mL tubes.
- 18 Put DNA LO BIND 1.5 mL tubes straight into the dry ice to snap freeze the tissue.
- 19 Repeat for however many tissues you want to convert into RNA.
- 20 The tissue can be stored at -80°C until ready to convert the tissue into RNA.

Qiagen Rneasy Mini Columns Protocol

- 21 From here on we follow the Qiagen Rneasy Mini Columns Protocol (Cat No 74104, Qiagen) for converting animal tissue into RNA.