

Feb 08, 2026

RNA Extraction, Library Preparation, Bulk RNA Sequencing, and Analysis of Mouse Proximal Large Intestine

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvy7wogmk/v1>

Anastasiya Moiseyenko¹

¹California Institute of Technology

Faecalibacterium prausn...



Anastasiya Moiseyenko

California Institute of Technology

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvy7wogmk/v1>

Protocol Citation: Anastasiya Moiseyenko 2026. RNA Extraction, Library Preparation, Bulk RNA Sequencing, and Analysis of Mouse Proximal Large Intestine. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvy7wogmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 08, 2026

Last Modified: February 10, 2026

Protocol Integer ID: 242882

Keywords: RNA extraction, library preparation, bulk RNA sequencing, RNA analysis, differentially expressed genes, proximal large intestine, mouse, large intestine this protocol detail, procedure for the rna extraction, proximal large intestine, bulk rna sequencing, rna extraction, large intestine, bulk rna, genes in the mouse, analysis of mouse proximal, expressed gene, mouse proximal

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: Grant ID: ASAP-020495

Aligning Science Across Parkinson's

Grant ID: Grant ID: ASAP-000375

Heritage Medical Research Institute

Grant ID: n/a

Abstract

This protocol details the procedure for the RNA extraction, library preparation, bulk RNA sequencing, and analysis used to assess differentially expressed genes in the mouse proximal large intestine.

Materials

- TRIzol Reagent (Zymo Research)
- Direct-zol RNA MiniPrep Plus Kit (Zymo Research)
- DNase I (included with kit)
- NanoDrop spectrophotometer
- Agilent 2100 Bioanalyzer (Agilent Technologies)
- ABclonal FAST RNA-seq Library Prep Kit
- Illumina NovaSeq X sequencing platform
- STAR aligner v2.7.1
- Mus musculus GRCm36 Primary Assembly reference genome (GENCODE Release M36)
- featureCounts
- R v4.3.0
- DESeq2 R package
- clusterProfiler R package

Tissue collection

- 1 Dissect approximately 1 cm of proximal large intestine from each mouse.
- 2 Immediately place tissue into TRIzol reagent (Zymo Research) and flash-freeze on dry ice.
- 3 Store samples at -80°C until RNA extraction.

Total RNA Extraction

- 4 Thaw samples on ice.
- 5 Extract total RNA using the Direct-zol RNA MiniPrep Plus Kit (Zymo Research) according to the manufacturer's protocol.
- 5.1 Perform on-column DNase I treatment as specified by the kit to remove genomic DNA.
- 6 Measure RNA concentration and purity using a NanoDrop spectrophotometer.
- 7 Assess RNA integrity using an Agilent 2100 Bioanalyzer.

Library Preparation and Sequencing

- 8 Prepare RNA-seq libraries using the ABclonal FAST RNA-seq Library Prep Kit according to the manufacturer's instructions.
- 9 Perform sequencing using 50 bp paired-end reads on an Illumina NovaSeq X platform with a target a sequencing depth of ~30 million reads per sample.

Bioinformatics Analysis



- 10 Align raw sequencing reads to the *Mus musculus* GRCm36 Primary Assembly reference genome (GENCODE Release M36) using STAR aligner v2.7.1 with default parameters.
- 11 Quantify gene-level read counts using featureCounts. Retain only protein-coding genes with a minimum total count ≥ 10 across all samples.
- 12 Import count data into R (v4.3.0). Normalize counts and perform differential expression analysis using the DESeq2 package.
- 12.1 Identify differentially expressed genes by an adjusted p-value < 0.05 and Log2 fold change threshold as defined per comparison.
- 13 Perform gene ontology and pathway enrichment analyses using the clusterProfiler package. Using the following criteria: Adjusted p-value < 0.1 , log2 fold change > 0.2