

Nov 21, 2024

RNA Extraction using Trizol

 Journal of Virology

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v922pml3e/v1>

Carolina Lopez¹

¹Washington University



Carolina Lopez

Washington University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v922pml3e/v1>

External link: <https://doi.org/10.1128/jvi.00894-25>

Protocol Citation: Carolina Lopez 2024. RNA Extraction using Trizol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v922pml3e/v1>

Manuscript citation:

Pye SE, Achouri E, Yang Y, Musa A, López CB Validation of diverse and previously untraceable Sendai virus copyback viral genomes by direct RNA sequencing. Journal of Virology 99(8). doi: [10.1128/jvi.00894-25](https://doi.org/10.1128/jvi.00894-25)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: May 30, 2024

Last Modified: July 18, 2025

Protocol Integer ID: 100946

Keywords: trizol protocol for rna extraction, rna extraction, using trizol protocol, rna, cell line

Abstract

Protocol for RNA extraction from cell lines and tissue

Materials

	A	B	C
	Reagent	Source	Identifier
	<i>TRIzol Reagent</i>	<i>Life Technologies</i>	<i>REF#15596018</i>
	<i>Linear Acrylamide</i>	<i>Ambion</i>	<i>CAT#AM9520</i>



RNA Extraction using Trizol

1h 15m 15s

1 From cell lines:

1h 15m 15s

1. Add  1 mL of TRIZOL for RNA isolation.
2. Add  200 μ L of chloroform and shake by  00:00:15 vigorously.
3. Leave for 3-10 min at RT to allow for phase separation.
- 4*. Centrifuge for  00:15:00 at Max speed,  4 °C *
5. Take 3/4 of the upper aqueous phase (very carefully).
6. Mix with the same volume of cold 100% isopropanol as the volume taken from the aqueous phase. (500 μ L)
- Optional step:
 - ♦ Add  7 μ L of acrylamide (if you expect too low amount of RNA from the sample, typically don't need if extracting from 6 well or larger)
7. Gently invert.
- 8*. Leave for  00:10:00 at RT or longer at  -20 °C (If expecting low yield, leave O/N at  -20 °C)
9. Centrifuge for  00:15:00 at Max at  4 °C .
10. Discard the supernatant with care.
11. Wash the pellet by  1 mL of 75% (or 70%) Cold EtOH.
12. Vortex (or pipette) the sample briefly to completely wash pellet.
13. Centrifuge for  00:10:00 at max speed at  4 °C .
14. Remove the supernatant.
 - Remove supernatant by first dumping into waste container, then draining onto paper towel
 - Spin tubes down in microfuge
 - Remove residual EtOH from the pellet.
15. Dry the pellet for about 10-20min.
 - Dry upside down on kimwipe
 - Some people like to dry on the air-intake of the hood
- 16*. Resuspend the pellet with  15 μ L of dH₂O
17. Incubate at  65 °C for  00:05:00 .

*steps different when extracting RNA from **in vivo samples**

4. Spinning for  00:20:00 can help remove "white floaties" from aqueous phase
8. Do not let RNA precipitate for longer than 10 min, usually can start centrifuge step right away



16. Resuspend in  50-60 μL of dH_2O

2

From tissue samples:

* Indicate step modification when extracting RNA from **tissue**:

4. Spinning for  00:20:00 can help remove "white floaties" from aqueous phase

8. Do not let RNA precipitate for longer than 10 min, usually can start centrifuge step right away

16. Resuspend in  50-60 μL of dH_2O