

Mar 07, 2019

## RNA extraction with PGTX



Forked from [RNA extraction with PGTX](#)

DOI

[dx.doi.org/10.17504/protocols.io.ywwfxfe](https://dx.doi.org/10.17504/protocols.io.ywwfxfe)

Lutz LCB Berwanger<sup>1</sup>

<sup>1</sup>Heinrich-Heine Universität Düsseldorf

M4455 - Synthetische B...



Sebastian ST Triesch

Institute for Synthetic Microbiology, HHU Düsseldorf

---

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.ywwfxfe](https://dx.doi.org/10.17504/protocols.io.ywwfxfe)

**Protocol Citation:** Lutz LCB Berwanger 2019. RNA extraction with PGTX. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.ywwfxfe>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** March 07, 2019

**Last Modified:** March 07, 2019

**Protocol Integer ID:** 21174



## Abstract


This protocol is used for RNA extraction in cyanobacteria after Pinto et al. 2009.

**Pinto, Fernando Lopes; Thapper, Anders; Sontheim, Wolfgang; Lindblad, Peter (2009): Analysis of current and alternative phenol based RNA extraction methodologies for cyanobacteria. In: *BMC molecular biology* 10, S. 79. DOI: 10.1186/1471-2199-10-79**




## Materials

### STEP MATERIALS

 Chloroform


 Isopropanol


 Ethyl alcohol, Pure 200 proof, for molecular biology Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023

## Protocol materials


 Chloroform

 Isopropanol

 Ethyl alcohol, Pure 200 proof, for molecular biology Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023

 Chloroform

 Isopropanol

 Ethyl alcohol, Pure 200 proof, for molecular biology Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023

## Safety warnings

### ! Additional Information

Always provide RNA work on ice and always wear Gloves.



### GHS Label elements, including precautionary statements

Signal word                      Danger

Hazard statement(s)

H227                      Combustible liquid.

H302                      Harmful if swallowed.

H319                      Causes serious eye irritation.

H311                      Toxic in contact with skin.

H314                      Causes severe skin burns and eye damage.

H330                      Fatal if inhaled.

H341                      Suspected of causing genetic defects.

H371                      May cause damage to organs.

H373                      May cause damage to organs through prolonged or repeated exposure.

H402                      Harmful to aquatic life.

Precaution Statement(s)

P260                      Do not breathe dust/ fume/ gas/ mist/ vapours/ spray.

P280                      Wear protective gloves/ protective clothing/ eye protection/ face protection.

P284                      Wear respiratory protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310                      Immediately call a POISON CENTER or doctor/ physician.

### Material list

Chemical	Volume/Mass	Company	Serial no.	Comments
Phenol	39.6g			See add. Info
Glycerol	6.9mL			
8-hydroxyquinoline	0.1g			
EDTA	0.58g			See add. Info



	Sodium acetate	0.8g			
	Guanidine thiocyanate	9.5g			
	Guanidine hydrochloride	4.6g			
	Chloroform				See add. Info
	Isopropanol				See add. Info
	EtOH	70%			
	DEPC treated Water		Roth		

### Buffer

Name	Ingredients		Comments
100 mL PGTX Solution:	Phenol Glycerol 8-hydroxyquinoline EDTA Sodium acetate Guanidine thiocyanate Guanidine hydrochloride	39.6 g 6.9 mL 0.1 g 0.58 g 0.8 g 9.5 g 4.6 g	

### Equipment list

Device	Company	Comments
Centrifuge	Eppendorf	4°C
Thermo Mix	Eppendorf	95°C
Centrifuge	Heraeus	4°C
Freezer	Liebherr	-20°C



1 Fill 50 mL Falcon tube with ice. Fill with bacterial liquid culture, up to a volume of ~45 mL.

2 Centrifuge for 3 minutes at 4°C and 4.000 g.

00:03:00

3 Discard supernatant. Resuspend cell pellet in residual water (~1 mL). Transfer to 2 mL 'safe lock' tube.

4 00:00:15

5 Resuspend the cell pellet in 1 ml of PGTX. Freeze in liquid nitrogen and store at -20 °C.

#### Safety information

Wear goggles, a lab coat and gloves when dealing with PGTX and liquid nitrogen.

6 Incubate for 5 min at 95 °C, shaking at 250 rpm in Thermomixer (Eppendorf)

00:05:00

7 Rapidly chill 5 min on ice.

00:05:00

8 Add 700 µl Chloroform.

#### Safety information

Wear goggles, a lab coat and gloves when dealing with phenol and chloroform.

Chloroform

9 Let the samples incubate for 10 min at room temperature in a Thermomixer. Vortex from time to time.

00:10:00

10 Centrifuge for 15 min at 14.000 g, 4 °C. Transfer the upper aqueous phase (**contains RNA**) to a fresh reaction tube and add the same volume (450 µL) of Aqua- P/C/I (Phenol/Chloroform/Isoamylalcohol).

00:15:00

**Safety information**

Wear goggles, a lab coat and gloves when dealing with phenol and chloroform.

- 11 Thoroughly mix by vortexing. Centrifuge for 15 min at high speed. Transfer the upper aqueous phase to a 1.5 mL reaction tube.

00:15:00

**Safety information**

Wear goggles, a lab coat and gloves when dealing with phenol and chloroform.

- 12 Add 1 volume of isopropanol.

Isopropanol

- 13 Mix and incubate for at least 30 min at -20 °C. (Can be left overnight.)

00:30:00

- 14 Centrifuge at 14.000 g for 30 min.

00:30:00

- 15 Discard supernatant.

- 16 Wash with 75% chilled ethanol. Avoid resuspending the pellet.

Ethyl alcohol, Pure 200 proof, for molecular biology **Merck MilliporeSigma**  
(Sigma-Aldrich) Catalog #E7023

- 17 Centrifuge at 14.000 g.

00:05:00

- 18 Repeat washing step with 75% chilled ethanol.

**Note**

Remove excess of EtOH by using a pipet.

- 19 Air dry pellet at RT. Do not overdry!

00:10:00



20 Resuspend the pellet with 30  $\mu\text{L}$  volume of ddH<sub>2</sub>O.

Note

Usually 40  $\mu\text{L}$  DECP-treated H<sub>2</sub>O