

Feb 04, 2026

Version 1

RNA Extraction with Trizol V.1



Forked from [RNA Extraction with Trizol](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvm1zxbv3p/v1>

Dakota Betz¹

¹ucsd



Tianyuan Zhou

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.yxmvm1zxbv3p/v1>

Protocol Citation: Dakota Betz 2026. RNA Extraction with Trizol. [protocols.io](#)
<https://dx.doi.org/10.17504/protocols.io.yxmvm1zxbv3p/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: February 04, 2026

Last Modified: February 04, 2026

Protocol Integer ID: 242560

Keywords: protocol for rna extraction, rna extraction, using trizol, trizol, rna, extraction, rna extraction with trizol

Abstract

This is a protocol for RNA extractions (using Trizol). A maximum of 3 samples can be extracted at once using this protocol. The entire protocol should be run in the fumehood. This protocol is time- and temperature-sensitive, so read all instructions before attempting and make sure you are adequately prepared.

Guidelines

A maximum of 3 samples can be extracted at once using this protocol.

Materials

- RNase Away
- RNA free tubes (1.6 and 2 ul)
- Trizol
- 100% ethanol (new every time)
- RNase free water (new every time)
- Direct-zol RNA miniprep kit (Zymo research) ← contains all other reagents needed

Protocol materials

 UltraPure™ Distilled Water

 Chloroform

Safety warnings

- ❗ Trizol is BAD--conduct entire extraction protocol in the fume hood! You will need to move the centrifuge, bead basher, and other equipment into the hood before starting. Make sure you have icepacks available as this protocol is time- and temperature-sensitive.

Before start

Clean all surfaces, tools, and gloves with RNase Away. Don't forget to clean pipettes, plastic pipettes, ice blocks, centrifuge, etc.

Make sure buffers have been prepped if it's a new kit (instructions in kit for which reagents/how much to add to each).



RNA purification & HAV

1h 5m 30s

- 1
Add  200 μL of virus susoension to be extracted to a  1.5 mL tube
- 2
Lysis
Add  800 μL of  UltraPure TM Distilled Water to tube (1 per sample). Mix thoroughly by brief pulse vortexing and incubate for  00:05:00 at room temperature. 5m
- 3
Phase Separation
Add  160 μL of  Chloroform to each tube. Shake vigorously for  00:00:30 , then incubate for  00:02:00 to  00:03:00 at  Room temperature . 5m 30s
- 4
Centrifugation
Centrifuge samples  12000 rpm, 4°C, 00:15:00 15m
- 5
Aqueous Phase Collection
After centrifugation, the mixture separates into lower red, phenol chloroform phase, an interphase, and a colorless upper aqueous phase. RNA remains exclusively in the aqueous phase. Carefully remove  480 μL of the clear upper aqueous phase and transfer it to a new, labeled microcentrifuge tub  0 °C e. **Please note: Extract 80% of the aqueous RNA layer leaving 20% behind in the tube to prevent contaminating the aqueous layer.**
- 6
RNA Precipitation
Add  400 μL of  UltraPure TM Distilled Water and  1 μL of  UltraPure TM Distilled Water co-precipitant to the aqueous phase. Mix by tilting 10 times and incub  800 μL ate for  00:20:00 at  -20 °C . 20m
- 7
RNA Pelleting
Centrifuge for  12000 rpm, 4°C, 00:10:00 10m
- 8
Ethanol Wash
A dark blue RNA pellet should be visible at the bottom of the tube if GlycoBlue was used.



Carefully discard the supernatant and wash the pellet with  800 μ L of

[M] 75 % (v/v) Ethanol

9 Wash Centrifugation

Centrifuge for  12000 rpm, 4°C, 00:05:00

5m



10 Ethanol Removal

Discard the supernatant. Briefly centrifuge again to collect residual ethanol, then use a P200 pipette tip to remove the remaining ethanol without disturbing the pellet. Vacuum or air dry the RNA pellet for  00:05:00 . It is important not to let the RNA pellet dry completely as this will greatly decrease its solubility.

5m

11 1. RNA Resuspension

Resuspend the RNA pellet in  50 μ L of  UltraPure TM Distilled Water by vortexing.

The dilution step was intentionally omitted to maximize sensitivity.

12 10) Store everything at -80 °C