



Jan 31, 2025

## RNA In situ hybridization (ISH) using RNAscope®

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygwx25ov8j/v1>

Yue Sun<sup>1</sup>, Jun B. Ding<sup>1</sup>

<sup>1</sup>Stanford University



Christine Weber-Schmidt

Allen Institute

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.kxygwx25ov8j/v1>

**Protocol Citation:** Yue Sun, Jun B. Ding 2025. RNA In situ hybridization (ISH) using RNAscope®. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kxygwx25ov8j/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 23, 2025

**Last Modified:** January 31, 2025

**Protocol Integer ID:** 118932

**Keywords:** ASAPCRN, In situ hybridization, rna in situ hybridization, procedures for rna, rna, situ hybridization, assay, frozen mouse brain section, acd

**Funders Acknowledgements:**

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-020551

## Abstract

This protocol outlines the procedures for RNA in situ hybridization to simultaneously detect up to three targets using the RNAscope® Multiplex Fluorescent v2 assay (ACD). The steps adhere to the manufacturer's guidelines (Doc. No: 323100-USM) with custom modifications optimized for fixed frozen mouse brain sections.

## Attachments



Ding\_Lab\_RNA In situ...

23KB

## Materials

- o RNAscope Target Probes: C1, C2, C3 probes
- o Probe dilute (Cat. No. 300041)
- o RNAscope Multiplex Fluorescent Reagent Kit v2 (Cat. No. 323270)
- o TSA plus fluorophores: Fluorescein (NEL741001KT), Cyanine 3 (NEL744001KT), Cyanine 5 (NEL745001KT)
- o 1X Phosphate-buffered saline (PBS)
- o 4% Paraformaldehyde (PFA): Freshly prepared in 1X PBS
- o 30% Sucrose Solution: Prepared in 1X PBS.
- o 5X SSC Buffer
- o Optimal Cutting Temperature (OCT) Embedding Medium
- o Ethanol: 50%, 70%, and 100% solutions
- o Superfrost Plus microscopic slides
- o Coverslips
- o Mounting Media: Fluorescence-compatible, anti-fade mounting media
- o Steamer
- o Hybridization incubator
- o Humidity Tray
- o ImmEdge Hydrophobic Barrier Pen



## Safety warnings

! Wear appropriate PPE as required by your institution.

## Ethics statement

Prior ethics approval (e.g. IACUC) should be obtained before performing these experiments. Approval was obtained by the Stanford University IACUC before any procedures were performed.

## Sample Preparation

- 1 Perfuse mice with freshly prepared 4% paraformaldehyde (PFA) in 1X PBS.
- 2 Dissect the brain and fix in 4% PFA at  4 °C overnight.
- 3 Immerse the tissue in 30% sucrose in 1X PBS at  4 °C until the tissue sinks to the bottom of the container (approximately 48 HRS for brain tissue).
- 4 Freeze tissue in OCT embedding media with dry ice or liquid nitrogen, and store it in an airtight container at  -80 °C .
- 5 Section tissue into 20 µm thick slices and mount on Superfrost Plus slides. Store at  -80 °C or proceed to pretreatment.

## Pretreatment (Day 1)

- 6 Wash slide in 1xPBS for 5 min
- 7 Bake slides at  60 °C for 30 min in a hybridization incubator.
- 8 Post-fix slides in 4% PFA in PBS for 15 min at  4 °C .
- 9 Dehydrate tissue by immersing slides sequentially in 50%, 70%, and 100% ethanol for 5 minutes each.  
**Note:** Handle PFA in a chemical hood and collect waste from step 8-9 as biohazard.
- 10 Air dry slide for 5 min at  Room temperature (RT).





- 11 Cover sections with RNAscope® Hydrogen Peroxide and incubate for 10 minutes at Room temperature .
- 12 Wash slides twice for 5 minutes each in 1X PBS
- 13 Preheat 1xPBS and RNAscope® 1X Target Retrieval Reagent in a steamer (> 99 °C ).
- 14 Briefly immerse slides in PBS for 10 seconds, then transfer to the Target Retrieval Reagent. Cover the steamer immediately.
- 14.1 Steam slides for 3 mins.
- 15 Transfer slides to fresh PBS and rinse for 15 seconds.
- 16 Immerse slides in 100% ethanol for 3 minutes.
- 17 Air dry slides at Room temperature .
- 18 Draw barriers around sections with ImmEdge pen. Allow to dry completely.
- 19 Place slides in a humidity tray and cover sections with Protease III  
**Note:** From this step forward, keep slides in the humidity tray for all incubations.
- 20 Incubate slides at 40 °C in a hybridization incubator for exactly 20 mins.
- 21 Wash slides **twice** for 5 minutes each in PBS.  
**Note:** From this step forward, ensure sections do not dry out.



## Probe Hybridization (Day 1)



22 Dilute and mix probes of target RNA (C1, C2 and/or C3 probes) using probe dilute (1:50~2000).

**Note:** Do not mix probes from the same channel.

23 Apply probe mixture on slides to cover all sections.

24 Incubate at  40 °C for 2 hours.

25 Wash slides twice for 2 minutes each in 1x RNAscope® Wash Buffer.

26 Store slides in 5x SSC buffer overnight at  Room temperature .



## Signal Amplification (Day 2)

27 Hybridize AMP1: Apply AMP1 to cover sections. Incubate at  40 °C for 30 minutes. Wash twice for 2 minutes each in Wash Buffer.

28 Hybridize AMP2: Apply AMP2 to cover sections. Incubate at  40 °C for 30 minutes. Wash twice for 2 minutes each in Wash Buffer.

29 Hybridize AMP3: Apply AMP3 to cover sections. Incubate at  40 °C for 15 minutes. Wash twice for 2 minutes each in Wash Buffer.

## Fluorescence Signal Development (Day 2)

30 **Note:** Fluorescence signals are developed sequentially for each channel. The steps below use HRP-C1 as an example. Select the appropriate HRP channel for your probes. Channels and fluorophores (Fluorescein, Cy3 and Cy5) can be paired as desired.

31 Apply HRP-C1 to cover sections. Incubate at  40 °C for 15 minutes. Wash twice for 2 minutes each in Wash Buffer.

32 Dilute TSA fluorophores at 1:1000 (Fluorescein, Cy3 or Cy5) in TSA buffer.



- 33 Apply diluted fluorophore to cover sections. Incubate at  40 °C for 30 minutes. Wash twice for 2 minutes each in Wash Buffer.
- 34 Apply HRP Blocker to cover sections. Incubate at  40 °C for 15 minutes. Wash twice for 2 minutes each in Wash Buffer.
- 35 Repeat steps 31–34 for C2 and/or C3 probes with the appropriate HRP and fluorophores.

## Counterstaining and mounting (Day 2)

- 36 Apply DAPI to cover sections and incubate for 30 seconds at  Room temperature .
- 37 Remove DAPI by tapping or flicking the slides, and immediately add 1–2 drops of mounting media on each slide.
- 38 Place coverslips over the tissue, avoiding air bubbles.
- 39 Dry slides in dark and store at  4 °C .