


Aug 14, 2019

## RNA Isolation

DOI

[dx.doi.org/10.17504/protocols.io.6hjhb4n](https://doi.org/10.17504/protocols.io.6hjhb4n)



Natalie Brown<sup>1</sup>

<sup>1</sup>Joslin



Natalie Brown

Joslin Diabetes Center

OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.6hjhb4n](https://doi.org/10.17504/protocols.io.6hjhb4n)

**Protocol Citation:** Natalie Brown 2019. RNA Isolation. [protocols.io https://dx.doi.org/10.17504/protocols.io.6hjhb4n](https://dx.doi.org/10.17504/protocols.io.6hjhb4n)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 14, 2019

**Last Modified:** August 14, 2019

**Protocol Integer ID:** 26891

## Abstract

Real Time Polymerase Chain Reaction

## Attachments



RT PCR.docx

9KB



- 1     •     Get plates, aspirate media
  
- 2     •     Add 2mL ice cold DPBS
  
- 3     ○     Aspirate DPBS
  
- 4     •     Add 500uL Trizol/well in chemical hood
  
- 5     •     Bring yellow tubes to hood, collect trizol and add to tubes
  
- 6     •     Add 100uL chloroform to tubes
  
- 7     ○     hand shake 15s
  
- 8     ○     Place on ice 5min until see layers (aq/lipids/bottom)
  
- 9     •     Spin full speed 20min cold centrifuge
  
- 10    ○     During, make 70% ethanol on ice (3mL RNase-free water + 7mL ethanol)
  
- 11    ○     Make new ep set 1-9
  
- 12    •     In hood, collect top layer after spinning and place into eppendorfs
  
- 13    ○     \*\*take note of first volume and make all the same



- 14     ○     Mix 1:1 ethanol with supernatant, use same tip to transfer from eppendorfs into pink spin tubes
- 15     ●     Spin 9000× 1min bench centrifuge
- 16     ○     Discard flowthrough in hood, keep tube
- 17     ●     Add 350uL RW1 buffer to tubes
- 18     ○     Spin 1 min
- 19     ○     During, prepare DNase soln (with syringe)
- 20     ●     Add 70uL DNase soln to middle of each column, leave 15min RT
- 21     ●     Repeat RW1 and spin
- 22     ○     Discard collection tube and get new
- 23     ●     Wash RPE 500uL
- 24     ○     Spin 1min
- 25     ○     Discard collection and keep tubes
- 26     ●     Repeat RPE, spin 2min
- 27     ○     Discard tubes



- 28     •     Spin 1min empty to dry
  
- 29     ○     Get new eppendorfs, label
  
- 30     •     Place pink tubes on top of eppendorfs
  
- 31     •     Add 30uL RNase free water, centrifuge 1min
  
- 32     ○     Close tubes as quickly as possible and keep on ice→ nanodrop