



Jul 07, 2023

RNA isolation and qRT-PCR

DOI

<https://dx.doi.org/10.17504/protocols.io.261ge3qojl47/v1>

Narayana Yadavalli^{1,2,3,4,5}, Shawn M. Ferguson^{1,2,3,4,6,5}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Berrak Ugur

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261ge3qojl47/v1>

Protocol Citation: Narayana Yadavalli, Shawn M. Ferguson 2023. RNA isolation and qRT-PCR. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261ge3qojl47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 25, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82430

Keywords: Primers, SYBR, PCR, ASAPCRN, rna isolation from cultured cell, rna isolation, quantitative rt pcr, pcr this protocol, pcr, cultured cell, qrt, isolation

Funders Acknowledgements:

ASAP

Grant ID: 000580

Abstract

This protocol describes the RNA isolation from cultured cells and quantitative RT PCR.

Attachments



[724-1818.docx](#)

136KB

Materials

Reagents require:

1. RNeasy Micro Plus kit (Qiagen)
2. iScript cDNA synthesis Kit (Bio-Rad)
3. Primers
4. SYBR Green Master Mix (BioRad)
5. 96 well plates (BioRad)
6. Optical Adhesive Covers



RNA isolation and qRT-PCR

- 1 Aspirate media from cells and rinse cells with PBS  On ice . 
- 2 Isolate RNA using RNeasy Micro Plus kit (Qiagen) according to manufacturer's protocol.
- 3 Generate cDNA from  1 μg purified RNA using iScript cDNA synthesis Kit (Bio-Rad) according to manufacturer's protocol.
- 4 The iScript cDNA was diluted 1:10 by using sterile water.
- 5 Combine  10 μL SYBR Green Master Mix (BioRad) with  6.78 μL Sterile Water per sample.
- 6 Combine  16.78 μL diluted SYBR Green Master Mix with  0.61 μL each of  10 micromolar (μM) forward and reverse primers per sample. Pipette this mixture into wells of 96-well qPCR plate 2 technical replicates were ran for each sample.
- 7 Pipette  2 μL of diluted RNA from step 4 in well with SYBR Green Master Mix.
- 8 Cover plate with Optical Adhesive Covers (Applied Biosystems).
- 9 Spin down plate in tabletop centrifuge.
- 10 Run qPCR in CFX96 Real-Time System (BioRad) using the following PCR steps.
- 11 PCR cycle. 

temperature	time	
95	5 minutes	39 cycles
95	15 seconds	
60	30 seconds	
75	10 seconds	Melt curve
95	20 seconds	

Analysis of qRT PCR data

- 12 Calculate the ΔCT values by subtracting respective gene CT value from housekeeping gene GAPDH value.
- 13 Followed by calculating $2^{\Delta CT}$ values, normalize the mRNA expression levels for genes of interest to GAPDH and represented relative to control.