

Apr 09, 2019 Version 2

## RNA Isolation for Tissue using TRIzol V.2

DOI

[dx.doi.org/10.17504/protocols.io.zvwf67e](https://dx.doi.org/10.17504/protocols.io.zvwf67e)

Sze-Xian Lim<sup>1</sup>

<sup>1</sup>Duke University



Sze-Xian Lim

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.zvwf67e](https://dx.doi.org/10.17504/protocols.io.zvwf67e)

**Protocol Citation:** Sze-Xian Lim 2019. RNA Isolation for Tissue using TRIzol. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.zvwf67e>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development













**We are still developing and optimizing this protocol**

**Created:** April 08, 2019

**Last Modified:** April 09, 2019

**Protocol Integer ID:** 22166



- 1 Add  1 mL TRIzol of TRIzol per  30 mg tissue of tissue and homogenize using handheld homogenizer.
- 2 Incubate at  25 °C for  00:05:00 to allow nucleoprotein complexes to dissociate.
- 3 Add  200 µL chloroform of chloroform to each ml of TRIzol carefully and vortex to mix well.
- 4 Centrifuge at max speed for  00:15:00 at  4 °C .
- 5 Carefully remove the top aqueous phase and transfer to a new Eppendorf tube. The interphase and bottom organic phase can be saved for DNA and protein respectively.
- 6 Add  500 µL 100% isopropanol of 100% isopropanol to the aqueous phase, mix by inversion and incubate at  -20 °C for a minimum of  02:00:00 .
- 7 Spin down at maximum speed for 30 mins to precipitate RNA.
- 8 Remove supernatant, and add  1 mL of  75 % volume ethanol ethanol to wash pellet.
- 9 Spin down at max speed for 15 minutes and remove supernatant.
- 10 Resuspend pellet in appropriate volume of nuclease free water.

