

Apr 07, 2019 Version 2

RNA isolation for tissue V.2

DOI

dx.doi.org/10.17504/protocols.io.zvdf626

Chin Yee Tan¹

¹Duke University



Chin Yee Tan

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.zvdf626

Protocol Citation: Chin Yee Tan 2019. RNA isolation for tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.zvdf626>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development












We are still developing and optimizing this protocol

Created: April 07, 2019

Last Modified: April 07, 2019

Protocol Integer ID: 22149



- 1 Homogenization in TriZol
Add  1 mL of Trizol reagent per  30 mg of tissue and homogenize using handheld homogenizer
- 2 Incubate at RT for  00:05:00 to allow nucleoprotein complexes to dissociate
- 3 Add 1/5 the volume of Trizol ( 0.2 mL) of chloroform carefully, and vortex to mix well
- 4 Spin down at max speed in a chilled centrifuge for  00:15:00
- 5 Carefully remove the top aqueous phase and transfer to a new Eppendorf tube
The interphase and bottom organic phase can be saved for DNA and protein respectively
- 6 To the aqueous phase, add  500 μ L of 100% isopropanol, mix by inversion and incubate at  -20 °C for a minimum for  02:00:00
- 7 Spin down at max speed for  00:30:00 s to precipitate RNA
- 8 Remove supernatant, and add  1 mL of 75% Ethanol to wash the pellet
- 9 Spin down at max speed for  00:15:00 and remove supernatant
- 10 Resuspend pellet in appropriate volume of nuclease free water