

Apr 08, 2019 Version 3

## RNA isolation for tissue V.3

DOI

[dx.doi.org/10.17504/protocols.io.zvqf65w](https://dx.doi.org/10.17504/protocols.io.zvqf65w)

Chin Yee Tan<sup>1</sup>

<sup>1</sup>Duke University



Chin Yee Tan

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.zvqf65w](https://dx.doi.org/10.17504/protocols.io.zvqf65w)

**Protocol Citation:** Chin Yee Tan 2019. RNA isolation for tissue. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.zvqf65w>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working












**We use this protocol and it's working**

**Created:** April 08, 2019

**Last Modified:** April 08, 2019

**Protocol Integer ID:** 22160



- 1 Homogenization in TriZol  
Add 1ml  1 mL Trizol per  30 mg tissue and homogenize using handheld homogenizer.
- 2 Incubate at RT for  00:05:00 to allow nucleoprotein complexes to dissociate.
- 3 Add  200  $\mu$ L Chloroform carefully, and vortex to mix well.
- 4 Spin down at max speed in a chilled centrifuge for  00:15:00 .
- 5 Carefully remove the top aqueous phase and transfer to a new Eppendorf tube. The interphase and bottom organic phase can be saved for DNA and protein respectively.
- 6 To the aqueous phase, add  500  $\mu$ L of 100% isopropanol, mix by inversion and incubate at  -20 °C for a minimum for  02:00:00 .
- 7 Spin down at max speed for  00:30:00 to precipitate RNA.
- 8 Remove supernatant, and add  1 mL 75% Ethanol to wash the pellet.
- 9 Spin down at max speed for  00:15:00 and remove supernatant.



10 Resuspend pellet in appropriate volume of nuclease free water.