

Aug 15, 2019


RNA Isolation from Plant Tissue

 [GigaScience](#)

 Peer-reviewed method

DOI

dx.doi.org/10.17504/protocols.io.439gyr6



Marc T. J. Johnson¹, Eric J. Carpenter², Zhijian Tian³, Richard Bruskiewich⁴, Jason N. Burris⁵, Charlotte T. Carrigan⁶, Mark W. Chase⁷, Neil D. Clarke⁸, Sarah Covshoff⁹, Claude W. dePamphilis¹⁰, Patrick P. Edger¹¹, Falicia Goh⁸, Sean Graham¹², Stephan Greiner¹³, Julian M. Hibberd⁹, Ingrid Jordon-Thaden¹⁴, Toni M. Kutchan¹⁵, James Leebens-Mack⁶, Michael Melkonian¹⁶, Nicholas Miles¹⁴, Henrietta Myburg¹⁷, Jordan Patterson², J. Chris Pires¹¹, Paula Ralph¹⁰, Megan Rolf¹⁵, Rowan F. Sage¹⁸, Douglas Soltis¹⁹, Pamela Soltis²⁰, Dennis Stevenson²¹, C. Neal Stewart Jr⁵, Barbara Surek¹⁶, Christina J. M. Thomsen¹, Juan Carlos Villarreal²², Xiaolei Wu³, Yong Zhang³, Michael K. Deyholos², Gane Ka-Shu Wong²³

¹Department of Biology, University of Toronto at Mississauga, Mississauga, Ontario, Canada;

²Department of Biological Sciences, University of Alberta, Edmonton, Alberta, Canada;

³BGI-Shenzhen, Bei Shan Industrial Zone, Yantian District, Shenzhen, China;

⁴International Rice Research Institute, Metro Manila, Philippines;

⁵Department of Plant Sciences, University of Tennessee, Knoxville, Tennessee, United States of America;

⁶Department of Plant Biology, University of Georgia, Athens, Georgia, United States of America;

⁷Jodrell Laboratory, Royal Botanic Gardens, Kew, Richmond, Surrey, United Kingdom;

⁸Genome Institute of Singapore, Singapore, Singapore;

⁹Department of Plant Sciences, University of Cambridge, Cambridge, United Kingdom;

¹⁰Department of Biology and Intercollege Graduate Program in Plant Biology, Huck Institutes of the Life Sciences, Pennsylvania State University, University Park, Pennsylvania, United States of America;

¹¹Division of Biological Sciences, University of Missouri, Columbia, Missouri, United States of America;

¹²Department of Botany and UBC Botanical Garden, University of British Columbia, Vancouver, British Columbia, Canada;

¹³Max Planck Institute for Molecular Plant Physiology, Wissenschaftspark Golm, Am Mühlenberg 1, Potsdam-Golm, Germany;

¹⁴Department of Biology, University of Florida, Gainesville, Florida, United States of America, Florida Museum of Natural History, University of Florida, Gainesville, Florida, United States of America;

¹⁵Donald Danforth Plant Science Center, St. Louis, Missouri, United States of America;

¹⁶Department of Botany, Cologne Biocenter, University of Cologne, Cologne, Germany;

¹⁷Department of Plant Biology, North Carolina State University, Raleigh, North Carolina, United States of America;

¹⁸Department of Ecology and Evolutionary Biology, University of Toronto, Toronto, Ontario, Canada;

¹⁹Department of Biology, University of Florida, Gainesville, Florida, United States of America;

²⁰Florida Museum of Natural History, University of Florida, Gainesville, Florida, United States of America;

²¹New York Botanical Garden, Bronx, New York, United States of America;

²²Department of Ecology and Evolutionary Biology, University of Connecticut, Storrs, Connecticut, United States of America;

²³Department of Biological Sciences, University of Alberta, Edmonton, Alberta, Canada, BGI-Shenzhen, Bei Shan Industrial Zone, Yantian District, Shenzhen, China, Department of Medicine, University of Alberta, Edmonton, Alberta, Canada

GigaScience Press



Eric Carpenter

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.439gyr6](https://doi.org/10.17504/protocols.io.439gyr6)

Collection Citation: Marc T. J. Johnson, Eric J. Carpenter, Zhijian Tian, Richard Bruskiewich, Jason N. Burris, Charlotte T. Carrigan, Mark W. Chase, Neil D. Clarke, Sarah Covshoff, Claude W. dePamphilis, Patrick P. Edger, Falicia Goh, Sean Graham, Stephan Greiner, Julian M. Hibberd, Ingrid Jordon-Thaden, Toni M. Kutchan, James Leebens-Mack, Michael Melkonian, Nicholas Miles, Henrietta Myburg, Jordan Patterson, J. Chris Pires, Paula Ralph, Megan Rolf, Rowan F. Sage, Douglas Soltis, Pamela Soltis, Dennis Stevenson, C. Neal Stewart Jr, Barbara Surek, Christina J. M. Thomsen, Juan Carlos Villarreal, Xiaolei Wu, Yong Zhang, Michael K. Deyholos, Gane Ka-Shu Wong 2019. RNA Isolation from Plant Tissue. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.439gyr6>

Manuscript citation:

Appendix S1 of "Evaluating Methods for Isolating Total RNA and Predicting the Success of Sequencing Phylogenetically Diverse Plant Transcriptomes" Marc T. J. Johnson et al. PLOS ONE, November 21, 2012. <https://doi.org/10.1371/journal.pone.0050226>

License: This is an open access collection distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

These protocols were used for RNA extraction from plant tissues in order to support the One Thousand Plants initiative's work to produce RNA-Seq transcriptomes from a diverse collection of plant samples.

Created: July 03, 2019

Last Modified: August 15, 2019

Collection Integer ID: 25441

Keywords: RNA, Plant RNA, RNA Isolation, Plant Tissue

Abstract

Methods for RNA isolation

These methods were originally included in Appendix S1 of "Evaluating Methods for Isolating Total RNA and Predicting the Success of Sequencing Phylogenetically Diverse Plant Transcriptomes" Marc T. J. Johnson et al. PLOS ONE, November 21, 2012. <https://doi.org/10.1371/journal.pone.0050226>

Many of the protocols share elements or combine components from several methods. For each method, we describe the reagents and procedures used, and identify the researchers or institute that implemented the protocol.

Due to the potential for contamination and degradation by RNase enzymes, as well as health concerns in handling some substances and chemicals, best practices in aseptic wet lab techniques must be practiced at all times during RNA isolation. Chief among these are the critical need to avoid contamination of samples by using extreme care when moving liquids and opening and closing tubes to avoid aerosols. Because of the risk of degradation by RNase enzymes, it is essential to use sterile RNase-free equipment, disposable plastics and solutions. RNase degradation and contamination can be avoided by keeping samples constantly frozen at low temperature (< -80 °C) prior to adding buffers that denature or immobilize RNase. Treating equipment with RNase denaturants (e.g. RNase Zap, Ambion, Austin, TX) and solutions with diethylpyrocarbonate (DEPC) can also prevent contamination and/or degradation of samples, but it can have some negative effects on samples¹. Many additional helpful tips for successful RNA isolation are available in Sambrook and Russell¹ and in Appendix A of [Qiagen's RNeasy Mini Handbook](#) downloadable from www.qiagen.com.

¹ Sambrook, J. & Russell, D.W. Molecular Cloning: A Laboratory Manual, 3rd ed. (Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY, 2001).

Attachments



[journal.pone.0050226...](#)

297KB



Materials

Standard equipment used

Although the reagents vary among protocols, most methods used the same basic equipment, which include:

- Ceramic or porcelain mortar and pestle for tissue homogenization, or some other equipment that can homogenize frozen tissue (e.g. bead mill)
- Water bath or heating block, capable of holding temperatures up to 70°C
- Non-refrigerated microcentrifuge (with rotor for 2 ml tubes)
- Refrigerated microcentrifuge (with rotor for 2 ml tubes) – only for protocols where noted
- Liquid nitrogen (and associated thermos or dewer)
- Sterile, RNase-free disposable tips with filter barriers
- Sterile, RNase-free disposable microcentrifuge tubes (2 ml and 15 ml were most commonly used, but sizes vary among protocols)
- Pipettors – (capable of pipetting various volumes between 0.1–1000µl)
- Stainless steel spatulas
- RNase-free water
- Sterile razor blades
- Analytical balance
- Glass beakers, flasks and graduated cylinders

Safety warnings

⚠ Please see SDS (Safety Data Sheet) for hazards and safety warnings.



Attachments

 [journal.pone.005022](#)

[journal.pone.0050226...](#)
297KB

Files

 SEARCH

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 1: Qiagen RNeasy Plant Mini Kit

VERSION 1

CREATED BY



Eric Carpenter

OPEN →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 2: McKenzie et al's Qiagen hybrid method

VERSION 1

CREATED BY



Eric Carpenter

OPEN →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 3: CTAB-PVP Method

VERSION 1

CREATED BY



Eric Carpenter

OPEN →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 4: CTAB-PVP-TRIzol Method

VERSION 1

CREATED BY

Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 5: pBIOZOL Method

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 6: pBIOZOL and Qiagen RNeasy Plant Mini Kit Method

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 7: pBIOZOL-LiCl Method

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 8: CTAB/Acid Phenol/Silica Membrane Method

VERSION 1



CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 9: CTAB/Acid Phenol/Silica Membrane Method

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 10: TRIzol LS Reagent Method

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 11: Tri Reagent Method

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



NAME

RNA Isolation from Plant Tissue Protocol 12: Hot Acid Phenol Method for Angiosperms

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



#13

NAME

RNA Isolation from Plant Tissue Protocol 13: Trizol/RNAqueous Midi-Kit

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



#14

NAME

RNA Isolation from Plant Tissue Protocol 14: Ambion Trizol RNA Extraction in Microcentrifuge Tubes with Turbo DNAfree Digestion

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



#15

NAME

RNA Isolation from Plant Tissue Protocol 15: Hot Acid Phenol Method for Algae

VERSION 1

CREATED BY



Eric Carpenter

[OPEN](#) →

Protocol



#16

NAME

RNA Isolation from Plant Tissue Protocol 16: CTAB-Hot Acid Phenol Method for Algae

VERSION 1

CREATED BY



Eric Carpenter

OPEN →

Protocol



#17

NAME

RNA Isolation from Plant Tissue Protocol 17: Invitrogen PureLink-Qiagen RNeasy Hybrid

VERSION 1

CREATED BY



Eric Carpenter

OPEN →

Protocol



#18

NAME

RNA Isolation from Plant Tissue Protocol 18: innuPREP Plant RNA Kit

VERSION 1

CREATED BY



Eric Carpenter

OPEN →