

Aug 20, 2019

🌐 RNA Isolation from Plant Tissue Protocol 15: Hot Acid Phenol Method for Algae

📖 [GigaScience](#)

✓ Peer-reviewed method

📁 In 1 collection


DOI

dx.doi.org/10.17504/protocols.io.4u3gwyn

GigaScience Press



Eric Carpenter



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4u3gwyn

External link: <https://doi.org/10.1093/gigascience/giz126>

Protocol Citation: Eric Carpenter: RNA Isolation from Plant Tissue Protocol 15: Hot Acid Phenol Method for Algae. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4u3gwyn>

Manuscript citation:

Carpenter EJ, Matasci N, Ayyampalayam S, Wu S, Sun J, Yu J, Jimenez Vieira FR, Bowler C, Dorrell RG, Gitzendanner MA, Li L, Du W, K Ullrich K, Wickett NJ, Barkmann TJ, Barker MS, Leebens-Mack JH, Wong GK. Access to RNA-sequencing data from 1,173 plant species: The 1000 Plant transcriptomes initiative (1KP). Gigascience. 2019 Oct 1;8(10):giz126. doi: 10.1093/gigascience/giz126.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

These protocols were used for RNA extraction from plant tissues in order to support the One Thousand Plants initiative's work to produce RNA-Seq transcriptomes from a diverse collection of plant samples.

Created: June 27, 2019

Last Modified: August 22, 2019

Protocol Integer ID: 25211

Abstract

Implemented by: Falcia Goh and Neil Clarke

This RNA isolation method is modified from that described by Köhrer and Domdey⁵.

This protocol is part of a collection of eighteen protocols used to isolate total RNA from plant tissue. (RNA Isolation from Plant Tissue Collection: <https://www.protocols.io/view/rna-isolation-from-plant-tissue-439gyr6>)

⁵ Kohrer, K. & Domdey, H. Preparation of high molecular weight RNA. Methods in Enzymology 194, 398-405 (1991).

Attachments



[journal.pone.0050226...](#)

291KB



Materials

Reagents

Extraction Buffer:

- 1 % SDS (v/v, starting from 10 % SDS stock solution)
- 51 mM sodium acetate pH 5.5
- 10 mM EDTA
- DEPC treated water

Note

The final reaction buffer was filter purified using Nalgene 0.22 μ M filter.

Other reagents:

- Acid phenol (pH 4.3)
- Phenol:chloroform (5:1) acid equilibrated to pH 4.7 from Sigma
- Isopropanol
- 70 % ethanol (diluted in DEPC treated water H₂O)
- 3 M Sodium acetate pH 5.5












Safety warnings


⚠ Please see SDS (Safety Data Sheet) for hazards and safety warnings.



- 1 Preheat phenol and phenol:chloroform to 65 °C . Heated phenol should not be re-used.
- 2 Collect algae cells via centrifugation for 00:10:00 at 16100 x g at Room temperature .
- 2.1 Flash freeze pellets with liquid nitrogen and keep at -80 °C until extractions are carried out.
- 3 Re-suspend frozen pellet in 800 µL of preheated **extraction buffer**.
- 4 Immediately add 800 µL of hot acid phenol.
- 4.1 Vortex the tubes for 00:00:15 .
- 5 Incubate at 65 °C for 00:10:00 . Vortex every 1 min for 00:00:10 .
- 6 Centrifuge at 16100 x g at 4 °C for 00:05:00 .
- 7 The aqueous phase was transferred to fresh 1.5 ml micro-centrifuge tube.
- 8 Repeat steps 5 - 7. Repeat 3x (depending on the amount of cells used). [go to step #5](#)
- 9 Extract with equal volume of phenol:chloroform (5:1).



- 9.1 Vortex for  00:01:00 at  Room temperature .
- 9.2 Spin for  00:05:00 in a microcentrifuge at top speed.
- 9.3 Repeat step 9 three times.  [go to step #9](#)
- 10 Transfer aqueous phase to a new 1.5 ml microfuge tube. Volume should be ~  700 μL .
- 11 Add 1/10 volume of  3 Molarity (M) sodium acetate, pH 5.5, and 1 volume of isopropanol.
- 11.1 Hold at  4 °C for  00:30:00 or more.
- 12 Spin in micro-centrifuge at  4 °C at top speed for  00:20:00 .
- 13 Remove the supernatant without dislodging the pellet.
- 14 Wash the pellet with 70 % ethanol.
- 15 Invert and air dry tubes at room temperature.
- 15.1 The pellet was re-suspended in  50 μL DEPC treated H_2O .

15.2 The RNAs were stored at  -20 °C .