

Aug 19, 2019

RNA Isolation from Plant Tissue Protocol 3: CTAB-PVP Method

 [GigaScience](#)

✓ Peer-reviewed method

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.4vygw7w>

Beijing Genomics Institute¹


¹Beijing Genomics Institute

GigaScience Press

BGI



Eric Carpenter



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4vygw7w>

External link: <https://doi.org/10.1093/gigascience/giz126>

Protocol Citation: Beijing Genomics Institute 2019. RNA Isolation from Plant Tissue Protocol 3: CTAB-PVP Method. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4vygw7w>

Manuscript citation:

Carpenter EJ, Matasci N, Ayyampalayam S, Wu S, Sun J, Yu J, Jimenez Vieira FR, Bowler C, Dorrell RG, Gitzendanner MA, Li L, Du W, K Ullrich K, Wickett NJ, Barkmann TJ, Barker MS, Leebens-Mack JH, Wong GK. Access to RNA-sequencing data from 1,173 plant species: The 1000 Plant transcriptomes initiative (1KP). Gigascience. 2019 Oct 1;8(10):giz126. doi: 10.1093/gigascience/giz126.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 28, 2019

Last Modified: August 19, 2019

Protocol Integer ID: 25240

Keywords: RNA, RNA isolation, RNA extraction, plant tissue, plant RNA, rna isolation from plant tissue protocol, rna isolation from plant tissue collection, rna isolation, total rna from plant tissue, plant tissue protocol, rna, total rna, plant tissue collection, beijing genomics institute, isolation, pvp method

Abstract

Implemented by: Beijing Genomics Institute

Note

A similar protocol was used by C. dePamphilis and P. Ralph.

This protocol is part of a collection of eighteen protocols used to isolate total RNA from plant tissue. (RNA Isolation from Plant Tissue Collection: <https://www.protocols.io/view/rna-isolation-from-plant-tissue-439gyr6>)

Attachments



[journal.pone.0050226...](#)

276KB



Materials

Reagents:

CTAB-PVP Buffer:

- CTAB (Hexadecyltrimethylammonium bromide; 2 % w/v)
- PVP-40 (2 % w/v)
- 100 mM Tris-HCl (pH 8.0)
- 25 mM EDTA
- 2 M NaCl (Warmed to 65 °C in a water bath to suspend in solution)
- Add β -ME to final concentration of 2 % before use

SSTE buffer:

- 1 M NaCl
- SDS (0.5 % w/v)
- 10 mM Tris-HCl (pH 8.0)
- EDTA (1 mM)

Other reagents:

- 75 % ethanol (treated with 0.1 % DEPC)
- 96–100 % ethanol
- Acid phenol (pH 4.5)
- Chloroform
- Isoamyl alcohol
- 10 M LiCl
- Glycogen (5 mg/ml)
- 3 M Sodium Acetate (pH 5.2)
- RNase free water

Safety warnings

 Please see SDS (Safety Data Sheet) for hazards and safety warnings.



- 1 Grind tissue to a powder in liquid nitrogen.
- 2 Add  200 mg –  500 mg of ground, frozen tissue to  3.0 mL of pre-heated **extraction buffer** in a 5 ml tube.
- 3 Vortex the tube until the tissue is mixed with the buffer.
- 4 Incubate the tube at  65 °C for  00:30:00 (min), vortexing briefly ( 00:00:15) every 2–3 min during the incubation.
- 5 Aliquot the mixture into four 2 ml RNase free tubes,  1 mL each tube.
- 6 Spin the tube at  12000 x g for  00:10:00 in a microcentrifuge.

Note

All of the insoluble matter should form a pellet at the bottom of the tube.

Note

When used for algae the centrifugation was performed at  3000 x g .

- 7 Pour the supernatant into a new 2 ml RNase free tube.

**Note**

When used for algae, 15 ml tubes were used.

8 Add equal volume of 24:1 chloroform:isoamyl alcohol to fill the tube.

9 Vortex tubes until the phases mix and appear cloudy.

9.1 Incubate at 20 °C – 25 °C for 00:05:00 .

10 Spin the tubes at 12000 x g for 00:10:00 in a microcentrifuge.

Note

When used for algae the centrifugation was performed at 3000 x g .

11 Transfer the upper, aqueous phase to new 2 ml RNase free tubes.

Note

When used for algae, 15 ml tubes were used.

11.1 Repeat step 7 to 9 one more time. [go to step #7](#)

12 Transfer the upper, aqueous phase to new 2 ml RNase free tubes.



- 12.1 Add 1/3 volume of **10 M** 10 Molarity (M) LiCl to each tube.
- 12.2 Mix and let stand at **4 °C** for **06:00:00** – **08:00:00** or overnight to precipitate RNA.
- 13 Spin the tubes at **18000 x g** for **00:20:00** in a microcentrifuge.
- 13.1 Decant the supernatant, taking care not to lose the pellet.
- 14 Add **1 mL** 75 % ethanol to the pellet.
- 15 Spin the tube at maximum speed (**> 11270 x g**) for **00:05:00** in a microcentrifuge.
- 15.1 Decant the supernatant carefully.
- 16 Repeat steps 14 and 15 one more time. [↩ go to step #14](#)
- 17 Open cap and air-dry the pellet.
- 18 Add **30 µL** RNase free water to dissolve the pellet.
- 18.1 Then add **70 µL** **SSTE buffer** to each tube.



19 Combine all 4 tubes into 1 tube.

20 Add equal volume of 25:24:1 acid phenol:chloroform:isoamyl alcohol to the tube.

21 Vortex the tubes until the phases mix and appear cloudy.

21.1 Then incubate at  20 °C for  00:05:00 .

22 Spin the tube at  12000 x g for  00:10:00 .

23 Transfer the upper, aqueous phase to a new 2 ml RNase free tube.

23.1 Add equal volume of 24:1 chloroform:isoamyl alcohol to the tube.

24 Vortex the tubes until the phases mix and appear cloudy.

24.1 Then incubate at  20 °C for  00:05:00 .

25 Spin the tube at  12000 x g for  00:10:00 in a microcentrifuge.

26 Transfer the upper, aqueous phase to a new 2 ml RNase free tube.

26.1 Add 2 volumes of cooled 100 % ethanol, 1/10 volumes of NaAc (pH 5.2) and  2 µL glycogen.



26.2 Mix and incubate at -20 °C for 02:00:00 .

27 Spin the tube at 18000 x g for 00:20:00 in a microcentrifuge.

27.1 Then decant the supernatant, taking care not to lose the pellet.

28 Add 1 mL 75 % cooled ethanol to the pellet.

28.1 Leave at 20 °C for 00:03:00 .

29 Centrifuge at 4 °C for 00:05:00 at 12000 x g .

29.1 Decant the liquid carefully, taking care not to lose the pellet.

29.2 Briefly centrifuge to collect the residual liquid and remove it with a pipette.

30 Wash the pellet twice with cooled 75 % DEPC-ethanol.

30.1 Open cap and air-dry the pellet.

31 Add 30 µL RNase-free water to dissolve the pellet.

32 Treat RNA with DNase I as per supplier's protocols.