


Aug 21, 2019

RNA Isolation from Plant Tissue Protocol 6: pBIOZOL and Qiagen RNeasy Plant Mini Kit Method

GigaScience

✓ Peer-reviewed method

In 1 collection



DOI

<https://dx.doi.org/10.17504/protocols.io.4rfgv3n>

Marc T. J. Johnson¹, Eric J. Carpenter², Zhijian Tian³, Richard Bruskiewich⁴, Jason N. Burris⁵, Charlotte T. Carrigan⁶, Mark W. Chase⁷, Neil D. Clarke⁸, Sarah Covshoff⁹, Claude W. dePamphilis¹⁰, Patrick P. Edger¹¹, Falicia Goh⁸, Sean Graham¹², Stephan Greiner¹³, Julian M. Hibberd⁹, Ingrid Jordon-Thaden^{14,15}, Toni M. Kutchan¹⁶, James Leebens-Mack⁶, Michael Melkonian¹⁷, Nicholas Miles^{14,15}, Henrietta Myburg¹⁸, Jordan Patterson², J. Chris Pires¹¹, Paula Ralph¹⁰, Megan Rolf¹⁶, Rowan F. Sage¹⁹, Douglas Soltis¹⁴, Pamela Soltis¹⁵, Dennis Stevenson²⁰, C. Neal Stewart Jr⁵, Barbara Surek¹⁷, Christina J. M. Thomsen¹, Juan Carlos Villarreal²¹, Xiaolei Wu³, Yong Zhang³, Michael K. Deyholos², Gane Ka-Shu Wong^{2,3,22}

¹Department of Biology, University of Toronto at Mississauga, Mississauga, Ontario, Canada;

²Department of Biological Sciences, University of Alberta, Edmonton, Alberta, Canada;

³BGI-Shenzhen, Bei Shan Industrial Zone, Yantian District, Shenzhen, China;

⁴International Rice Research Institute, Metro Manila, Philippines;

⁵Department of Plant Sciences, University of Tennessee, Knoxville, Tennessee, United States of America;

⁶Department of Plant Biology, University of Georgia, Athens, Georgia, United States of America;

⁷Jodrell Laboratory, Royal Botanic Gardens, Kew, Richmond, Surrey, United Kingdom;

⁸Genome Institute of Singapore, Singapore, Singapore;

⁹Department of Plant Sciences, University of Cambridge, Cambridge, United Kingdom;

¹⁰Department of Biology and Intercollege Graduate Program in Plant Biology, Huck Institutes of the Life Sciences, Pennsylvania State University, University Park, Pennsylvania, United States of America;

¹¹Division of Biological Sciences, University of Missouri, Columbia, Missouri, United States of America;

¹²Department of Botany and UBC Botanical Garden, University of British Columbia, Vancouver, British Columbia, Canada;

¹³Max Planck Institute for Molecular Plant Physiology, Wissenschaftspark Golm, Am Mühlenberg 1, Potsdam-Golm, Germany;

¹⁴Department of Biology, University of Florida, Gainesville, Florida, United States of America;

¹⁵Florida Museum of Natural History, University of Florida, Gainesville, Florida, United States of America;

¹⁶Donald Danforth Plant Science Center, St. Louis, Missouri, United States of America;

¹⁷Department of Botany, Cologne Biocenter, University of Cologne, Cologne, Germany;

¹⁸Department of Plant Biology, North Carolina State University, Raleigh, North Carolina, United States of America;

¹⁹Department of Ecology and Evolutionary Biology, University of Toronto, Toronto, Ontario, Canada;

²⁰New York Botanical Garden, Bronx, New York, United States of America;

²¹Department of Ecology and Evolutionary Biology, University of Connecticut, Storrs, Connecticut, United States of America;

²²Department of Medicine, University of Alberta, Edmonton, Alberta, Canada



Eric Carpenter

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4rfgv3n>

Protocol Citation: Marc T. J. Johnson, Eric J. Carpenter, Zhijian Tian, Richard Bruskiewich, Jason N. Burris, Charlotte T. Carrigan, Mark W. Chase, Neil D. Clarke, Sarah Covshoff, Claude W. dePamphilis, Patrick P. Edger, Falicia Goh, Sean Graham, Stephan Greiner, Julian M. Hibberd, Ingrid Jordon-Thaden, Toni M. Kutchan, James Leebens-Mack, Michael Melkonian, Nicholas Miles, Henrietta Myburg, Jordan Patterson, J. Chris Pires, Paula Ralph, Megan Rolf, Rowan F. Sage, Douglas Soltis, Pamela Soltis, Dennis Stevenson, C. Neal Stewart Jr, Barbara Surek, Christina J. M. Thomsen, Juan Carlos Villarreal, Xiaolei Wu, Yong Zhang, Michael K. Deyholos, Gane Ka-Shu Wong 2019. RNA Isolation from Plant Tissue Protocol 6: pBIOZOL and Qiagen RNeasy Plant Mini Kit Method. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4rfgv3n>

Manuscript citation:

Appendix S1 of "Evaluating Methods for Isolating Total RNA and Predicting the Success of Sequencing Phylogenetically Diverse Plant Transcriptomes" Marc T. J. Johnson et al. PLOS ONE, November 21, 2012. <https://doi.org/10.1371/journal.pone.0050226>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

These protocols were used for RNA extraction from plant tissues in order to support the One Thousand Plants initiative's work to produce RNA-Seq transcriptomes from a diverse collection of plant samples.

Created: June 26, 2019

Last Modified: August 21, 2019

Protocol Integer ID: 25095

Keywords: RNA, RNA isolation, RNA extraction, plant tissue, rna isolation from plant tissue protocol, rna isolation from plant tissue collection, total rna from plant tissue, rna isolation, diverse plant transcriptome, isolating total rna, qiagen rneasy plant mini kit method this protocol, qiagen rneasy plant mini kit method, plant tissue protocol, total rna, rneasy plant mini kit, plant tissue collection, rna, pbiozol, isolation

Abstract

This protocol is part of a collection of eighteen protocols used to isolate total RNA from plant tissue. (RNA Isolation from Plant Tissue Collection: <https://www.protocols.io/view/rna-isolation-from-plant-tissue-439gyr6>) and was originally published as part of Appendix S1 of "Evaluating Methods for Isolating Total RNA and Predicting the Success of Sequencing Phylogenetically Diverse Plant Transcriptomes" Marc T. J. Johnson et al. PLOS ONE, November 21, 2012. <https://doi.org/10.1371/journal.pone.0050226>

Attachments



[journal.pone.0050226...](#)

283KB

Materials

MATERIALS

 RNeasy Plant Mini Kit **Qiagen Catalog #74904**

Safety warnings

 Please see SDS (Safety Data Sheet) for hazards and safety warnings.



- 1 Grind tissue to a powder in liquid nitrogen.
- 2 Add  1.3 mL of cold ( 4 °C) pBIOZOL reagent for up to  100 mg of frozen ground tissue.
 - 2.1 Mix by briefly vortexing or flicking the bottom of the tube until the sample is thoroughly re-suspended.
- 3 Incubate the tube for  00:05:00 at  Room temperature .

Note

Lay the tube down horizontally to maximize surface area during RNA extraction.

- 4 Centrifuge for  00:10:00 at  12000 x g in a microcentrifuge at  Room temperature .
 - 4.1 Transfer the supernatant to a new  1.5 mL RNase-free tube.
- 5 Add  100 µL of  5 Molarity (M) NaCl and  300 µL chloroform.
 - 5.1 Vortex vigorously.
- 6 Centrifuge at  12000 x g for  00:10:00 .



- 7 Transfer the top aqueous phase to a new 1.5 ml RNase-free tube.
- 7.1 Add an equal volume of 5:1 acid phenol:chloroform to the tube.
- 8 Vortex the tube until the phases mix and appear cloudy.
- 8.1 Incubate at  20 °C for  00:05:00 .
- 9 Centrifuge at  12000 x g for  00:10:00 .
- 10 Transfer the top aqueous phase to a new 1.5 ml RNase-free tube.
- 10.1 Add to the aqueous phase equal volume of 24:1 chloroform:isoamyl alcohol.
- 10.2 Vortex the tube until the phases mix and appear cloudy.
- 10.3 Then incubate at  Room temperature for  00:05:00 .
- 11 Centrifuge at  12000 x g for  00:10:00 .
- 12 Transfer the top aqueous phase to a new 1.5 ml RNase-free tube.
- 12.1 Add 1/2 volume of 100 % ethanol.



- 13 Pour the contents of the tube into a Qiagen mini RNA spin column (pink), until the column is almost filled with liquid.
- 14 Cap the tube.

14.1 Centrifuge at  12000 x g for  00:00:15 .

Note

The column should be empty at the end of this spin.

- 15 Discard the flow-through from the collection tube.
- 16 Repeat the previous two steps with the same mini RNA spin column, until all of the liquid in the tube(s) has been passed through the column.

Note

The nucleic acid is now bound to the silica membrane in the spin column.

17 Apply  700 µL of solution RW1 to the spin column.

18 Cap the tube.

18.1 Centrifuge at  12000 x g for  00:00:15 .

Note

The column should be empty at the end of this spin.



19 Discard the flow-through from the collection tube.

20 Apply  500 µL of solution RPE to the spin column.

20.1 Cap the tube

20.2 Centrifuge at  12000 x g for  00:00:15 .

Note

The column should be empty at the end of this spin.

21 Discard the flow-through.

22 Repeat previous two steps one time.  [go to step #20.2](#)

23 Spin at maximum speed for  00:02:00 to remove remaining liquid from the silica membrane.

24 Transfer the spin column to a new 1.5 ml conical bottom microcentrifuge tube.

25 Add  30 µL –  50 µL of RNase-free water to the column.

25.1 Then let tube incubate at  20 °C for  00:03:00 .



26 Spin at maximum speed for  00:01:00 to collect RNA solution.