

Aug 21, 2019 Version 1

🌐 RNA Isolation from Plant Tissue Protocol 9: CTAB/Acid Phenol/Silica Membrane Method V.1

📖 [GigaScience](#)

✓ Peer-reviewed method

📁 In 1 collection


DOI

dx.doi.org/10.17504/protocols.io.4yegxte

GigaScience Press



Eric Carpenter



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4yegxte

External link: <https://doi.org/10.1093/gigascience/giz126>

Protocol Citation: Eric Carpenter: RNA Isolation from Plant Tissue Protocol 9: CTAB/Acid Phenol/Silica Membrane Method. [protocols.io](https://dx.doi.org/10.17504/protocols.io.4yegxte) <https://dx.doi.org/10.17504/protocols.io.4yegxte>

Manuscript citation:

Carpenter EJ, Matasci N, Ayyampalayam S, Wu S, Sun J, Yu J, Jimenez Vieira FR, Bowler C, Dorrell RG, Gitzendanner MA, Li L, Du W, K Ullrich K, Wickett NJ, Barkmann TJ, Barker MS, Leebens-Mack JH, Wong GK. Access to RNA-sequencing data from 1,173 plant species: The 1000 Plant transcriptomes initiative (1KP). Gigascience. 2019 Oct 1;8(10):giz126. doi: 10.1093/gigascience/giz126.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

These protocols were used for RNA extraction from plant tissues in order to support the One Thousand Plants initiative's work to produce RNA-Seq transcriptomes from a diverse collection of plant samples.

Created: June 30, 2019

Last Modified: August 23, 2019

Protocol Integer ID: 25318

Keywords: RNA, RNA Extraction, RNA Isolation, plant tissue

Abstract

Implemented by: Henrietta Myburg and Marc Johnson

This protocol is a modification of protocol 8. It was developed after protocol 8 and several commercially available plant RNA isolation kits failed to produce a sufficient yield and quality of RNA from *Oenothera* spp. (Onagraceae) for next-generation sequencing. *Oenothera* are rich in polysaccharides, oils, flavonoids and complex ellagitannins that likely interfere with isolation.

We suspect that this protocol will be most useful for species and tissues with complex secondary chemistry and rich in oils (e.g. some Rosaceae and Pinaceae).

The most important modifications to this protocol versus Protocol 8 is the use of less plant tissue, more extraction buffer, and repeating solvent extractions until the interphase is clean of debris. The protocol is regrettably longer and more involved than Protocol 8. We attempted to remove or reduce the replication of the solvent extractions steps without success (i.e. yield and quality are always decreased in *Oenothera* when any steps are removed).

This protocol is part of a collection of eighteen protocols used to isolate total RNA from plant tissue. (RNA Isolation from Plant Tissue Collection: <https://www.protocols.io/view/rna-isolation-from-plant-tissue-439gyr6>)

Attachments




[journal.pone.0050226...](#)

276KB

Materials

MATERIALS



 RNaseZap™ RNase Decontamination Solution Thermo Fisher Scientific Catalog #AM9780

Safety warnings

 Please see SDS (Safety Data Sheet) for hazards and safety warnings.





- 1 Chill mortar and pestle with liquid nitrogen (fill the mortar carefully with liquid nitrogen, let it evaporate, repeat one more time. Keep the pestle in the mortar when doing this).
- 2 Fill the mortar again.

2.1 Add  50 mg –  100 mg of tissue into the mortar while the liquid nitrogen is still “bubbling”.

- 3 Carefully grind tissue to a powder in liquid Nitrogen.

Note

Add more nitrogen if needed, but it is important to be careful when you pour the liquid nitrogen. Especially if you are grinding a small amount of tissue in a small sized mortar. The liquid nitrogen tends to “blow” the tissue out of the mortar.

- 4 Add  2 mL to  5 mL (depending on the amount of tissue you start with) of prepared extraction buffer on top of the tissue powder and mix in with the pestle.

Note

The extraction buffer will freeze because the mortar is still very cold from the liquid NITROGEN.
Do not worry about this – start grinding your next sample.
The tissue/CTAB mix will slowly thaw and at some point you will be able to mix it well using the pestle.

Note

Keep in mind that if you use more extraction buffer you will be working in more 2 ml tubes per sample and that the number of tubes per sample will increase once you reach the Qiagen RNeasy part of the protocol. Thus, there will be a lot of pipetting into the Qiagen spin column. One can do the solvent extraction steps in 15 ml Falcon tubes, but keep in mind that you will still have to pass a significantly larger amount of liquid through the Qiagen spin column.



- 5 Once thawed and well mixed, pipette or pour the slurry into a 2 ml tube.


Note




If the slurry is too thick, cut the tip of the 1 ml pipette tip with a clean razor blade or scissors.


If you start with a lot of tissue, you might have to rinse the side of the mortar with an additional ml of extraction buffer (or more if needed.)

- 6 Vortex each tube until the tissue is mixed with the buffer.

Note

To facilitate mixing, you may have to invert the tube on the vortex, and/or heat it briefly in a  65 °C water bath.

- 7 Incubate tubes at  65 °C for  00:10:00 –  00:15:00 , vortexing briefly (15 seconds) twice during the incubation.

- 8 Spin the tubes at maximum speed ( 14000 x g) for  00:03:00 in a microcentrifuge.

Note

All of the insoluble matter should form a pellet at the bottom of the tube.

- 9 Pour or pipette the supernatant into a new 2 ml tube.

**Note**


Make sure that the total volume of the supernatant per tube is ≤ 1 ml. (You will be adding approximately 1:1 volume ratios in the following steps)



Solvent Extraction #1

- 10 Add an equal volume (compared to sample) of 24:1 chloroform:isoamyl alcohol to each tube.

Note

We attempted to remove solvent extraction #1 and #2, but doing so lowers the quality and quantity of total RNA in evening primrose (*Oenothera*).

- 11 Vortex the tubes for  00:00:15 or until the phases mix and appear cloudy.


- 12 Spin the tubes at maximum speed ( 14000 x g) for  00:03:00 in a microcentrifuge.
Most of the chlorophyll will be dissolved in the lower, organic phase.

- 13 Pipette the upper, aqueous phase to a new 2 ml tube.

Note

Avoid transferring any of the material (usually a white precipitate) from the boundary between the phases.

Solvent Extraction #2

- 14 Repeat steps 10 to 13 one more time and proceed to step 15.  [go to step #10](#)

Solvent Extraction #3



15 Add equal volume (compared to the sample) of 5:1 phenol:chloroform (pH 4.5) to the tubes containing the aqueous phase.


16 Vortex each tube for  00:00:15 or until the phases mix and appear cloudy.

17 Spin the tube at maximum speed ($> 14000 \times g$) for  00:03:00 in a microcentrifuge.

18 Pipette the upper aqueous phase to a new 2 ml tube, using a disposable pipette.

Note

Avoid transferring any of the material (usually a white precipitate) from the boundary between the phases.

19 Repeat steps 15 to 18 until the interphase is clean (2-3x for most *Oenothera*) and then proceed to step 20.  [go to step #15](#)

Note

This is likely the most critical step. Do not proceed to 20 until the interphase appears clean.

Solvent Extraction #4

20 Add an equal volume (compared to sample) of 24:1 chloroform:isoamyl alcohol to each tube containing the aqueous phase.

21 Vortex each tube for  00:00:15 or until the phases mix and appears cloudy.

22 Spin the tubes at maximum speed ($> 14000 \times g$) for  00:03:00 in a microcentrifuge.



23 Pipette the upper, aqueous phase to a new 2 ml tube.

Note

Avoid transferring any of the material (usually a white precipitate) from the boundary between the phases.

Note

At this point the interphase should be relatively clean, especially if you repeated the phenol/chloroform steps to the point where the interphase appeared to be clean.

24 Repeat steps 20 to 23 one more time and then proceed to step 25. [➡ go to step #20](#)

25 Estimate the volume of the aqueous phase based on the markings on the tube.

25.1 Add at least 0.5 volumes of buffer RLT and mix well by briefly shaking.

Note

We attempted excluding this step but doing so decreased the RNA yield and quality

26 Estimate the new total volume in the tube.


26.1 Add 0.5 volumes of 95–100 % ethanol.

26.2 Mix by briefly shaking.

Note

If you hold the tube against the light while mixing, you sometimes see the RNA swirl while it precipitates out.



27 Pipette the contents of the tube into a Qiagen miniRNA spin column (pink), until the column is almost filled with liquid (\leq  750 μL).

28 Cap the tube and spin for  00:00:15 at $>$  8000 x g .

Note


The column should be empty at the end of this spin.

29 Discard the flow-through from the collection tube.

30 Repeat the previous two steps with the same spin column, until all of the liquid in all the tube(s) from the same sample have passed through the column.

Note

The nucleic acid is now bound to the silica membrane in the spin column.


31 Apply  350 μL of buffer RW1 to the spin column.

32 Cap the tube and spin for  00:00:15 at $>$  8000 x g .

Note

The column should be empty at the end of this spin.

33 Discard the flow-through from the collection tube.

34 Apply  80 μL of DNase digestion solution to the membrane of the spin column.



35 Incubate at Room temperature for 00:15:00 .

36 Apply 350 μL of solution RW1 to the spin column.

37 Cap the tube and spin for 00:00:15 at > 6300 x g .

Note

The column should be empty at the end of this spin.

38 Discard the flow-through from the collection tube.

39 Apply 500 μL of solution RPE to the spin column.

40 Cap the tube and spin for 00:00:15 at > 6300 x g .

Note

The column should be empty at the end of this spin.

41 Discard the flow-through from the collection tube.

42 Apply 500 μL of solution RPE to the spin column.


43 Cap the tube and spin for 00:00:15 at > 6300 x g .

**Note**


The column should be empty at the end of this spin.



44 Discard the flow-through from the collection tube.

45 Transfer the spin column to a new collection tube.

45.1 Spin at maximum speed for  00:03:00 to remove remaining liquid from the silica membrane.

46 Transfer the spin column to a new 1.5 ml conical bottom microcentrifuge tube

47 Add  44 μL of RNase-free water onto the column.

47.1 Incubate for  00:01:00 to elute at  Room temperature and spin to elute.