

Jan 30, 2023

RNA isolation protocol from unwilling plant tissue

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8j6n7g2w/v1>

sa.parra¹

¹Universidad de Chile



Samuel Parra

Universidad de Chile

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8j6n7g2w/v1>

Protocol Citation: sa.parra 2023. RNA isolation protocol from unwilling plant tissue. protocols.io
<https://dx.doi.org/10.17504/protocols.io.5jyl8j6n7g2w/v1>

Manuscript citation:

Modified from Chang, S., Puryear, J., & Cairney, J. (1993). A simple and efficient method for isolating RNA from pine trees. Plant Molecular Biology Reporter, 11(2), 113–116. doi:10.1007/bf02670468

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: January 30, 2023

Last Modified: January 30, 2023

Protocol Integer ID: 76076

Keywords: rna isolation protocol from unwilling plant tissue, rna isolation protocol, high quality rna extraction from plan tissue, high quality rna extraction, rna, qiagen rneasy plant kit, unwilling plant tissue, extraction, isolation, plant

Abstract

This protocol is intended for HIGH quality RNA extraction from plan tissues

We used this protocol after several attempts with Trizol, Qiagen RNeasy plant kit and other protocols with no results.

This procedure redereed us between 300 to 1000 ng with good RNA

Materials

Extraction buffer:

2% CTAB (hexadecyltrimethylammonium bromide)

2% PVP (polyvinylpyrrolidinone K 30)

100 mM Tris-HCl (pH 8.0)

25 mM EDTA

2.0 M NaCl

0.5g/L spermidine

(mix and autoclave) Alt. microweave heat

2% [B-mercaptoethanol] (add just before)

Ethanol 100%

LiCl2 10M

Chloroform:IAA 24:1

1.5 mL tubes

2 mL tubes

SSTE buffer:

1M NaCl

0.5% SDS

10mM TrisHCl (pH 8)

1mM EDTA (pH 8)



RNA isolation protocol

- 1 This protocol is an adaptation of <https://doi.org/10.1007/BF02670468> designed for extraction of RNA from pine needles with high polyphenolic compounds or high sugar content.

For all buffers use DEPC treated water or Nuclease free ddH₂O

- 2 **Extraction buffer:**
2% CTAB (hexadecyltrimethylammonium bromide)
2% PVP (polyvinylpyrrolidinone K 30)
100 mM Tris-HCl (pH 8.0)
25 mM EDTA
2.0 M NaCl
0.5g/L spermidine
(mix and autoclave) *Alt. microwave heat*
2% [B-mercaptoethanol] (add just before)

Ethanol 100%

LiCl2 10M

Chloroform:IAA 24:1

1.5 mL tubes

2 mL tubes

SSTE buffer:

1M NaCl

0.5% SDS

10mM TrisHCl (pH 8)

1mM EDTA (pH 8)

Procedure, Extraction

- 3 For  200 mg sample use  3 mL of extraction buffer.

10m

Grind the sample with liquid nitrogen until to a fine powder. Add the 3 ml of Extraction buffer, the frozen paste will become softer in time, use the cold paste to finish grinding any remanent plant tissue.

Separate the mixture in 3 separated 2mL tubes

Keep the tubes on ice



On ice

- 4 Add an equal volumen of ChCl:IAA and vortex x 30 sec

10m

Room temperature

Extract the WHOLE aqueose phase, do not leave any of it from each tube.

ignore any coloration of the aqueus phase

**Optional : Gycogen can be added at this point*

Transfer the phase to 1.5 ml tubes and add 1/4 volume of 10M LiCl and mix by inverting

- 5 Precipitate overnight at 4°C 4 °C

12h

Phases can be formed due to cold SDS in the buffers, **ignore** every phase or precipitation and go on.

RNA cleaning

2h 30m

- 6 13000 rpm, Room temperature, 00:10:00

10m

Centrifuge 10m at 13000 at RT, discard the liquid, mind the pellet!

- 7 **If the SSTE buffer is below 25°C SDS on it will form flecks, warm until sds is completely dissolved*

Dissolve the pellet with 500 µL of SSTE buffer.

Join the pellet from the 3 initial tubes using the same 500 ul of SSTE.

- 8 Extract with an equal volumen of CHCl:IAA by vortexing for 30 sec

10m

13000 rpm, Room temperature, 00:10:00 Centrifuge at 13000 rpm x 10m at RT

- 9 Recover the upper phase completely, do not leave ANY traces, a *white smudge between the phases must be retrieved. dont worry if SOME chcl2 is tranfered.*

2h

Add 2 volumes of cold ETOH100% and precipitate for at least 02:00:00 at



-20 °C

- 10 After precipitation centrifuge at 13000 rpm, 4°C, 00:10:00

10m

Discard the ethanol, DRY the rna pellet (might be invisible!)

Resuspend the pellet on 30-50 µL of ddH2O nuclease free or TE buffer