



May 04, 2020

Version 1

RNA Purification from Buccal Swabs, Nasopharyngeal Samples (swab or aspirate) and Saliva using the Monarch Total RNA Miniprep Kit V.1



DOI

dx.doi.org/10.17504/protocols.io.be93jh8n

New England Biolabs¹

¹New England Biolabs

New England Biolabs (N...

Coronavirus Method De...



Isabel Gautreau

New England Biolabs

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.be93jh8n>

External link: <https://www.neb.com/protocols/2020/03/11/rna-purification-from-buccal-swabs-using-the-monarch-total-rna-miniprep-kit-neb-t2010>



Protocol Citation: New England Biolabs 2020. RNA Purification from Buccal Swabs, Nasopharyngeal Samples (swab or aspirate) and Saliva using the Monarch Total RNA Miniprep Kit. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.be93jh8n>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: April 18, 2020

Last Modified: May 14, 2021

Protocol Integer ID: 35867

Keywords: buccal swab, saliva swab, nasopharyngeal, Viral RNA, RNA Extraction, Saliva RNA, Buccal Swab, Nasopharyngeal swab, viral rna isolation, viral RNA extraction, rna purification from buccal swab, monarch total rna miniprep kit this protocol, rna from buccal swab, monarch total rna miniprep kit, rna purification, rna, purification, nasopharyngeal sample, swab, saliva

Abstract

This protocol utilizes the Monarch Total RNA Miniprep Kit to purify RNA from buccal swabs, nasopharyngeal samples, and saliva.

Guidelines

This protocol is to be used for research use only.

Materials

MATERIALS

☒ Nuclease-free Water

☒ Microcentrifuge

☒ Monarch Total RNA Miniprep Kit **New England Biolabs Catalog #T2010S**

☒ RNase-free Microfuge Tubes (0.5 mL) **Thermo Fisher Catalog #AM12300**

☒ RNase-free Microfuge Tubes (1.5 mL) **Thermo Fisher Catalog #AM12400**

Additional Materials:

- isopropanol
- ≥95% ethanol
- 2X Monarch DNA/RNA Protection Reagent
- collection tubes (additional)
- Monarch Proteinase K
- Monarch Lysis Buffer
- 1X Monarch DNA/RNA Protection Reagent



Troubleshooting

Safety warnings

⚠ Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Before start

- Monarch DNA/RNA Protection Reagent is supplied as a 2X concentrate. Dilute with nuclease-free water only as needed, as some sample types require resuspension in the 2X concentrate, while others require a 1X solution. If purifying samples stored in Monarch DNA/RNA Protection Reagent, please review the related guidance.
- For the 50 prep kit, add 275 µl nuclease-free water to the lyophilized DNase I vial and resuspend by gentle inversion. We suggest making aliquots of DNase I, sized to your processing needs, and storing at -20°C to minimize freeze-thaw cycles (3 F/T cycles maximum)
- For the 50 prep kit, add 1,040 µl Proteinase K Resuspension Buffer to the lyophilized Proteinase K (Prot K) vial and vortex to resuspend. Store at -20°C.
- For the 50 prep kit, add 100 ml ethanol ≥ 95% (not included) to the 25 ml RNA Wash Buffer concentrate and store at room temperature
- Addition of RNA Lysis Buffer and all subsequent steps should be performed at room temperature (this will prevent precipitation of detergent in the lysis buffer). If samples are accidentally placed on ice and precipitate forms, allow the samples to return to room temperature to resolubilize before loading onto the column.











1 Sample Disruption and Homogenization

STEP CASE

Buccal + nasopharyngeal with transport medium

1 step

Buccal swabs and nasopharyngeal samples with transport medium


- 2 Place swab into a tube containing  300 µL of 1X Monarch DNA/RNA Protection Reagent to an aliquot of transport medium and vortex briefly.
- 3 For every  300 µL of DNA/RNA Protection Reagent/Sample Mixture , add  15 µL Monarch Proteinase K.
- 4 Vortex briefly and incubate at  Room temperature for  00:30:00 . 
- 5 Vortex sample briefly and spin for  00:02:00 (16,000 x g) to pellet debris. 
- 6 Transfer supernatant to an RNase-free microfuge tube.
- 7 Add an equal volume of **Monarch RNA Lysis Buffer** and vortex briefly.
- 8 Proceed to **Step 1** of **Part 2: RNA Binding and Elution**

Part 2: RNA Binding and Elution

9

Note

All centrifugation steps should be carried out at 16,000 x g.

Transfer up to  800 μ L of the sample from Part 1 to a gDNA Removal Column (light blue) fitted with a collection tube.

Note

For sample identification, label collection tubes, as gDNA removal columns will be discarded after spinning.

- 10 Spin for  00:00:30 to remove most of the gDNA.

Note

SAVE THE FLOW-THROUGH (RNA partitions here).

- 11 Discard the gDNA Removal Column.

- 12 Add an equal volume of ethanol ($\geq 95\%$) to the flow-through and mix thoroughly by pipetting.

Note

To exclude RNA ≤ 200 nt, add only 1/2 volume ethanol to flow-through. The addition of ethanol creates favorable conditions for RNA to bind to the RNA Purification column.

- 13 Transfer mixture to an RNA Purification Column (dark blue) fitted with a collection tube.



- 14 Spin for  00:00:30 .

- 15 Discard flow-through.

Note

If further gDNA removal is essential for downstream applications, proceed to on-column DNase I treatment, Step 10.1–10.3 (recommended). If not, proceed to Step 5.





- 15.1 Add  500 μ L RNA Wash Buffer and spin for  00:00:30 and discard flow-through.



Note



This ensures all salts are removed prior to the addition of DNase I.

Note

If using a vacuum manifold, add 500 μ l of RNA Wash Buffer and switch the vacuum on. Allow the solution to pass through the column, then switch the vacuum source off.

- 15.2 In an RNase-free microfuge tube (not included), combine  5 μ L DNase I with  75 μ L DNase I Reaction Buffer and pipet mixture directly to the top of the matrix.



- 15.3 Incubate for  00:15:00 at  Room temperature .

- 16 Add  500 μ L RNA Priming Buffer and spin for  00:00:30 .

- 17 Discard flow-through.

Note



If using a vacuum manifold, add 500 μ l of RNA Priming Buffer and switch the vacuum on. Allow the solution to pass through the column, then switch the vacuum source off.

- 18 Add  500 μ L RNA Wash Buffer and spin for  00:00:30 .

- 19 Discard flow-through.

**Note**

If using a vacuum manifold, add 500 µl of RNA Priming Buffer and switch the vacuum on. Allow the solution to pass through the column, then switch the vacuum source off.

- 20 Add another  500 µL RNA Wash Buffer and spin for  00:02:00 .






- 21 Transfer column to an RNase-free microfuge tube.

Note

Use care to ensure the tip of the column does not contact the flow-through. If in doubt, re-spin for 1 minute to ensure no ethanol is carried over.

Note

If using a vacuum manifold, add 500 µl of RNA Wash Buffer and switch the vacuum on. Allow the solution to pass through the column, then switch the vacuum source off.


- 22 Add  30 µL to  100 µL Nuclease-free Water directly to the center of the column matrix and spin for  00:00:30 .

**Note**

For best results, elute with at least 50 µl, which is the minimum volume needed to wet the membrane. Lower volumes can be used but will result in lower recovery (elution in 30 µl results in > 80% recovery and 100 µl provides maximum recovery). For spectrophotometric analysis of eluted RNA, it may be necessary to re-spin eluted samples and pipet aliquot from top of the liquid to ensure that the A 260/230 is unaffected by possible elution of silica particles.

- 23 Place RNA on ice if being used for downstream steps at:

 -20 °C short-term storage (less than one week)

 -80 °C long-term storage .



Note

Addition of EDTA to 0.1–1.0 mM may reduce the activity of any contaminating RNases.