

Nov 01, 2024

## RNA-seq for human iPSC-CM following GSK3 inhibition

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg32rb1v25/v1>

Taylor Anglen<sup>1</sup>, Irene M. Kaplow<sup>1</sup>, Baekgyu Choi<sup>2</sup>, Kevin Hagy<sup>1</sup>, Duc Tran<sup>1</sup>, Magan E. Ramaker<sup>1</sup>, Alejandro Barrera<sup>1</sup>, Svati Shah<sup>1</sup>, Inkyung Jung<sup>2</sup>, Ravi Karra<sup>1</sup>, Yarui Diao<sup>1</sup>, Charles A. Gersbach<sup>1</sup>

<sup>1</sup>Duke University; <sup>2</sup>Korea Advanced Institute of Science and Technology

Gersbach Lab



Taylor Anglen

Duke University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg32rb1v25/v1>

**Protocol Citation:** Taylor Anglen, Irene M. Kaplow, Baekgyu Choi, Kevin Hagy, Duc Tran, Magan E. Ramaker, Alejandro Barrera, Svati Shah, Inkyung Jung, Ravi Karra, Yarui Diao, Charles A. Gersbach 2024. RNA-seq for human iPSC-CM following GSK3 inhibition. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg32rb1v25/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 19, 2024

**Last Modified:** November 01, 2024

**Protocol Integer ID:** 110344

**Keywords:** gsk3 inhibition this protocol, gsk3 inhibition, following gsk3 inhibition, seq for human ipsc, seq methods in human ipsc, rna, human ipsc, seq, using chir99021

**Funders Acknowledgements:**

NIH

Grant ID: HG012053

## Abstract

This protocol describes RNA-seq methods in human iPSC-CM following growth with or without GSK3 inhibition using CHIR99021.

## Troubleshooting

## Matrigel coating tissue culture plates

- 1 Prepare Matrigel-coated vessels as needed, following the manufacturer's instructions. Use coated vessels within two weeks.
- 2 Thaw Matrigel aliquot and dilute 1:30 in cold DMEM/F12.
- 3 Add Matrigel solution to plates quickly (2 ml for a 6-well, 12 ml for a 10-cm, 32 ml for a 15-cm).
- 4 Incubate at 37°C for a minimum of 1 hour. Use plates before the medium evaporates or within two weeks.

## iPSC to cardiomyocyte (iPSC-CM) differentiation with WTC11 cell line

- 5 Aspirate media and add dPBS to wash. Then, aspirate dPBS. The differentiation protocol is based on previously published work. **(1)**
- 6 Add Accutase<sup>TM</sup> to cells to dissociate them (1 ml for a 6-well, 3 ml for a 10-cm). Incubate at 37°C for 3 to 5 minutes.
- 7 Add 1:1 warm DMEM/F12 to the Accutase<sup>TM</sup> and pipette gently 2 to 3 times to promote a single-cell suspension. Transfer cells to a conical tube and centrifuge at 300g for 5 minutes.
- 8 Aspirate medium and resuspend cells in mTeSR+<sup>TM</sup> (mTeSR+<sup>TM</sup> basal medium + supplement) + Rock inhibitor (10uM final concentration).
- 9 Seed cells in Matrigel coated plates quickly so that in 72hrs cells are 70% to 80% confluent. (2ml for 6-well, 12ml for a 10-cm of mTeSR+<sup>TM</sup>)
- 10 Change medium to RB- (RPMI 1640 + 50x B-27<sup>TM</sup> supplement - insulin) + CHIR99021 (10 µM final concentration) (3 ml for a 6-well, 18 ml for a 10-cm). Volumes per tissue culture plate will remain consistent during differentiation.
- 11 After 48 hours, change medium to RB- + IWP2 (7.5 µM final concentration). Change media as close to 48 hours as possible, within +/- 5 minutes.



- 12 After another 48 hours, change medium to RB-. Change media as close to 48 hours as possible, within +/- 5 minutes.
- 13 After another 48 hours, change medium to RB+ (RPMI 1640 + 50x B-27™ supplement with 100 U/ml penicillin and 100 µg/ml streptomycin). Change media as close to 48 hours as possible, within +/- 5 minutes.
- 14 After another 48 hours, change medium to RB+. Change media as close to 48 hours as possible, within +/- 5 minutes. Check for cell beating.
- 15 After another 48 hours, check for cell beating. If beating has not occurred, replace with RB+ until beating occurs. If beating has occurred, change medium to NG+ (RPMI 1640 - glucose + 50x B-27™ supplement with 100 U/ml penicillin, 100 µg/ml streptomycin, 500 µg/ml recombinant human albumin, 213 ng/ml ascorbic acid, and 0.748 µl/ml 60% w/w sodium lactate solution), based on CDM3 medium. **(2)**
- 16 After another 48 hours, change medium to NG+.
- 17 After 48 hours, replate cells in an appropriate dish and density ( $3 \times 10^6$  cells for a 6-well,  $1.8 \times 10^7$  cells for a 10-cm). To do this, aspirate media, add dPBS to wash, and then aspirate dPBS. Differentiation is complete at this point. Cells are now referred to as iPSC-CMs.
- 18 Add trypsin-EDTA (0.05%) to cells to dissociate them (1 ml for a 6-well, 3 ml for a 10-cm). Incubate at 37°C for 3 minutes.
- 19 Add a 1:1 mixture of stop media (DMEM + 5% FBS + DNase at 20 µg/ml final concentration) to quench the trypsin-EDTA (0.05%). Pipette gently to dislodge cells from the tissue culture plate and form a single-cell suspension, limiting to 3 to 5 repetitions. Transfer cells to a conical tube and centrifuge at 300g for 5 minutes.
- 20 Resuspend cells in RB+ + Rock inhibitor (10 µM final concentration). Seed cells in a 6-well plate over six wells ( $3 \times 10^6$  cells per well).
- 21 After 24 hours, change medium to RB+ (3 ml for a 6-well). Repeat every 48 hours for 144 hours.

## Induce proliferation and developmental phenotypes in iPSC-CMs

- 22 Change the medium to RB+ for half the wells. For the remaining wells, change to RB+ + CHIR99021 (4 µM final concentration). Repeat every 48 hours for 168 hours.



## RNA-seq

- 23 iPSC-CMs were cultured for 28 days following the initiation of differentiation.
- 24 Isolate RNA using the Total RNA Purification Plus Kit (Norgen) and submit to Azenta for standard RNA-seq with ribosomal RNA depletion.
- 25 Trim reads using Trimmomatic v0.32 to remove adapters, and then align to GRCh38 using the STAR aligner v2.4.1a. **(3,4)**
- 26 Obtain gene counts with featureCounts from the subread package (version 1.4.6-p4) using the comprehensive gene annotation in Gencode v22. **(5)**
- 27 Conduct differential expression analysis with DESeq2, where gene counts are fitted to a negative binomial generalized linear model, and a Wald test identifies significantly differentially expressed genes. **(6)**

## Protocol references

- (1)** Jiang, L., Yin, M., Wei, X., Liu, J., Wang, X., Niu, C., Kang, X., Xu, J., Zhou, Z., Sun, S., Wang, X., Zheng, X., Duan, S., Yao, K., Qian, R., Sun, N., Chen, A., Wang, R., Zhang, J., ... Meng, D. (2015). Bach1 represses Wnt/ $\beta$ -catenin signaling and angiogenesis. *Circulation Research*, *117*(4), 364–375. <https://doi.org/10.1161/CIRCRESAHA.115.306829>
- (2)** Burridge, P. W., Matsa, E., Shukla, P., Lin, Z. C., Churko, J. M., Ebert, A. D., Lan, F., Diecke, S., Huber, B., Mordwinkin, N. M., Plews, J. R., Abilez, O. J., Cui, B., Gold, J. D., & Wu, J. C. (2014). Chemically defined generation of human cardiomyocytes. *Nature Methods*, *11*(8), 855–860. <https://doi.org/10.1038/nMeth.2999>
- (3)** Bolger, A. M., Lohse, M., & Usadel, B. (2014). Trimmomatic: A flexible trimmer for Illumina sequence data. *Bioinformatics*, *30*(15), 2114–2120. <https://doi.org/10.1093/bioinformatics/btu170>
- (4)** Dobin, A., Davis, C. A., Schlesinger, F., Drenkow, J., Zaleski, C., Jha, S., Batut, P., Chaisson, M., & Gingeras, T. R. (2013). STAR: Ultrafast universal RNA-seq aligner. *Bioinformatics*, *29*(1), 15–21. <https://doi.org/10.1093/bioinformatics/bts635>
- (5)** Liao, Y., Smyth, G. K., & Shi, W. (2014). FeatureCounts: An efficient general purpose program for assigning sequence reads to genomic features. *Bioinformatics*, *30*(7), 923–930. <https://doi.org/10.1093/bioinformatics/btt656>
- (6)** Love, M. I., Huber, W., & Anders, S. (2014). Moderated estimation of fold change and dispersion for RNA-seq data with DESeq2. *Genome Biology*, *15*(12). <https://doi.org/10.1186/s13059-014-0550-8>