

Oct 29, 2025

RNA Sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvw42k7lmk/v1>

Dahlia Rohm^{1,2}, Joshua Black^{1,2}, Sean R. McCutcheon^{1,2}, Alejandro Barrera^{1,2}, Shante Berry^{1,2}, Daniel Morone^{1,2}, Xander Nuttle^{3,4}, Celine de Esch^{3,4}, Derek Tai^{3,4}, Michael Talkowski^{3,4}, Nahid Iglesias¹, Charles A. Gersbach^{1,2}

¹Duke University; ²Center for Advanced Genomic Technologies; ³Massachusetts General Hospital;

⁴Broad Institute of Massachusetts Institute of Technology and Harvard University

Gersbach Lab



Andrea R Daniel

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvw42k7lmk/v1>

Protocol Citation: Dahlia Rohm, Joshua Black, Sean R. McCutcheon, Alejandro Barrera, Shante Berry, Daniel Morone, Xander Nuttle, Celine de Esch, Derek Tai, Michael Talkowski, Nahid Iglesias, Charles A. Gersbach 2025. RNA Sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvw42k7lmk/v1>



Manuscript citation:

Activation of the imprinted Prader-Willi Syndrome locus by CRISPR-based epigenome editing.

Rohm D, Black JB, McCutcheon SR, Barrera A, Berry SS, Morone DJ, Nuttle X, de Esch CE, Tai DJC, Talkowski ME, Iglesias N, Gersbach CA. *Cell Genomics*. 2025 Feb 12;5(2):100770. doi: 10.1016/j.xgen.2025.100770. PMID: 39947136

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 08, 2025

Last Modified: October 29, 2025

Protocol Integer ID: 229362

Keywords: RNA Sequencing, iPSC, Prader-Willi Syndrome, induced pluripotent stem cell, pluripotent stem cell, methods for rna, rna

Funders Acknowledgements:

NIH

Grant ID: HG012053

Disclaimer

This protocol was adapted from the work of Dahlia Rohm and colleagues in the Gersbach lab at Duke University.

Abstract

This protocols describes methods for RNA sequencing in human induced pluripotent stem cells.

Materials

Norgen Total RNA Purification Plus Micro kit (Norgen, 48500); Qubit RNA Broad Range quantification assay (Thermo, Q10211); Agilent TapeStation (Agilent TapeStation 4200) RNA ScreenTape (Agilent, 5067-5576); Genewiz (Azenta Life Sciences); Trimmomatic v0.32; STAR aligner; featureCounts; DESeq2.



Cell preparation and RNA preparation

- 1 Harvest cells, protocol defined by user. Pellet ~100,000 cells (between 1000 and 500,000) and proceed to RNA isolation immediately or flash-freeze in liquid nitrogen and store pellets at -80°C until processing.
- 0.2 Isolate total RNA using Norgen Total RNA Purification Plus Micro kit (Norgen, 48500).

RNA sequencing

- 0.3 Quantify RNA with the Qubit RNA Broad Range quantification assay (Thermo, Q10211) and optionally validate RNA integrity with Agilent TapeStation (Agilent TapeStation 4200) RNA ScreenTape (Agilent, 5067-5576) or a similar automated electrophoresis system.
- 0.4 Submit total RNA to Genewiz (Azenta Life Sciences) for total RNA sequencing for 2×150bp paired-end short read sequencing. Genewiz verifies quality of samples and libraries and performed DNase treatment as necessary.
- 0.5 Subject resulting reads to adapter trimming using Trimmomatic v0.32, align to GRCh37 or GRCh38 with the STAR aligner, and retrieve gene counts using featureCounts.
- 0.6 Perform differential expression analysis using DESeq2. Determine genes with significant differential expression using a threshold of $FDR < 0.01$ and absolute value of $\text{Log}_2(FC) > 1$.