

Jun 28, 2023

RNAseH-based ribodepletion

DOI

dx.doi.org/10.17504/protocols.io.36wgqj9n5vk5/v1

Nicole Schonrock¹, Helaine Grazielle Santos Vieira¹, Oguzhan Begik², Eva Maria Novoa³

¹Garvan Institute for Medical Research; ²Centre for Genomic Regulation; ³Center for Genomic Regulation (CRG)



Eva Maria Novoa

Center for Genomic Regulation (CRG)

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36wgqj9n5vk5/v1>

Protocol Citation: Nicole Schonrock, Helaine Grazielle Santos Vieira, Oguzhan Begik, Eva Maria Novoa 2023. RNAseH-based ribodepletion. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqj9n5vk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 12, 2022



Last Modified: June 28, 2023

Protocol Integer ID: 69866

Keywords: using rnaseh nuclease, rnaseh enzyme, rnaseh nuclease, based ribodepletion goal, rnaseh, ribodepletion goal, rna sample, rna molecule, dna oligonucleotide, complementary to human rrna, rna, approaches as nuclease, dna oligo, samples from rrna, rrna molecule, ribodepletion, human rrna, many rounds of ribodepletion, ribodepleted material, nuclease, rrna, 28s ribodepletion, dna, enzyme, dna duplex

Disclaimer

Do not use for long read approaches.

This protocol is adapted from the protocol used in Adiconis et al., Nature Methods 2013
(<https://doi.org/10.1038/nmeth.2483>)

Abstract

Goal: Method to ribodeplete samples from rRNAs using RNaseH nuclease.

Summary: DNA oligos are complementary to human rRNAs. They will anneal to the rRNA molecules. RNaseH will degrade RNA-DNA duplexes, leaving the RNA samples depleted from rRNAs

Tested on : Human, Mouse and Yeast (for 18s and 28s ribodepletion).

Cost: DNA oligonucleotides will be expensive the first time ordered, but will last for many rounds of ribodepletion. Once the DNA oligonucleotides have been purchased, the major cost comes from the RNaseH enzyme. In our hands, the cost is much lower than commercial kits, thus being scalable for high input samples where the ribodepleted material needs to be high.

Recommended for: NGS-based approaches. Do not use for Nanopore-based approaches as nuclease-based methods degrade RNA molecules.

Troubleshooting

