

Oct 10, 2024

## ROS Measurement Using CellROX

DOI

<https://dx.doi.org/10.17504/protocols.io.yxmvme82og3p/v1>

Ali Ghoochani<sup>1,2,3,4</sup>, Monther Abu-Remaileh<sup>1,2,3,4</sup>

<sup>1</sup>Department of Chemical Engineering, Stanford University, Stanford, CA 94305, USA;

<sup>2</sup>Department of Genetics, Stanford University, Stanford, CA 94305, USA;

<sup>3</sup>The Institute for Chemistry, Engineering and Medicine for Human Health (Sarafan ChEM-H), Stanford University, Stanford, CA 94305, USA;

<sup>4</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA

asap



Ali Ghoochani

Stanford University

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.yxmvme82og3p/v1>

**External link:** <https://www.thermofisher.com/document-connect/document-connect.html?url=https://assets.thermofisher.com/TFS-Assets%2FLSG%2Fmanuals%2Fmp10422.pdf>



**Protocol Citation:** Ali Ghoochani, Monther Abu-Remaileh 2024. ROS Measurement Using CellROX. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.yxmvme82og3p/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 10, 2024

**Last Modified:** October 10, 2024

**Protocol Integer ID:** 109560

**Keywords:** CellRox, Reactive oxygen species measurement, oxidative stress levels in cell, fluorescent probe cellrox, oxidative stress level, using cellrox, flow cytometry, cellrox this protocol, ros measurement, cell, cytometry

## Abstract

This protocol describes a method to assess oxidative stress levels in cells using the fluorescent probe CellROX using flow cytometry

## Materials

- **CellROX-Green reagent (5  $\mu$ M)** (Invitrogen, Cat. No. C10444)
- **PBS**
- **Trypsin**
- **Flow cytometer**



## Preparation for Experiment

- 1 Once cells are ready for the experiment, aspirate the media and add fresh media

## CellROX Treatment

30m

- 2 Treat the cells with 5  $\mu$ M CellROX-Green (Invitrogen, Cat No. C10444) for  00:30:00 at  37 °C

30m

- 2.1 Wash the cells once with PBS.

## Cell Detachment

3m

- 3 Trypsinize the cells to detach them from the culture surface.  00:03:00  37 °C

3m

- 3.1 Centrifuge the cells at  1000 rpm  00:03:00 and discard the supernatant.

3m

## Flow Cytometry Analysis

- 4 Resuspend the cells in PBS.
- 5 Measure the intensity of CellROX fluorescence using flow cytometry to assess oxidative stress levels.

## Protocol references

### CellROX