

Sep 10, 2020

Version 2

## Rosbash/Janelia StickLAMP Protocol V.2

DOI

[dx.doi.org/10.17504/protocols.io.bk89kzz6](https://dx.doi.org/10.17504/protocols.io.bk89kzz6)



Albert Yu<sup>1</sup>, Timothy A Brown<sup>2</sup>, Jasmine Quynh Le<sup>1</sup>, Kristina Galatsis<sup>1</sup>, Michael Rosbash<sup>1</sup>

<sup>1</sup>Brandeis University; <sup>2</sup>Janelia Research Campus

XPRIZE Rapid Covid Tes...



Albert Yu

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bk89kzz6>

**Protocol Citation:** Albert Yu, Timothy A Brown, Jasmine Quynh Le, Kristina Galatsis, Michael Rosbash 2020. Rosbash/Janelia StickLAMP Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bk89kzz6> Version created by **Albert Yu**

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 10, 2020

**Last Modified:** September 10, 2020

**Protocol Integer ID:** 41985

**Keywords:** saliva sample, saliva contaminant, carryover of saliva contaminant, saliva,  $\mu$ l of saliva, rapid purification step, nucleic acid amplification product, rna, faithful detection of sar, selective separation of bead, novel high contrast dye, colorimetric readout


## Abstract


A protocol for the detection of SARS-CoV-2 from saliva samples featuring a rapid purification step and a high-contrast colorimetric readout. Saliva is first inactivated using a 100x inactivation reagent consisting of 2.5M TCEP, 100 mM EDTA, 1.2N NaOH solution diluted to approximately 1x final concentration and heated to 95C for 5 minutes. RNA is rapidly purified and concentrated with magnetic beads in a PEG/NaCl-based buffer using a 3D-printed magnetic stick that enables selective separation of beads without carryover of saliva contaminants. Beads are eluted directly into an RT-LAMP reaction mix, which uses a novel high contrast dye that turns from purple to clear when acidified by nucleic acid amplification products that enables unambiguous identification of successful amplification. This protocol is sensitive down to 1 copy/ $\mu$ l of SARS-CoV-2 in 300  $\mu$ l of saliva. This degree of sensitivity enables faithful detection of SARS-CoV-2 even in pooled samples.


## Materials


### MATERIALS


 NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**

 Twist synthetic SARS-CoV-2 RNA control **Twist Bioscience Catalog #Mt007544.1**

 SARS-CoV-2 Master Mix

 Actin Master Mix


 100x Inactivation Reagent

 Bead Mix

 Magnetic Tips

 Heat Block at 65C


 Heat Block at 95C

 Magnetic Stick

### STEP MATERIALS

-

-

 Magnetic Tips

 NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**

 SARS-CoV-2 Master Mix

 Actin Master Mix

 Water

 Twist synthetic SARS-CoV-2 RNA control **Twist Bioscience Catalog #Mt007544.1**

100x Inactivation Reagent

2.5M TCEP

150mM EDTA

1.2N NaOH

SARS-CoV-2/Actin Master Mix

12.5µl SARS-CoV-2/Actin Buffer/Dye/Primer Mix (Currently only available from us)

0.5µl WarmStart RTx NEB M0380L

1µl Bst2.0 NEB M0537L

11µl H2O


Primers used


See [https://docs.google.com/spreadsheets/d/11n-9754VqtsXszTC2tUxFq-\\_gKIGgjPL-KtKqXVevH4/edit#gid=0](https://docs.google.com/spreadsheets/d/11n-9754VqtsXszTC2tUxFq-_gKIGgjPL-KtKqXVevH4/edit#gid=0)


Bead Mix


See [https://ethanomics.files.wordpress.com/2012/08/serapure\\_v2-2.pdf](https://ethanomics.files.wordpress.com/2012/08/serapure_v2-2.pdf) with 300µl beads instead of 1000µl


## Protocol materials


 Actin Master Mix

 Twist synthetic SARS-CoV-2 RNA control **Twist Bioscience Catalog #Mt007544.1**

 SARS-CoV-2 Master Mix

 Twist synthetic SARS-CoV-2 RNA control **Twist Bioscience Catalog #Mt007544.1**

 Magnetic Tips


 100x Inactivation Reagent

 Heat Block at 65C


 NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**

 SARS-CoV-2 Master Mix

 Water

 Actin Master Mix

 NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**

 Bead Mix

 Magnetic Tips


 Heat Block at 95C

 Magnetic Stick

 100x Inactivation Reagent

 SARS-CoV-2 Master Mix

 Actin Master Mix

 Bead Mix

 Magnetic Tips

 NaCl **Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014**

 Water

 Twist synthetic SARS-CoV-2 RNA control **Twist Bioscience Catalog #Mt007544.1**

## Troubleshooting



## Safety warnings

! Do not open up PCR tubes after amplification.

## Before start

Prepare:

Saliva collection kit (2.0ml Tube and funnel provided by us, or your own saliva collection device from standard labware, such as 1.5ml, 5ml, 15ml, or 50ml tubes. Saliva samples >1ml will likely have to be subsampled)

Magnetic stick

1 magnetic tip per sample

Bead mix: Let bead mix come to room temperature for 20 minutes prior to use, and ensure beads are suspended in solution by vortexing or pipetting up and down



130mM NaCl

SARS-CoV-2 Master mix

Actin Master mix







## Saliva Collection

- 1 Instruct patient to avoid food, drink, toothbrushing, and nasal sprays for a minimum of  00:30:00 prior to sample collection
- 2 Begin pooling saliva in your mouth. Saliva production can be stimulated by thinking about food, or about the saliva collection itself.
- 3 Gently expel saliva into the funnel, tapping to collect in the tube, until amount of saliva is approximately flush with the base of the funnel  750  $\mu$ L Approximately

## Inactivation

5m

- 4 Add inactivation reagent to approximately 1x final concentration. Reaction is tolerant of between 0.7x to 2x final concentration.  7.5  $\mu$ L Approximately  
 100x Inactivation Reagent
- 5 Invert 40 times to mix
- 6 Heat tube to approximately  95  $^{\circ}$ C for  00:05:00 . Viral RNA release is similar between 93-98C. Use tube clip to prevent popping.

5m

## Equipment

ThermoMixer

NAME

Benchtop Incubator

TYPE

Eppendorf

BRAND

5382000023

SKU

<https://online-shop.eppendorf.us/US-en/Temperature-Control-and-Mixing-44518/Instruments-44519/Eppendorf-ThermoMixerC-PF-19703.html>


LINK

Any heat block will suffice


SPECIFICATIONS




7 Remove tube from heat and let rest at  Room temperature for at least



 00:03:00


OR


 On ice

for at least


 00:00:30 .

8 While tube is resting, aliquot  25  $\mu$ L SARS-CoV-2 mastermix and  25  $\mu$ L Actin mastermix to separate wells of PCR strip tube, 96-well plate, or 1.5ml tube per sample

 On ice .

Per run, prepare two additional  25  $\mu$ L SARS-CoV-2 mastermixes for positive and negative controls.

 SARS-CoV-2 Master Mix


 Actin Master Mix


## STEP CASE

## If pooling

From 1 to 10 steps

Prepare one 25ul SARS-CoV-2 reaction and one 25ul Actin reaction per 5 samples

- 9 Add approximately 0.7x volumes of bead mix -  525  $\mu$ L Approximately . Sample is tolerant of between 0.7x-1.2x volumes of bead mix. Pipette up and down to mix.






 Bead Mix

## STEP CASE

## If pooling

9 steps

Remove 60ul of inactivated saliva from 5 samples and add to a single tube, for a total of 300 $\mu$ l. Add 210 $\mu$ l bead mix to pooled tube.

- 10 Let stand at  Room temperature for  00:03:00 10m
- 11 Cap magnetic stick with a clean tip and dip in bead/sample mix for  00:02:00 , 2m  
dipping up and down 5 times every  00:00:30 . Meanwhile, prepare  500  $\mu$ L  
130mM NaCl in a separate 1.5ml or 2ml tube.





## Equipment

### Magnetic Stick

NAME

Rosbash/Brown

BRAND

None

SKU

Magnetic stick used for bead purifications

SPECIFICATIONS

Magnetic Tips

NaCl Merck MilliporeSigma (Sigma-Aldrich) Catalog #53014

- 12 Remove magnetic stick from sample and swirl in clean 130mM NaCl solution for 00:00:05 . Discard NaCl solution.

5s

- 13 Remove magnetic stick from wash sample and place in SARS-CoV-2 mix for 00:00:30

30s

- 14 Remove magnetic stick from SARS-CoV-2 mix and place in Actin mix for 00:00:30

30s

- 15 Add 5  $\mu$ L water to additional SARS-CoV-2 Mix (negative control) and 5  $\mu$ L synthetic Twist SARS-CoV-2 positive RNA control to additional SARS-CoV-2 Mix, prepared in Step 8.

Water

Twist synthetic SARS-CoV-2 RNA control Twist Bioscience Catalog #Mt007544.1

16 Cap tubes and place on 🌡️ 65 °C heating apparatus for 🕒 00:40:00 .

40m

If using a thermal cycler, run with the following program:

65C for 40 minutes

4C indefinitely

#### Equipment

ThermoMixer

NAME

Benchtop Incubator

TYPE

Eppendorf

BRAND

5382000023

SKU

<https://online-shop.eppendorf.us/US-en/Temperature-Control-and-Mixing-44518/Instruments-44519/Eppendorf-ThermoMixerC-PF-19703.html>

LINK

Any heat block will suffice

SPECIFICATIONS



17 Remove tubes from heating apparatus and examine color change.

Positive		Negative		Inconclusive			
SARS-CoV-2	Actin	SARS-CoV-2	Actin	SARS-CoV-2	Actin		
							
SARS-CoV-2	Actin			Negative Control			
							
				Positive Control			
							

- 18 If a positive sample is found when pooling, re-test pooled samples individually.