

Nov 19, 2025

Version 6

RSVAB WGS and GF protocols V.6

DOI

<https://dx.doi.org/10.17504/protocols.io.kqdg3xbzqg25/v6>

Maria Iglesias¹

¹Laboratory of Reference and Research in Respiratory Viruses, National Centre for Microbiology, Instituto de Salud Carlos III, 28220 Majadahonda, Spain



Maria Iglesias

National Institute of Health Carlos III

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kqdg3xbzqg25/v6>

Protocol Citation: Maria Iglesias 2025. RSVAB WGS and GF protocols. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.kqdg3xbzqg25/v6> Version created by [Maria Iglesias](#)

Manuscript citation:

Generic novel system for genomic characterization of Respiratory Syncytial Virus obtaining whole genome sequencing and a full-length G and F sequences.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 19, 2025

Last Modified: November 19, 2025

Protocol Integer ID: 232920

Keywords: novel rsv amplicon, viral genome, cdna from rsv, rsvab wg, viral nucleic acid extract, rsv, specific sequences of the main antigen, systems for genomic characterization, generating cdna, producing amplicon, genomic characterization, main antigen

Abstract

This SOP describes the procedure for generating cDNA from RSV viral nucleic acid extracts and subsequently producing amplicons tiling the viral genome using. We propose two systems for genomic characterization of RSV. First, a novel RSV amplicon-based system for WGS, and second, a method focused on obtaining the specific sequences of the main antigens, G and F.



Sample management

- 1 No prior treatment is necessary. However, in samples with high Cts, the use of DNaseI in the sample as a pretreatment can improve coverage.

dsCDNA generation:

- 2 During this step three master mixes will be prepared: MMI, MMII and MMIII.

Materials: **Kit Superscript III First Strand (Invitrogen)**
100% DMSO
RNAseH (Invitrogen)
Klenow fragment 3' → 5' exo (New England Biolabs)
Primer FR26RV-N : 5'GCC GGA GCT CTG CAG ATA TCNNNNNN 3'

Note

This step must be performed in a RNase free, pre-PCR environment in which post PCR RSV amplicons are not present, to minimise risk of sample contamination.

Citation

Díez-Fuertes F, Iglesias-Caballero M, García-Pérez J, Monzón S, Jiménez P, Varona S, Cuesta I, Zaballos Á, Jiménez M, Checa L, Pozo F, Pérez-Olmeda M, Thomson MM, Alcamí J, Casas I
(2021). A Founder Effect Led Early SARS-CoV-2 Transmission in Spain.

<https://doi.org/10.1128/JVI.01583-20>

LINK

- 3 **MMI Preparation:**

	A	B
	FR26RV-N (10uM)	2
	DMSO	0,5
	Total	2,5 ul

Mix thoroughly by vortexing.

4 **MMII Preparation:**

	A	B
	10x First Strand Buffer	2
	DTT 100 mM	2
	MgCl ₂ 25mM	4
	dNTPs	1
	RNaseOUT	0,5
	SSIII RT	0,5
	Total	10 ul

Kit Superscript III First Strand (Invitrogen)

5 **MMIII Preparation:**

	A	B
	Klenow 5'-3'	1
	RNaseH	0,5
	Total	1,5 ul

6 Defrost extracted RNA.



Maintain **on ice** the MMI,MMII and MMIII mixes.

7 **MMI Amplification:**

Add  5 μL  Sample in MMI mix

Place the tube on a thermocycler and run the following program:

A	B
65°C	5 min
4°C	2 min

Briefly tube centrifugation

8 **MMII Amplification:**

Addition of  10 μL from MMII in the tube with the MMI and the viral extraction.

Place the tube on a thermocycler and run the following program:

A	B
25°C	10 min
50°C	50 min
85°C	10 min
4°C	∞

Briefly tube centrifugation

9 **MMIII Amplification:**

Addition of  1.5 μL of MMIII into the tube with the previous mixes and the viral extraction

Place the tube on a thermocycler and run the following program:

A	B
37°C	60 min

A	B
75°C	15 min

Briefly tube centrifugation

- 10 **STOP POINT:** cDNA can be stored at 4°C (same day) or -20°C (up to a week).

RSVAB WGS protocol

11

Citation

Iglesias-Caballero M, Camarero-Serrano S, Varona S, Mas V, Calvo C, García ML, García-Costa J, Vázquez-Morón S, Monzón S, Campoy A, Cuesta I, Pozo F, Casas I (2023)

. Genomic characterisation of respiratory syncytial virus: a novel system for whole genome sequencing and full-length G and F gene sequences.

<https://doi.org/10.2807/1560-7917.ES.2023.28.49.2300637>

LINK

Materials:

2x MyTaqRed mix (Bioline)

Primers:

A	B
Primer ID	Sequence (5'-3')
Mix A	
RSVCombinitia	ACGC GAAA AAAT GCGT



A	B
	ACWA CA
RSVWGS4R	CATG WTG WYTT ATTT GCCC C
RSVWGS2F	CACT WACA ATAT GGGT GCC
RSVWGS2R.2	CRTT YCTT AARG TRGG CC
RSVGF-R	TTCG YGAC ATATT TGCC CC
RSVWGS3.2F	ACAT GGAA AGAY ATYA GCC
RSVWGS3.2R	TTGC ATCT GTAG CAGG AATG G
RSVWGS6F	TATAY AGAT AYCA YATG GGTG G
RSVWGS14R	CGAT ATAA CAAR TTRG GATC ACC
RSVWGS13F	AAAA GATT GGGG



A	B
	AGAG GG
RSVWGS13R	GGRC CTAT TGTA AGGA CTAR GT
Mix B	
RSVWGS9F	GARC AACT CAAA GAAA ATGG
RSVWGS9R	AYTG RAAC ATRG GCAC CC
OG1-21	GGGG CAAA TGCA ACCA TGTC C
RSVWGS5F	CATA ATTA YTTT GAAT GGC
RSVWGS5R.2	CAAA CATT TAAT CTRC TRAG GCTR A
RSVWGS14F. 2	TYAA CAAC YCKA ATCA YGTG G
RSVWGS6R	CCCT CTCC CCAA TCTT TTTC

A	B
RSVWGS13F	AAAA GATT GGGG AGAG GG
RSVWGS3.2R	TTGC ATCT GTAG CAGG AATG G
RSVCombEnd ing	ACGA GAAA AAAA GTGT CAAA AACT AA

Primer description by mix

	1	NC_038235.1	1	22	RSVC ombin ital	1	+
	2	NC_038235.1	2202	2221	RSVW GS9F	1	+
	3	NC_038235.1	2329	2348	RSVW GS4R	2	-
	4	NC_038235.1	3352	3370	RSVW GS2F	1	+
	5	NC_038235.1	3364	3381	RSVW GS9R	2	-
	6	NC_038235.1	4672	4692	OG1- 21	1	+



	7	NC_038235.1	7593	7611	RSVG F-R	2	-
	8	NC_038235.1	7674	7692	RSVW GS5F	1	+
	9	NC_038235.1	9304	9322	RSVW GS3.2 F	1	+
	10	NC_038235.1	9317	9341	RSVW GS5R.2	2	-
	11	NC_038235.1	9932	9949	RSVW GS2R.2	2	-
	12	NC_038235.1	10313	10333	RSVW GS14F.2	1	+
	13	NC_038235.1	10799	10820	RSVW GS6F	1	+
	14	NC_038235.1	13036	13055	RSVW GS6R	2	-
	15	NC_038235.1	13038	13055	RSVW GS13F	1	+
	16	NC_038235.1	11455	11477	RSVW GS14 R	2	-
	17	NC_038235.1	14213	14233	RSVW GS3.2 R	2	-
	18	NC_038235.1	14508	14531	RSVW GS13R	2	-
	19	NC_038235.1	15192	15217	RSVC ombE nding	2	-
	1	NC_001781.1	1	22	RSVC ombin itial	1	+
	2	NC_001781.1	2202	2221	RSVW GS9F	1	+



	3	NC_001781.1	2331	2350	RSVWGS4R	2	-
	4	NC_001781.1	3354	3372	RSVWGS2F	1	+
	5	NC_001781.1	3366	3383	RSVWGS9R	2	-
	6	NC_001781.1	4673	4693	OG1-21	1	+
	7	NC_001781.1	7607	7624	RSVG F-R	2	-
	8	NC_001781.1	7687	7705	RSVWGS5F	1	+
	9	NC_001781.1	9316	9334	RSVWGS3.2F	1	+
	10	NC_001781.1	9329	9353	RSVWGS5R.2	2	-
	11	NC_001781.1	9944	9961	RSVWGS2R.2	2	-
	12	NC_001781.1	10325	10345	RSVWGS14F.2	1	+
	13	NC_001781.1	10811	10832	RSVWGS6F	1	+
	14	NC_001781.1	13048	13067	RSVWGS6R	2	-
	15	NC_001781.1	13050	13067	RSVWGS13F	1	+
	16	NC_001781.1	11467	11489	RSVWGS14R	2	-
	17	NC_001781.1	14225	14245	RSVWGS3.2R	2	-

18	NC_001781.1	14522	14543	RSVWGS13R	2	-	
19	NC_001781.1	15197	15222	RSVC ombE nding	2	-	

Primer scheme

Note

The protocol is based in the RSV genome amplification in two separate mixes with an unique amplification program. The mixes that will be mixed at the end of cycling.

12 Preparation of RSV Amplification Mix A:

A	B
Mezcla A	1x (μl)
MyTaqRed 2x	15
H2O	7,8
RSV CombInit (10μM)	0,2
WGS 4R (10μM)	0,2
WGS 2F (10μM)	0,2
WGS 2R (10μM)	0,2
RSV GF_R (10μM)	0,2



	A	B
	WGS 3.2F (10µM)	0,2
	WGS 3.2R (10µM)	0,2
	WGS 6F (10µM)	0,2
	WGS 14R (10µM)	0,2
	WGS 13F (10µM)	0,2
	WGS 13R (10µM)	0,2
	Volumen total	25

13 Preparation of RSV Amplification Mix B:

	A	B
	Mezcla B	1x (µl)
	MyTaqRed 2x	15
	H2O	8
	WGS 9F (10µM)	0,2
	WGS 9R (10µM)	0,2
	RSV OG1_21 (10µM)	0,2




A	B
WGS 5F (10μM)	0,2
WGS 5R.2 (10μM)	0,2
WGS 14F.2 (10μM)	0,2
WGS 6R (10μM)	0,2
WGS 13F (10μM)	0,2
WGS 3.2R (10μM)	0,2
CombEnd (10μM)	0,2
Volumen total	25

14 Addition of  5 μL of the previous prepared double stranded cDNA on each mix.

15 **Amplification protocol:**

A	B	C
95°C	1 min	
95°C	30 seg	x45
55°C	7 min	
72°C	2 min	
72°C	5 min	
12°C	∞	

16 To assess PCR performance, the amplicons can be loaded onto a 1% agarose gel for electrophoresis.



1-4: <22
5-8: 22-28
9-12: >30

- 17 Finally, mix in one single tube both mixes and proceed to purification and library preparation.

RSVAB GF protocol starting from viral extraction

- 18 Materials:
Qiagen OneStep RT-PCR kit.

- 19 **Preparation of GF amplification mix:**

	A	B
	H2O	10
	5xQ PCR MM	6
	dNTPs	1
	RSVWGS_2F (10uM)	1
	RSVGF-R (10 uM)	1
	RT-PCR mix	1
	Total	20 ul

20 Addition of  10 µL of the viral extraction

21 **Amplification protocol GF:**

	A	B	C
	48°C	60 min	
	95°C	15 min	
	95°C	30 seg	x 45
	55°C	2 min	
	72°C	1 min	
	72°C	5 min	
	12°C	∞	

Citations

Step 2

Díez-Fuertes F, Iglesias-Caballero M, García-Pérez J, Monzón S, Jiménez P, Varona S, Cuesta I, Zaballos Á, Jiménez M, Checa L, Pozo F, Pérez-Olmeda M, Thomson MM, Alcamí J, Casas I. A Founder Effect Led Early SARS-CoV-2 Transmission in Spain.

<https://doi.org/10.1128/JVI.01583-20>