


May 07, 2025

RT-LAMP for BYV detection

DOI

<https://dx.doi.org/10.17504/protocols.io.261geeo1og47/v1>



Marta Niedzicka¹, Stephen Byrne¹

¹Teagasc



Marta Niedzicka

Teagasc

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261geeo1og47/v1>

Protocol Citation: Marta Niedzicka, Stephen Byrne 2025. RT-LAMP for BYV detection. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261geeo1og47/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 02, 2025

Last Modified: May 07, 2025

Protocol Integer ID: 125978

Keywords: beet yellows virus, BYV, RT-LAMP, plant virus, detecting beet yellows virus, symptomatic beet samples from ireland, byv genome, lamp for byv detection, symptomatic beet sample, byv detection, virus, byv, rna, easy colorimetric assay, yellow color, total rna, assembled genome, genome

Funders Acknowledgements:

Horizon Europe - Marie Skłodowska-Curie Actions

Grant ID: 101106728

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This protocol describes the process of detecting beet yellows virus (BYV) in plant material utilising a quick and easy colorimetric assay. The starting point is total RNA extracted from the material. The positive results are shown in yellow color, while negative samples stay bright pink - an example of these results is provided at the end of the protocol. The basis for development of the assay was assembly of a BYV genome from symptomatic beet samples from Ireland - the sequence data and assembled genome have been deposited in a public repository.

Materials

 WarmStart® Colorimetric LAMP 2X Master Mix **New England Biolabs Catalog #M1800**

 Nuclease-free water **Merck MilliporeSigma (Sigma-Aldrich)**



Protocol materials

☒ WarmStart® Colorimetric LAMP 2X Master Mix **New England Biolabs Catalog #M1800**

☒ Nuclease-free water **Merck MilliporeSigma (Sigma-Aldrich)**

Safety warnings

⚠ Please read SDS associated with various consumables and kits used in this protocol and wear appropriate PPE. A site specific procedural risk assessment should be carried out prior to introducing this protocol to the lab.

Before start

This protocol assumes that you have already isolated Total RNA from plant material.

Primers

- 1 Primers were designed using NEB LAMP Primer design tool: <https://lamp.neb.com/#/>
The BYV genome sequence used for the design is deposited in GenBank (accession no. PV593736)

	A	B
	BYV_F3_P243L1	AGTCGATGAAAGAGTGATCG
	BYV_B3_P243L1	AGCAGTGTCAAACCTTCTCG
	BYV_FIP_P243L1	ATGAATAGCACACAGGGCGGTGACGCATTTGGTTCACAG
	BYV_BIP_P243L1	TTGTACCCTAAAGTTAAGGGTTGGAACAGTTCGCGAATTTGC
	BYV_LF_P243L1	TCACCACTTTCGTACCCGTACTTA
	BYV_LB_P243L1	TACGGTAAACTGAAATTTGTGCTGC

Designed primers within ORF1a/b gene of BYV

gBlock sequence ordered and used for testing the primers:

CAAAGTTTTTCGTTATGACCGGTCACGACGTTTTCTTCGTGCCTGACCCTTATAAACTTTTA
GTGAAATTGGGAGCTTCTAAGGATGAAGTGGACGATGAGTTTCTGTTTGAAGTGTTCACCTC
TTTTCGCGATTTAACGAAAGATTTAGTCGATGAAAGAGTGATCGAACTCTTGACGCATTTGG
TTCACAGTAAGTACGGGTACGAAAGTGGTGACACGTACGCCGCCCTGTGTGCTATTCATTGT
ATTCGTTCAAACTTTTTCATCGTTCAAGAAATTGTACCCTAAAGTTAAGGGTTGGGTCGTTCA
CTACGGTAAACTGAAATTTGTGCTGCGCAAATTCGCGAACTGTTTTCGCGAGAAGTTTGAC
ACTGCTTTCGGCGAAGCGTACTTTCTTACTTACGACGAAGCTTGAGACTGTGTGGTAACTC
GGTTGTTGTTTTTGTGTTGTTGTGTGTCTGTAGTCCTACTTCGCGGTGATGGACTGTGTACT
CCGCTC

Primer mix

- 2 We recommend making stockings of the LAMP primers at a useable concentrations.
Here, we suggest a 10X Primer Mix containing all LAMP primers.



	A	B	C
	Primer	10X Concentration (Stock)	1X Concentration (Final)
	FIP	16 µl	1.6 µl
	BIP	16 µl	1.6 µl
	F3	2 µl	0.2 µl
	B3	2 µl	0.2 µl
	Loop F (LF)	4 µl	0.4 µl
	Loop B (LB)	4 µl	0.4 µl

LAMP Primer Mix

- 3 Prepare primer stocks in nuclease free water and store at -20 °C for up to 2 years.

Note: Make primer stock in molecular biology graded water rather than other buffer to avoid carryover of additional buffer to the LAMP reaction

Reaction mix

40m

- 4 Thaw all components to be used at Room temperature and place on ice. Vortex briefly to mix thoroughly. Centrifuge to collect material and place On ice .

Components:

- Target nucleic acid
- LAMP Primer 10X
- WarmStart® Colorimetric LAMP 2X Master Mix **New England Biolabs Catalog #M1800**
- Nuclease-free water **Merck MilliporeSigma (Sigma-Aldrich)**

- 5 Prepare reaction mix as described below using **Colorimetric LAMP MasterMix**, **LAMP primers** and **nuclease-free water** (1-3 positions from the table below) for your number of samples. The final volume with target included per sample is 25 µL



A	B
WarmStart Colorimetric LAMP 2X Master Mix	12.5 µl
LAMP Primer Mix 10X	2.5 µl
dH ₂ O	9 µl
Target (DNA, RNA, or water for NTC)	1 µl
Total volume	25 µl

Reaction mix per 1 sample - final concentration

- 6 Vortex reaction mix and centrifuge to collect material.
- 7 Pipet  24 µL per reaction into desired reaction vessels and add  1 µL of sample.
Mix by vortexing or by pipetting if using a plate or similar vessel, centrifuge to collect if necessary.
Check that reaction solutions have a bright pink color, which indicates initial high pH required for successful pH-LAMP reaction.
- 8 Seal reaction vessels.
- 9 Incubate at  65 °C for  00:30:00 .

30m
- 10 Remove tubes or vessels from incubation and examine by eye. Positive reactions will have turned yellow while negative controls should remain pink. If color change is not robust, e.g. an orange color is visible, return reactions to  65 °C for an additional  00:10:00 . Reactions can be examined earlier if desired, and high copy of input reactions can exhibit full color change in as little as 10–15 minutes.

10m


Color will be visible directly on removal from incubation temperature, but can be intensified by allowing reaction to cool to room temperature.

Results

40m

- 11 The result can be photographed or scanned to record the colorimetric results, or simply kept at room temperature in the reaction vessel (the colour change remains stable).

Example result after  00:30:00 incubation



Example RT-LAMP results: 6 positive samples, negative sample (RNA from barley infected with BYDV) and NTC