


Apr 17, 2020

RTPCR Amplification of SARS-CoV2 Whole Genome for Illumina NGS

DOI

<https://dx.doi.org/10.17504/protocols.io.bew8jfhw>



Monica Galiano¹, Shahjahan Miah², Angie Lackenby², Omolola Akinbami², Joanna Ellis², Maria Zambon²

¹WHO Influenza Collaborating Centre, The Francis Crick Institute;

²Respiratory Virus Unit, Microbiology Services Colindale, Public Health England

Coronavirus Method De...



Monica Galiano

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bew8jfhw>

Protocol Citation: Monica Galiano, Shahjahan Miah, Angie Lackenby, Omolola Akinbami, Joanna Ellis, Maria Zambon 2020. RTPCR Amplification of SARS-CoV2 Whole Genome for Illumina NGS . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bew8jfhw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 10, 2020

Last Modified: April 17, 2020

Protocol Integer ID: 35520

Keywords: whole genome, Illumina, SARS-COV2, WGS SARS-COV2, COVID19, nCOV, nCOV-2019, cov2 whole genome for illumina ng, cov2 whole genome, whole genomes from sar, cov genome, rtPCR amplification of sar, obtaining full genome sequence, cov2 positive clinical sample, rtPCR amplification, performing rtPCR amplification, full genome sequence, PCR amplification, reverse transcription, respiratory virus unit, cov2 positive sample, phe with rdrp gene detection, whole genome, different locations on the genome, rdrp gene detection, viral nucleic acid, genome, contamination of the sequence, end of the genome, using random hexamer primer, random hexamer primer, full genome, transcription, redundancy in whole genome, separate PCR reaction, integrity of viral nucleic acid, bioinformatic assembly pipeline, PCR reaction, specific step on the bioinformatic assembly pipeline, overlapping amplicon, cov2, cov2 detection protocol, first sar, low viral load, location of primer, PCR, direct comparison of the sequence, rtPCR, positive clinic,

Abstract

SUMMARY

This document describes the procedure for performing RTPCR amplification for next-generation sequencing (NGS) of whole genomes from SARS-CoV2 positive clinical samples. This methodology was employed at the Respiratory Virus Unit, Microbiology Services Colindale, Public Health England (RVU-PHE) to sequence the first SARS-CoV2 positive samples.

Reverse transcription (RT) is performed using random hexamer primers. PCR amplification is done using sequence-specific primers which amplify two sets of 30-31 overlapping amplicons (size 1.0 to 1.8 kb), each set independently covering the entire length of the MERS-CoV genome. Once amplicons are obtained, rough equimolar mixtures, clean-up, DNA quantitation and final dilutions are performed to prepare samples for library preparation for Illumina MiSeq next-generation sequencing (NGS).

Primers from set A do not overlap location of primers from set B, except for those binding at the end of the genome (see figure 1). Amplicon mixtures from each set are treated as separate samples throughout the procedure, including assembly. Merging (or comparison) of consensus sequences of both sets is performed at the end.

IMPORTANT: Please note that every amplicon is obtained in a separate PCR reaction and, if using both primer sets, there is a total of ~61 PCR reactions; this protocol is NOT a multiplex approach. We recommend using a 96-well plate per sample.

ADVANTAGES AND LIMITATIONS

The use of both primer sets allows for redundancy in whole genome (WG) coverage, which increases the chances of obtaining full genome sequences when samples have low viral load. Although we haven't formally determined the sensitivity of this assay (or limit of amplification of WG), we have obtained full coverage with samples with CT values of up to 38 (CT value obtained at RVU-PHE with RdRp gene detection within the SARS-CoV2 detection protocol as described by Corman et al. [1]). Variation in storage or transport conditions of samples may affect the integrity of viral nucleic acids.

Another advantage of this protocol is that, as primers from set A and B match different locations on the genome, direct comparison of the sequences obtained with each set allows for curation of primer-induced contamination of the sequences, although we deal with this with a specific step on the bioinformatic assembly pipeline (not described here).

Alternatively, when samples have a relatively good CT (<30) only set A may be used, simplifying the amplification and library steps.

The main limitation of this protocol is its unsuitability for high throughput processing of large batches of samples, although current efforts are being made to adapt the primers to a multiplex approach.

Materials

EQUIPMENT

PCR thermocycler for 96-well plates
Vortex
Pipettes with disposable filter tips
Ice bucket
Microcentrifuge and picofuge
Plate centrifuge
Mother E-base
96-well plates with adhesive seals
1.5ml and 0.2ml polypropylene tubes/strip tubes
Benchtop UV transilluminator (UVP)
Qubit Fluorometer (Life Technologies)
Qubit Assay Tubes (Life Technologies)
FrameStar FastPlate 96 (4titude)
Aluminum StarSeal (PCR) E2796-9792, StarLab

REAGENTS

SuperScript™ VILO™ cDNA Synthesis Kit (ThermoFisher, Invitrogen cat 11754050)
dNTP set 100mM (ThermoFisher, Invitrogen cat 10297018)
Oligonucleotide primers 100pmol/μl stocks (Eurofins) (see Appendix 3)
Platinum Taq DNA polymerase HiFi (ThermoFisher, Invitrogen cat 11304029)
Nuclease-free water (Severn Biotech)
1% or 2% E-gels 48-wells (Invitrogen)
10x TBE buffer
DNA 1Kb plus ladder (Invitrogen cat 10787018)
BlueJuice 10X Gel loading buffer (ThermoFisher, Invitrogen cat10816015)
Qubit dsDNA HS Assay Kit (ThermoFisher, Invitrogen Q32854)
QIAquick PCR purification kit (cat 28106)

Safety warnings



Sample inactivation in lysis buffer should be performed at BSL3



- 1 A total amount of 140 μl of RNA (includes some excess) are needed to set up the necessary volume of RT reactions. RNA should be kept on ice or stored at $-80\text{ }^{\circ}\text{C}$ if not immediately used.
- 2 The RT reaction is as follows (Table 1):

RT mix	x 1 (μl)	Bulk mix (x 14)
Water	4.0	56
5X VILO Reaction mix	4.0	56
10X Superscript Enzyme mix	2.0	28
RNA	10.0	140
Final volume	20.0	280

Table 1

- 3 To ensure uniform heat distribution, we aliquot this mix into 7×0.2 (strip) tubes.
Cycling program:

$25\text{ }^{\circ}\text{C}$ 00:10:00

$42\text{ }^{\circ}\text{C}$ 01:00:00

$85\text{ }^{\circ}\text{C}$ 00:05:00

Hold at $4\text{ }^{\circ}\text{C}$

Keep cDNA on ice or store at $-80\text{ }^{\circ}\text{C}$ if not immediately used.

- 4 Combine primers according to Tables 1 and 2 (4.1) and dilute as follows to get a working primer dilution of $5\text{ pmol}/\mu\text{l}$

Note

Preparation of primer mix pairs (working dilution 5 pmol/μl)

Use different colours of strip tubes, e.g. yellow strip tubes for SET A and purple strip tubes for SET B. Add 180 μl of water + 10 μl of each Forward and Reverse primers at 100 pmol/μl, according to Appendix 2. Vortex to mix. Store at -20°C.

Thaw both sets of primer mixes, briefly vortex them and spin down. Using a multichannel pipette (volume 0.5 or 1 to 10 μl) load 4 μl of primer mixes to the corresponding wells of a 96-well PCR plate, according to Figure 2. Keep the PCR plate on a cooler.

Figure 2. SET A (yellow) and SET B (purple)

	1	2	3	4	5	6	7	8	9	10	11	12
A	1A	2A	3A	4A	5A	6A	7A	8A	9A	10A	11A	12A
B	13A	14A	15A	16A	17A	18A	19A	20A	21A	22A	23A	24A
C	25A	26A	27A	28A	29A	30A	31A					
D												
E	1B	2B	3B	4B	5B	6B	7B	8B	9B	10B	11B	12B
F	13B	14B	15B	16B	17B	18B	19B	20B	21B	22B	23B	24B
G	25B	26B	27B	28B	29B	30B						
H												

4.1

The position of primers on the nCoV genome sequence has been based on the sequence of strain hCoV-19/Wuhan/IVDC-HB-01/2019 (GISAID accession EPI_ISL_402119). We gratefully acknowledge the authors (Wenjie Tan et al), originating and submitting laboratory (National Institute for Viral Disease Control and Prevention, China CDC) of this sequence from GISAID's EpiCOV™ Database on which this protocol is based.

	Primer pair	Forward primer	Sequence	Position of 5'end	Reverse primer	Sequence	Position of 5'end	Amplification size
--	-------------	----------------	----------	-------------------	----------------	----------	-------------------	--------------------



1A	ncov-1F	ATTAAAGGTTTAT ACCTTCCCAGGT AAC	1	ncov-1032R	TCAAAAGGTG TCTGCAATTC AT	1032	1031
2A	ncov-15F	CCTTCCCAGGTA ACAAACCA	15	ncov-1032R	TCAAAAGGTG TCTGCAATTC AT	1032	1017
3A	ncov-693F	TTGACTTAGGCG ACGAGCTT	693	ncov-2030R	CAGATGTGAA CATCATAGCAT CAA	2030	133 7
4A	ncov-1684F	CGCCATTATTTT GGCATCTT	168 4	ncov-3019R	TGAAGCCAAT TTAAACTCAC CA	3019	133 5
5A	ncov-2665F	TGCCCTTGCACC TAATATGAT	266 5	ncov-4010R	GGAAGTCTAGT TTCTTCCAGA GTTG	4010	134 5
6A	ncov-3686F	CACGAAGTTCTA CTTGACCA	368 6	ncov-5015R	CAAGTTGCGT GTGGAGGTTA	5015	132 9
7A	ncov-4667F	TCTCTCAAAGTG CCAGCTACAG	466 7	ncov-6037R	GCTTGCCTTT GGATATGGTT	6037	137 0
8A	ncov-5687F	TCAGCACCACT GCTCAGTA	568 7	ncov-7021R	AACACCTAAA GCAGCGGTTG	7021	133 4
9A	ncov-6688F	TGCTAAGCCTTT TCTTAACAAAGT T	668 8	ncov-8027R	CTGCAACTTC CGCACTATCA	8027	133 9
10A	ncov-7663F	TGATGAAGTTGC GAGAGACTTG	766 3	ncov-9019R	TTCAGCAGCC AAAACACAAG	9019	135 6
11A	ncov-8700F	TTGATGGTGGTG TCACTCGT	870 0	ncov-10033R	AGAGGTTTGT GGTGGTTGGT	1003 3	133 3
12A	ncov-9699F	TCATTTGTATTT CACAAAGCA	969 9	ncov-11071R	CAAAGACCAT TGAGTACTCT GGA	11071	137 2
13A	ncov-10696F	TGGAGACAGGTG GTTTCTCA	106 96	ncov-12030R	TGCATGGAAA GCAAAACAGA	1203 0	133 4
14A	ncov-11668F	CCGCTACTTTAG ACTGACTCTTGG	1166 8	ncov-13050R	GCAGGCACTT CTGTTGCAT	1305 0	138 2
15A	ncov-12701F	CCTGTTGCACTA CGACAGATG	1270 1	ncov-14043R	AATACCAGCAT TTCGCATGG	1404 3	134 2



16 A	ncov-13663 F	CACACTTTCTCT AACTACCAACAT GAA	136 63	ncov-15028 R	ATGCGAAAAG TGCATCTTGA	1502 8	136 5
17A	ncov-14663 F	AAACTGTCAAAC CCGGTAATTTT	146 63	ncov-16023 R	AGACACGAAC CGTTCAATCA	1602 3	136 0
18 A	ncov-15686 F	TGCGTAAACATT TCTCAATGATG	1568 6	ncov-17038 R	TTGCAACATT GCTAGAAAAC TCA	1703 8	135 2
19 A	ncov-16689 F	TGCTACTGTACG TGAAGTGCTG	166 89	ncov-18034 R	AAGTTGCCAC ATTCCTACGT G	1803 4	134 5
20 A	ncov-17694 F	TGCAATTAACAG GCCACAAA	1769 4	ncov-19031 R	TCGTGAAGAA CTGGGAATTT G	1903 1	133 7
21A	ncov-18683 F	CATGCTTTTCCA CTGCTTCA	186 83	ncov-20043 R	AACACCATTA CGGGCATTTC	2004 3	136 0
22 A	ncov-19664 F	TTGATGGACAAC AGGGTGAA	196 64	ncov-21231 R	GTCCACCATG CGAAGTGTC	2123 1	156 7
23 A	ncov-20863 F	CATTTTGGTGCT GGTTCTGA	208 63	ncov-22226 R	CGAAAAACCC TGAGGGAGAT	2222 6	136 3
24 A	ncov-21896 F	TTCGAAGACCCA GTCCCTAC	2189 6	ncov-23214 R	CACCTGTGCC TGTAAACCA	2321 4	1318
25 A	ncov-22883 F	TCTTGATTCTAA GGTTGGTGGT	228 83	ncov-24231 R	CACCAAAGGT CCAACCAGAA	2423 1	134 8
26 A	ncov-23850 F	TTAAACCGTGCT TTAACTGGAATA	238 50	ncov-25243 R	ATGGCAATCA AGCCAGCTAT	2524 3	139 3
27 A	ncov-24858 F	GGCACACACTGG TTTGTAACAC	248 58	ncov-26224 R	TGCTTACAAA GGCACGCTAG T	2622 4	136 6
28 A	ncov-25886 F	TCTTCAATTGTC ATTACTTCAGGT G	258 86	ncov-27227 R	CCTGAAAGTC AACGAGATGA AA	2722 7	1341
29 A	ncov-26889 F	GCCACTCCATGG CACTATTC	268 89	ncov-28191 R	TTCATAGAAC GAACAACGCA CT	2819 1	130 2
30 A	ncov-27876 F	TTGTCACGCCTA AACGAACA	278 76	ncov-29226 R	ACATTCCGAA GAACGCTGAA	2922 6	135 0

31 A	ncov- 28871 F	AGGCAGCAGTAG GGGAAGCTT	288 71	ncov- 29848 R	AAAATCACAT GGGGATAGCA	2984 8	977
---------	---------------------	---------------------------	-----------	---------------------	--------------------------	-----------	-----

Table 2. Set A: Primer sequences & combinations

Primer pair	Forward primer	Sequence	Position of 5'end	Reverse primer	Sequence	Position of 5'end	Amplification size
1B	ncov-193F	CTTACGGTTT CGTCCGTGTT	193	ncov-1541R	GTGGAACCCAAT AGGCACAC	1541	1348
2B	ncov-1185F	CGTCACCAAA TGAATGCAAC	1185	ncov-2495R	TTCTGTGGGAAG TGTTTCTCC	2495	1310
3B	ncov-2193F	GAGACGGTTG GGAAATTGTT	2193	ncov-3523R	TTCAACTTGCAT GGCATTGT	3523	1330
4B	ncov-3172F	TCAACCTGAA GAAGAGCAAG AA	3172	ncov-4525R	AGCACCATAATC AACCACACC	4525	1353
5B	ncov-4173F	CTAAAAAGGC TGGTGGCACT	4173	ncov-5550R	AGGGTTGTCTGC TGTTGTCC	5550	1377
6B	ncov-5192F	CTGGGTAGGT ACATGTCAGC A	5192	ncov-6548R	GATCTGTGTGGC CAACCTCT	6548	1356
7B	ncov-6192F	CACCCTCTTT TAAGAAAGGA GCTA	6192	ncov-7537R	TGTACATTGAC TCTTGTTGCTC	7537	1345
8B	ncov-7194F	CCATTTTCATC TTTTAAATGG GATT	7194	ncov-8554R	CCACCCTTAAGT GCTATCTTTG	8554	1360
9B	ncov-8179F	GCAAGGGTTT GTTGATTGAG A	8179	ncov-9506R	TAAAGGCAACTA CATGACTGTATT CAC	9506	1327
10B	ncov-9169F	CCTTGAAGGT TCTGTTAGAG TGG	9169	ncov-10456R	GAAATTGGGCCT CATAGCAC	10456	1287
11B	ncov-10080F	CATCTGGTAA AGTTGAGGGT TGT	10080	ncov-11533R	TCTGGCCAAAA ACATGACAG	11533	1453



12B	ncov-11189F	CCTTCTCTTG CCACTGTAGC TT	11189	ncov-12547R	TGCTGATGCATA AGTAAATGTTG	12547	1358
13B	ncov-12163F	GGCTGTTGCT AATGGTGATT C	12163	ncov-13547R	TCAAAAGCCCT GTATACGACA	13547	1384
14B	ncov-13191F	CAGTTACACC GGAAGCCAAT	13191	ncov-14513R	TCCTGATTATGT ACAACACCTAGC TC	14513	1322
15B	ncov-14183F	CCAGGGCTTT AACTGCAGAG	14183	ncov-15544R	CCGTGACAGCTT GACAAATG	15544	1361
16B	ncov-15167F	TGAAATCAATA GCCGCCACT	15167	ncov-16511R	AAACCAAAAAC TTGTCCATTAGC	16511	1344
17B	ncov-16176F	TTCAAGGTAT TGGGAACCTG A	16176	ncov-17551R	CGAGGAACATGT CTGGACCTA	17551	1375
18B	ncov-17171F	CCGCTGTTGA TGCACTATGT	17171	ncov-18533R	TTTATACGCACT ACATTCCAAGG	18533	1362
19B	ncov-18165F	ACCTGGCATA CCTAAGGACA	18165	ncov-19542R	CATCATGTTATAA GCATCGAGATAC A	19542	1377
20B	ncov-19191F	TTGGAATTGC AATGTGCGATA G	19191	ncov-20478R	ACCTGTTTGCGC ATCTGTTA	20478	1287
21B	ncov-19966F	TGTGCACCAC TCACTGTCTT	19966	ncov-21730R	AAGAACAAGTCC TGAGTTGAATG	21730	1764
22B	ncov-21360F	CAAATCCAAT TCAGTTGTCT TCC	21360	ncov-22705R	CCATAACACTTA AAAGTGGAATA GA	22705	1345
23B	ncov-22363F	TGGGTTATCT TCAACCTAGG ACTT	22363	ncov-23713R	TTTGTGGGTATG GCAATAGAGTT	23713	1350
24B	ncov-23384F	GGTTGCTGTT CTTTATCAGG ATG	23384	ncov-24735R	GAGGTGCTGACT GAGGGAAG	24735	1351
25B	ncov-24369F	GACTCACTTT CTTCCACAGC AA	24369	ncov-25712R	AGAGAAAAGGG GCTTCAAGG	25712	1343
26B	ncov-25354F	TGCTCAAAGG AGTCAAATTA CATT	25354	ncov-26734R	AAACAGCAGCA AGCACAAAA	26734	1380

27 B	ncov- 26361 F	TGTGTGCGTA CTGCTGCAA	2636 1	ncov- 27712 R	TGCCGCAACAAT AAGAAAAA	2771 2	1351
28 B	ncov- 27361F	TGAAGAGCAA CCAATGGAGA	2736 1	ncov- 28719 R	GGGTGCCAATGT GATCTTTT	2871 9	135 8
29 B	ncov- 28397 F	GCCCCAAGGT TTACCCAATA	2839 7	ncov- 29734 R	GGTGAATATGTG GTGGCTCT	2973 4	133 7
30 B	ncov- 28397 F	GCCCCAAGGT TTACCCAATA	2839 7	ncov- 29871 R	TGATTTTAATAG CTTCTTAGGAGA ATGAC	2987 1	147 4

Table 3. Set B: Primer sequences & combinations

- 5 Prepare a PCR mix following *Table 4*. Thaw and briefly vortex all the required reagents except enzymes. Prepare a mix for ~70 samples (enough for both sets A and B).

Note

Preparation of working dilution of dNTPs (100 mM dNTP Set, Invitrogen, cat 10297018)
Mix the four tubes of dNTPs in one bijou tube. Add 1.5 ml of water to reach concentration of 10 mM (Final volume 2.5 ml). Aliquot in 2 ml tubes with screw cap.

PCR mix	X 1 (μl)	X 70
Water	35.8	2,506
10X High Fidelity Buffer (Invitrogen)	5.0	350
10 mM dNTPs mix (Invitrogen) See Appendix 1	1.0	70
50mM MgSO ₄ (Invitrogen)	2.0	140
Platinum Taq DNA Polymerase 5 U/μl (Invitrogen)	0.2	14
Final volume	44.0	3,080

Table 4



- 6 Mix with a gentle vortex. Keep on ice. Dispense 44µl of the PCR mix into every well with primers of rows A, B, C, E, F and G.
- 7 Pull together all 7 cDNA strip tubes into one and load 4 µl of cDNA in each well.
- 8 Cover the 96-well plate with an adhesive lid or strip caps. Spin briefly the plate.

- 9 Transfer the 96-well plate to a thermal cycler.
Amplify using the following cycling conditions:

95 °C 00:10:00

followed by 35 cycles of:

94 °C 00:00:30

52 °C 00:00:30

68 °C 00:05:00

Hold at 15 °C

10 Working in a post-PCR area

Use 1% or 2% E-gels 48 wells (Invitrogen) for visualisation of the bands, following manufacturer's instructions. For each sample prepare 1× 96-well microtiter plate to mix Loading buffer

- 11 Prepare diluted Blue Juice loading buffer. Add 10 µl of diluted Blue Juice to every well of rows A, B, C, E, F and G (mirroring the PCR plate) Use multichannel pipettes.

Note

Preparation of Blue Juice diluted 1/25 or 1/50

Mix 10 ml of buffer TBE with 400 µl (1/25) or 200 µl (1/50) µl of undiluted Blue Juice.

- 12 Add 5 µl of sample from the PCR plate to the microtiter plate. Use multichannel pipettes.
- 13 Separately, prepare enough mixture of diluted Blue Juice (14µl diluted Blue Juice + 1µl of 1 Kb Plus Ladder) to add to wells marked with M.



- 14 Transfer 15 μ l of samples and ladder dilutions into the E-gels. Load 15 μ l of diluted Blue Juice in empty wells (no wells should be empty). Run for 20 minutes and view under shortwave U.V.


Note

Alternative using multichannel pipette (tips will fit into every other well of 48 wells E-gels): if two samples are run in parallel, load amplicons from sample A in odd wells, and amplicons from sample 2 in even wells. Repeat the procedure with rows B to F. Load 15 μ l of the mixture of Blue Juice + Ladder into the wells marked with M at the beginning and end of each gel row. In that way, amplicon 1 of sample A will run next to amplicon 1 of sample B and similarly for each amplicon.

- 14.1 Some examples of gel pictures are shown here:

Some gel pictures (4 samples)

CT values were obtained with the RdRp gene detection described by Corman et al. (1)



- 15 To prepare equimolar mixes, ideally PCR reactions should be first cleaned up using e.g. silica columns or magnetic beads and then quantified, but the amount of PCR reactions makes this very impractical, unless using automated methods. Instead, rough equimolar mixes are prepared based on brightness of gel bands: 3 µl of a strong band, 5 µl of a

normal band, and 10 µl of a weak or non-visible band. Mix the appropriate amount of each amplicon from SET A in an Eppendorf tube (final volume 60 to 300µl). Do not include the negative control. Repeat the procedure for SET B amplicons in a separate tube.

- 16 Perform PCR product purification using QIAquick PCR purification kit columns, following manufacturer's instructions. Perform repeated additions and spins of the sample until all of it has gone through the column. Elute the purified DNA in 50 µl of molecular grade water.
- 17 Run quantitation of DNA using the Qubit dsDNA HS reagent (Life Technologies) and the Qubit fluorometer, according to manufacturer's instructions.
- 18 Dilute the samples to the appropriate concentration for Illumina library preparation using Nextera XT (not described here).